


# Total PFAS analysis using Combustion Ion Chromatography

*Theresa Steurer, Metrohm AG, Switzerland*

# Fluorine in focus



Lanciki, A. Metrohm AG (2021): Adsorbable organic fluorine (AOF), WP-078EN

- US EPA lists >12'000 chemicals as per- and polyfluorinated alkyl substances (PFAS)  
(Source: [CompTox Chemicals Dashboard \(epa.gov\)](https://www.epa.gov/comp-tox-chemicals) )
- PFAS are persistent, ubiquitous, and bioaccumulative
- PFAS are used in industrial and consumer products
- Harmful to human health
- Monitoring individual substances difficult to impossible
- Sum parameters to detect organic fluorine compounds

# What to expect...

Why to use CIC for AOF analysis and why a dedicated sample preparation for AOF is needed.

Introduction to IC and Combustion IC.

Sample preparation and analysis of AOF.

# Why use CIC for AOF analysis?

- AOX (adsorbable organic halogens, Cl, Br, I) combustion + microcoulometric titration after adsorption step (DIN 9562, EPA 1650)
- DIN 38409-59: AOF, AOCl, AOBr, and AOI after adsorption by CIC
- Improved and for AOF **optimized sample preparation** (neutralization)
- Superiority
  - Adsorption = preconcentration
  - Dedicated adsorption columns
  - IC – sensitive ion analysis
  - Beside TOP assay and EOF – 3<sup>rd</sup> approach for total PFASs assessment

## Why not AOF?

- Titration with  $\text{AgNO}_3$  – no precipitation

## What got improved?

- Sample preparation for AOF
- Detection of individual halogen-species




# Why dedicated sample preparation for AOF?

## REMOVAL OF INORGANIC FLUORINE


- AOX: Sample and washing solution at  $\text{pH} < 2$ 
  - Mostly non-dissociated HF
  - can adsorb to activated carbon
- AOF: Sample and washing solution at neutral pH
  - Inorganic fluorine can be removed






# Introduction to IC and to Combustion IC

# Ion chromatography in brief



- Analytical technique
- Separation of (chemically) similar substances in one run
- Specialized for anions, cations and polar substances
- Automated sample preparation as valuable feature

# Theory of CIC – a sample preparation technique





# Whole combustion process

## Setup for solid samples



930 IC Compact Flex  
(IC analysis)



920 Absorber Module  
(absorption of combustion  
gases and liquid handling)



Combustion Oven  
(hydropyrolysis of samples)



ABD with MMS 5000  
(sampling of solids)

- Sample is weighted in on boat
- Boat is picked up from rack by gripper and transferred into CIC tube
- Combustion happens according to light dispersion or fixed boat program
- Absorption in water or hydrogen peroxide after combustion
- Measurement of absorber solution by IC

# Samples for Combustion IC

- Combustible samples
- Samples can be in liquid, solid or gaseous form
- Samples can be analyzed for **halogens** (F, Cl, Br, I) and **sulfur** simultaneously
- Up and coming technique in norms and standards, especially in petroleum industry
- Breaks up organics which usually cannot be directly analyzed by IC






# Sample preparation and analysis of AOF


# Don't be sour - neutral is the solution

Sample preparation and preconcentration for AOF analysis (EDIN 38409-59)



Sample/blank  
neutralization  
with sodium  
nitrate

1



Adsorption step and  
rinsing  
(activated carbon  
columns for AOF)


2






Semi automation with *APU sim*  
(Analytik Jena)




# Complexity of sample preparation and sample analysis for AOF determination



# CIC analysis of a surface water sample




# Calculation of the AOF concentration from IC analysis



fluoride sample  
( $\mu\text{g/L}$ )

column 2  
 $\sum$   
column 1

AOF ( $\mu\text{g/L}$ )



fluoride blank  
( $\mu\text{g/L}$ )


Volume  
absorption  
solution (mL)

Volume for  
adsorption on  
activated carbon  
(100 mL)

$$\sum_{\text{column 1}}^{\text{column 2}} c(\text{AOF}) = \left( c(\text{F}^-)_{\text{IC}} * \frac{V_{\text{Abs}}}{V_{\text{Smpl}}} \right) - \left( c(\text{F}_{\text{BW}}^-)_{\text{IC}} * \frac{V_{\text{AbsBW}}}{V_{\text{SmplBW}}} \right)$$

# Replicate AOF measures in “natural” samples:

Excellent repeatability and performance for AOF analysis by CIC despite the complexity of this method



| AOF #1<br>( $\mu\text{g/L}$ ) | AOF #2<br>( $\mu\text{g/L}$ ) | AOF #3<br>( $\mu\text{g/L}$ ) | AOF #4<br>( $\mu\text{g/L}$ ) | AOF <sub>AVG</sub> $\pm$<br>SD ( $\mu\text{g/L}$ ) | RSD<br>(%) |
|-------------------------------|-------------------------------|-------------------------------|-------------------------------|--|------------|
| 6.3                           | 6.3                           | 6.8                           | 6.8                           | 6.5 $\pm$ 0.3                                      | 5          |
| 10.2                          | 10.0                          | 9.3                           | 9.2                           | 9.7 $\pm$ 0.5                                      | 5          |
| 7.4                           | 7.0                           | 7.6                           | 7.2                           | 7.3 $\pm$ 0.3                                      | 4          |

Excellent repeatability with RSDs $\leq$ 5%

Data  
evaluation

Combustion, absorption, and  
IC analysis

AOF  
adsorption



A person wearing a full-body white protective suit, a white face mask, and blue gloves is standing in a shallow, polluted stream. They are holding a clear plastic vial in their right hand and a silver metal case in their left hand. The background shows a forested area with a large, rusted metal pipe or structure partially submerged in the water. The water is murky and greenish-brown, indicating pollution. The overall scene suggests an environmental investigation or sampling of contaminated water.

WILL AOF ALSO BE YOUR  
FUTURE SOLUTION FOR PFAS  
SCREENING?