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Supersedes EN ISO 10304-1:1995, EN ISO 10304-2:1996

English Version

Water quality - Determination of dissolved anions by liquid chromatography of ions - Part 1: Determination of bromide, chloride, fluoride, nitrate, nitrite, phosphate and sulfate (ISO 10304-1:2007)

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METHOD 300.0

DETERMINATION OF INORGANIC ANIONS BY ION CHROMATOGRAPHY

1.0 SCOPE AND APPLICATION

1.1 This method covers the determination of the following inorganic anions:

PART A.	
Bromide	Nitrite
Chloride	Ortho-Phosphate-P
Fluoride	Sulfate
Nitrate	

PART B.

Bromate	Chlorite
Chlorate	

1.2 The matrices applicable to each method are shown below:



1.2.1 Drinking water, surface water, mixed domestic and industrial wastewaters, groundwater, reagent waters, solids (after extraction 11.7), leachates (when no acetic acid is used).

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






ICS 13.060.50

English Version

Water quality - Determination of dissolved anions by liquid chromatography of ions - Part 1: Determination of bromide, chloride, fluoride, nitrite, phosphate and sulfate (ISO 10304-1:2007)

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This European Standard was approved by CEN on 1 March 2008.

CEN members are bound to comply with the CEN/CENELEC Internal Regulations which stipulate the conditions for giving this European Standard the status of a national standard without any alteration. Up-to-date lists and bibliographical references concerning such national standards may be obtained on application to the CEN/CENELEC Management Centre or to any CEN member.

This European Standard exists in three official versions (English, French, German). A version in any other language made by translation under the responsibility of a CEN member into its own language and notified to the CEN/CENELEC Management Centre has the same status as the official versions.

CEN members are the national standards bodies of Austria, Belgium, Bulgaria, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland and United Kingdom.

Management Centre: Avenue Marnix 17, B-1000 Brussels

© 2009 CEN All rights of exploitation in any form and by any means reserved worldwide for CEN national Members.

Ref. No. EN ISO 10304-1

EN ISO 10304-1

March 2009

Supersedes EN ISO 10304-1:1995, EN ISO 10304-2:1996

EUROPEAN STANDARD
NORME EUROPÉENNE
EUROPÄISCHE NORM

ICS 13.060.50

English Version

Water quality - Determination of dissolved anions by liquid chromatography of ions - Part 4: Determination of chlorate, chlorite and chlorite in water with low contamination (ISO 10304-4:2022)

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This European Standard was approved by CEN on 8 January 2022.

CEN members are bound to comply with the CEN/CENELEC Internal Regulations which stipulate the conditions for giving this European Standard the status of a national standard without any alteration. Up-to-date lists and bibliographical references concerning such national standards may be obtained on application to the CEN/CENELEC Management Centre or to any CEN member.

This European Standard exists in three official versions (English, French, German). A version in any other language made by translation under the responsibility of a CEN member into its own language and notified to the CEN/CENELEC Management Centre has the same status as the official versions.

CEN members are the national standards bodies of Austria, Belgium, Bulgaria, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Republic of North Macedonia, Romania, Serbia, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and United Kingdom.

CEN-CENELEC Management Centre: Rue de la Science 23, B-1040 Brussels

© 2022 CEN All rights of exploitation in any form and by any means reserved worldwide for CEN national Members.

Ref. No. EN ISO 10304-4

EN ISO 10304-4

March 2022

Supersedes EN ISO 10304-4:1999

METHOD 300.1

DETERMINATION OF INORGANIC ANIONS IN DRINKING WATER BY ION CHROMATOGRAPHY

1. SCOPE AND APPLICATION

1.1 This method covers the determination of the following inorganic anions in reagent water, surface water, ground water, and finished drinking water. As a result of different specified injection volumes (See conditions in Tables 1A and 1B), these anions are divided between the common anions listed in Part A and the inorganic disinfection by-products listed in Part B. These different injection volumes are required in order to compensate for the relative concentrations of these anions in drinking water and maintain good chromatographic peak shape throughout the expected dynamic range of the detector. Bromide is included in both Part A, due to its importance as a common anion, as well as Part B due to its critical role as a disinfection by-product precursor.


PART A.-- Common Anions	
Bromide	Nitrite
Chloride	ortho-Phosphate-P
Fluoride	Sulfate
Nitrate	
PART B.-- Inorganic Disinfection By-products	
Bromate	Chlorite
Bromide	Chlorate

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



* HAA9 sont MCAA, MBAA, DCAA, BCAA, DBAA, TCAA, BDCAA, CDBAA, TBAA

METHOD 557

DETERMINATION OF HALOACETIC ACIDS, BROMATE, AND DALAPON IN DRINKING WATER BY ION CHROMATOGRAPHY ELECTROSPRAY IONIZATION TANDEM MASS SPECTROMETRY (IC-ESI-MS/MS)

1. SCOPE AND APPLICATION

- 1.1 Method 557 is a direct-inject, ion chromatography, negative-ion electrospray ionization, tandem mass spectrometry (IC-ESI-MS/MS) method for the determination of haloacetic acids in finished drinking water. Bromate and dalapon (2,2-dichloropropionic acid) may be measured concurrently with the haloacetic acids. Real time, chromatographic separation of common anions in drinking water (matrix elimination) is a key feature of this method. Acceptable method performance has been demonstrated for matrix ion concentrations of 320 mg/L chloride, 250 mg/L sulfate, 150 mg/L bicarbonate and 20 mg/L nitrate. Method 557 requires the use of MS/MS in Multiple Reaction Monitoring (MRM) mode to enhance selectivity. Precision and accuracy data have been generated for the detection of nine haloacetic acids, bromate, and dalapon in reagent water, synthetic sample matrix, and finished drinking water from both ground water and surface water sources. The single laboratory Lowest Concentration Minimum Reporting Level (LCMRL) has also been determined in reagent water. Method 557 is applicable for the measurement of the following analytes:


Analyte	Chemical Abstracts Services Registry Number (CASRN)
Bromate (BrO_3^-)	15541-45-4 (BrO_3^- anion)
Bromochloroacetic acid (BCAA)	5589-96-8
Bromodichloroacetic acid (BDCAA)	71133-14-7
Chlorodibromoacetic acid (CDBAA)	5278-95-5
Dalapon	75-99-0
Dibromoacetic acid (DBAA)	631-64-1
Dichloroacetic acid (DCAA)	79-43-6
Monobromoacetic acid (MBAA)	79-08-3
Monochloroacetic acid (MCAA)	79-11-8
Tribromoacetic acid (TBAA)	75-96-7
Trichloroacetic acid (TCAA)	76-03-9

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

Polyphosphates

Phosphate, pyrophosphate, trimetaphosphate, tripolyphosphate

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