












# Electrodes for Titration

Which electrode for which application?

Application	Specifics	Electrode	Order no.
<b>Aqueous acid/ base titrations</b>	General	Ecotrode Plus	6.0262.100
		Ecotrode Gel	6.0221.100
	Alkaline samples, Bayer liquors. Titrations performed at high temperatures	Unitrode	6.0259.100
	Acidity of alcoholic beverages	Unitrode <i>easyClean</i> with Pt1000	6.0260.010
	Titration with small sample volumes	Flat membrane pH electrode	6.0256.100
	Titration in water having a low conductivity	Aquatrode Plus	6.0253.100
	Carbonate hardness, acid capacity of water, p- & m-value	Aquatrode Plus with Pt1000	6.0257.600
	Electroplating baths, etching baths, phosphate coating	Profitrode	6.0255.100
	Electroplating and etching baths containing fluoride or hydrofluoric acid	Solitrode HF	6.0223.100
<b>Non-aqueous acid/base titrations</b>	Protein-containing samples	Porotrode	6.0235.200
	Determination of the base number (TBN) of mineral oils, titrations with perchloric acid, cyclohexylamine, alcoholic HCl	Solvotrode <i>easyClean</i> with LiCl <sub>sat</sub> in ethanol	6.0229.010
<b>Redox titrations</b> Arsenite, cerium sulfate, iron(III), iodine, potassium bromate, sodium nitrite, oxalic acid, permanganate, thiosulfate, titan(III), Hg(NO <sub>3</sub> ) <sub>2</sub>	Determination of the acid number (TAN) of mineral oils, titrations with alcoholic KOH, NaOH and TBAOH, potassiummethylate, free acids; hydroxyl number of oils and fats	Solvotrode <i>easyClean</i> c(TEABr) = 0.4 mol/L in ethylene glycol	6.0229.010 (+ 6.2320.000)
	Titration at constant pH value	Pt-Titrode	6.0431.100
	Redox titrations at varying pH value	Combined Pt-ring-electrode	6.0451.100
	Chemical oxygen demand in water (ferrometry)	Combined Au-ring-electrode	6.0452.100
	Penicillin, ampicillin	Combined Au-ring-electrode	6.0452.100
	Titration in I <sub>pot</sub> mode	Double Pt sheet electrode	6.0309.100
	Bromatometry, iodometry, cerimetry according to Ph. Eur. & USP	Pt-Titrode	6.0431.100
<b>Karl Fischer Titrations</b>	Water content according to Karl Fischer	Double Pt-wire electrode (HF-resistant version on request)	6.0338.100
		Combined Ca <sup>2+</sup> -ISE polymer membrane	6.0510.100
<b>Complexometric</b> Titrations with EDTA, Complexon <sup>®</sup> III and IV	Back-titration of the Ba <sup>2+</sup> excess with EDTA	Combined Ca <sup>2+</sup> -ISE polymer membrane	6.0510.100
	Determination of Ca <sup>2+</sup> , Mg <sup>2+</sup> (Application Bulletin 125)	Cu <sup>2+</sup> ISE – crystal membrane	6.0502.140
	Determination of Al, Ba, Bi, Ca, Cd, Co, Fe, Mg, Ni, Pb, Zn (Application Bulletin 101)	Ag-Titrode	6.00430.100
<b>Precipitation titrations</b> Titrations with silver nitrate	Chloride in general, sodium chloride in food	Ag-Titrode with Ag <sub>2</sub> S-coating	6.00430.100S
	Chloride in dialysis and infusion solutions	Ag-Titrode with Ag <sub>2</sub> S-coating	6.00430.100S
	Titration according to Ph. Eur. & USP	Ag-Titrode with Ag <sub>2</sub> S-coating	6.00430.100S
	Determination of hydrogen sulfide, mercaptans, carbonyl sulfide, sulfides	Ag-Titrode with Ag <sub>2</sub> S-coating	6.00430.100S
	Chloride, bromide, iodide and cyanide in electroplating baths	F-ISE crystal membrane	6.0502.150
<b>Photometric titrations</b>	Fluoride/hydrofluoric acid in etching baths	Optrode	6.1115.000
	Titration in aqueous and non-aqueous media; choice of eight wave lengths (470, 502, 520, 574, 590, 610, 640, 660 nm)	Surfactrode Resistant	6.0507.130
<b>Surfactants in non-aqueous media</b> Aromatic/aliphatic, hydrocarbons, ketones, chloroethanes, methylisobutylketone	Titration of anionic and cationic surfactants, titrations in chloroform, surfactant formulations containing oils, pH <10	Surfactrode Refill	6.0507.140
	Titration of anionic and cationic surfactants, titration of washing agents, soaps, pH >10	Cationic Surfactant	6.0507.150
<b>Surfactants in aqueous media</b>	Titration of cationic surfactants	Ionic Surfactant	6.0507.120
	Titration of anionic surfactants	NIO electrode	6.0507.010
	Titration of non-ionic surfactants, titration of pharmaceutical ingredients with sodium tetraphenylborate	Thermoprobe	6.9011.020
<b>Thermometric titrations</b>	Titration in aqueous and non-aqueous solutions without HF	Thermoprobe HF	6.9011.040
	Titration in aqueous solutions containing HF		

# Practical tips, care and maintenance of electrodes for titration

Ecotrode Plus	Ecotrode Gel	Unitrode	Aquatrode Plus	Profitrode	Solvotrode easyClean	Ag/Pt/Au-Titrodes	Combined Ag/Pt/Au electrodes	Surfactrodes	Surfactant electrodes	Optrode
<p>Fixed ground-joint diaphragm</p> <ul style="list-style-type: none"> <li>• Insensitive to contamination</li> <li>• Attractive price/performance ratio</li> </ul> <p><b>Ordering Info:</b> 6.0262.100<sup>1)</sup></p> <p>iEcotrode Plus: 6.0280.300<sup>2)</sup></p> <p>dEcotrode Plus: 6.00201.300<sup>3)</sup></p>	<p>Twin-pore diaphragma</p> <ul style="list-style-type: none"> <li>• For routine measurements in similar samples</li> <li>• With lifetime indicator</li> <li>• Maintenance-free reference electrolyte (gel)</li> </ul> <p><b>Ordering Info:</b> 6.0221.100<sup>1)</sup></p> <p>6.00221.600 (with Pt-1000)<sup>4)</sup></p>	<p>Fixed ground-joint diaphragm</p> <ul style="list-style-type: none"> <li>• Very low alkali error</li> <li>• Insensitive to contamination</li> <li>• High-temperature resistance</li> </ul> <p><b>Ordering Info:</b> 6.0258.600 (with Pt1000)<sup>4)</sup></p> <p>iUnitrode: 6.0278.300<sup>2)</sup></p> <p>dUnitrode: 6.00200.300<sup>3)</sup></p> <p>Unitrode <i>easyClean</i>: 6.0260.010 (with Pt1000, 1.2 m fixed cable)</p>	<p>Fixed ground-joint diaphragm</p> <ul style="list-style-type: none"> <li>• Rapid response in iondeficient or weakly buffered solutions</li> <li>• Insensitive to contamination</li> </ul> <p><b>Ordering Info:</b> 6.0253.100<sup>1)</sup></p> <p>6.0257.600 (with Pt1000)<sup>4)</sup></p> <p>iAquatrode Plus: 6.0277.300<sup>2)</sup></p> <p>dAquatrode Plus: 6.00202.300<sup>3)</sup></p>	<p>Ground-joint diaphragm</p> <ul style="list-style-type: none"> <li>• For difficult matrices</li> <li>• Double-junction construction</li> </ul> <p><b>Ordering Info:</b> 6.0255.100<sup>1)</sup> (Fitting length = 113 mm)</p> <p>6.0255.110<sup>1)</sup> (Fitting length = 170 mm)</p> <p>6.0255.120<sup>1)</sup> (Fitting length = 310 mm)</p> <p>dProfitrode: 6.00204.300<sup>3)</sup></p>	<p><i>easyClean</i> diaphragm</p> <ul style="list-style-type: none"> <li>• For titrations in nonaqueous matrices</li> <li>• Rapid response in organic solvents</li> <li>• Electrically shielded</li> </ul> <p><b>Ordering Info:</b> 6.0229.010 (1 m fixed cable)</p> <p>6.0229.020 (2 m fixed cable)</p> <p>dSolvotrode: 6.00203.300<sup>3)</sup></p>	<p>Maintenance-free pH glass reference system</p> <ul style="list-style-type: none"> <li>• Also available as micro-electrodes with 6.4 mm diameter</li> </ul> <p><b>Ordering Info:</b> Ag-Titrode: 6.00430.100<sup>1)</sup></p> <p>Pt-Titrode: 6.0431.100<sup>1)</sup></p> <p>Micro-Au Titrode: 6.0435.110<sup>1)</sup></p> <p>iAg-Titrode: 6.00470.300<sup>2)</sup></p> <p>iPt-Titrode: 6.0471.300<sup>2)</sup></p> <p>dAg-Titrode: 6.00404.300<sup>3)</sup></p> <p>dPt-Titrode: 6.00401.300<sup>3)</sup></p>	<p>Ceramic-pin diaphragm</p> <ul style="list-style-type: none"> <li>• For precipitation or redox titrations with change of pH value</li> </ul> <p><b>Ordering Info:</b> Ag: 6.00450.100<sup>1)</sup></p> <p>Pt: 6.0451.100<sup>1)</sup></p> <p>Au: 6.0452.100<sup>1)</sup></p> <p>iAg: 6.00450.300<sup>2)</sup></p> <p>iPt: 6.0451.300<sup>2)</sup></p> <p>dAg: 6.00402.300<sup>3)</sup></p> <p>dPt: 6.00403.300<sup>3)</sup></p>	<p>For surfactant titrations in non-aqueous media</p> <p><b>Surfactrode Resistant (6.0507.130)<sup>1)</sup></b> Resistant to chloroform and many other solvents.</p> <p><b>Surfactrode Refill (6.0507.140)<sup>1)</sup></b> Refillable surfactant electrode, thus practically unlimited working life. Not resistant to chloroform.</p> <p><b>Surfactrode refill paste: 6.2319.000</b></p>	<p>For surfactant titrations in aqueous media</p> <p><b>NIO electrode (6.0507.010)<sup>1)</sup></b> For titration of non-ionic surfactants</p> <p><b>Ionic Surfactant (6.0507.120)<sup>1)</sup></b> Optimized for anionic surfactants</p> <p><b>Cationic Surfactant (6.0507.150)<sup>1)</sup></b> Optimized for cationic surfactants</p>	<p>Power supplied by titrator via USB, no separate power adapter needed</p> <p><b>Ordering Info:</b> 6.1115.000</p> <p>Optional: 6.2166.000 USB power supply unit for titrators without USB connector.</p>
										
<p>Store in 6.2323.000 storage solution only.</p> <p>Do not wipe electrode.</p> <p>For cleaning/care 6.2325.000 pHit kit is recommended.</p>	<p>Store in 6.2308.000 KCl sat.</p> <p>Rinse with water or ethanol to remove contamination.</p>	<p>Use 6.2308.040 Idrolyte for titrations at temperatures 80...100 °C or to suppress protein precipitation by KCl.</p> <p>Rinse with water or ethanol.</p> <p>Do not wipe electrode.</p> <p>Unitrodes filled with c(KCl)=3 mol/L should be stored in 6.2323.000 storage solution.</p> <p>For cleaning/care 6.2325.000 pHit kit is recommended.</p>	<p>Store in 6.2323.000 storage solution only.</p> <p>Do not wipe electrode.</p> <p>For cleaning/care 6.2325.000 pHit kit is recommended.</p> <p>Lift sleeve ring for cleaning of ground-joint diaphragm. If the sleeve ring is blocked, immerse in hot water for a few minutes to dissolve adhering material.</p> <p>Spare ground-joint diaphragm for Profitrodes 6.0255.1X0: Order no. 6.1243.020</p>	<p>Store in bridge electrolyte.</p> <p>Do not wipe electrode.</p> <p>For cleaning/care 6.2325.000 pHit kit is recommended.</p> <p>Condition glass membrane in dist. water before next titration.</p> <p>Alternative electrolyte: c(TEABr) = 0.4 mol/L in ethylene glycol (6.2320.000).</p>	<p>Store in reference electrolyte.</p> <p>Do not wipe electrode.</p> <p>In case of contamination with organic residues, immerse sensor in an appropriate solvent (for 30 min.).</p> <p>Condition glass membrane in dist. water before next titration.</p> <p>Alternative electrolyte: c(TEABr) = 0.4 mol/L in ethylene glycol (6.2320.000).</p>	<p>Store in distilled water only.</p> <p>Testing of the electrodes according to AB-048.</p> <p>Ag Titrodes also available with Ag<sub>2</sub>S or Ag-halide coating.</p>	<p>Store in reference electrolyte.</p> <p>Testing of the electrodes according to AB-048.</p> <p>Ag ring electrodes also available with Ag<sub>2</sub>S or Ag-halide coating.</p>	<p>Store dry.</p> <p>Allow a few titrations for conditioning of the Surfactrodes.</p> <p>Reactivate Surfactrode Resistant with fine-grained sandpaper if response is poor.</p> <p>Testing of the electrodes according to AB-305.</p>	<p>Store dry.</p> <p>Rinse with dist. water or 20% methanol in water.</p> <p>Wipe carefully with a methanol-soaked tissue to remove any adhering contamination.</p> <p>Do not use in organic matrices or at temperatures &gt;40 °C.</p> <p>Testing of the electrodes according to AB-305.</p>	<p>Allow 5 minutes for warming up the LED.</p>

<sup>1)</sup> Electrodes with no cable and no temperature probe have plug-in head G

<sup>2)</sup> An iConnect 2.854.0010 is required to connect an iTrode to the instrument.

<sup>3)</sup> dTrodes can only be used with a digital measuring module at an OMNIS Titrator / Titration Module.

<sup>4)</sup> Electrodes with no cable and with an integrated temperature sensor have plug-in head U