6.0913.130 Measuring vessel cover
6.0913.130 Messgefäß-Deckel
6.0913.130 Couvercle de bécher de mesure

EN

6.0913.130 Measuring vessel cover

6.0913.130 Measuring vessel cover with integrated conductometric measuring cell.

General

Check the measuring cell for transportation damage immediately upon arrival.

Storage

The measuring cell has to be stored in a dry place. It should be thoroughly rinsed with distilled H₂O before usage.

Measurement

The two metal pins of the conductivity measuring cell must entirely be immersed in the measuring solution during the whole measurement. For this, there should not be more than 20 mL of the measuring solution in the measuring vessel. The evaporation loss during the measurement is usu-

ally in the range of 5-7 mL per day and should be compensated by a correspondingly higher volume at the beginning of the measurement or by regular addition of small volumes for very long measurements (> 3 days).

Ensure that no air or nitrogen bubbles directly sweep over the electrode during the measurement. Should the situation arise, adjust the PTFE tubing for the gas inlet to the measuring vessel differently.

Cleaning

For easier cleaning of the two metal pins of the conductivity measuring cell, the protection ring can be removed. It should however be remounted after the cleaning because it is important for the geometry of the measuring cell.

Clean the measuring cell first with an organic solvent (e.g. acetone, acetone:ethanol mixture (1:1) or 2-propanol). After this, thoroughly rinse with distilled H₂O.

With strong and oily contaminations the measuring cell can first be cleaned mechanically with the aid of a brush, dishwashing agent and hot water and then with organic solvent as well as with distilled H₂O.

Calibration

For the determination of the induction time a calibration of the cell constant is not necessary because the absolutely measured conductivity does not matter in this case - only the relative changing is being evaluated.

For the determination of the stability time or if the GLP test "Conductivity" is to be carried out, the cell constant must however be determined. Use a standard solution with a known conductivity, as e.g. the 6.2301.060 Metrohm conductivity standard (must be diluted 1:100 with distilled H₂O before usage) or 6.2324.010 (can be used directly). The exact instructions for the procedure of the calibration can be found in the manual of your stability measurement device.

Specifications

The conductivity measuring cell fulfills the following specifications at the time of delivery:

Cell constant	1.1 cm ⁻¹
Tolerance of the cell constants	± 0.055 cm ⁻¹

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