Platinum Electrode Tips for CVS Platin Elektrodentips für CVS

6.1204.170 / 6.1204.190 / 6.1204.610

6.1204.170 Platinum electrode tip 3 mm for CVS
6.1204.190 Platinum electrode tip 1 mm for CVS
6.1204.610 Platinum electrode tip 2 mm for CVS

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EN

Platinum Electrode Tips for CVS

The Platinum electrode tip (6.1204.170 or 6.1204.190 or 6.1204.610) together with a driving axle (6.1204.210 or 6.1204.220) constitute the Rotating Disk Electrode (RDE) for VA Stands.

Caution

- Remove the protective cap before using the electrode.
- Protect the electrode against mechanical load!
- PEEK electrode tips may not be exposed to concentrated mineral acids as they will otherwise be irreversibly damaged!
- Avoid extensive use of the electrode tips in organic solvents.
- PEEK electrode tips may be used at a maximum temperature of 40 °C.
- Glass electrode tips may be used at a maximum temperature of 50 °C.
- Do not use the electrode outside the recommended temperature range.
- It is recommended to use an individual electrode kit, consisting of working electrode (WE), reference electrode (RE), auxiliary electrode (AE), and measuring vessel for each application.
- The electrodes are polished extensively ex works, we do not recommend to repolish them mechanically. Repolishing the electrode will impair the signal-to-noise ratio and damage the electrode irreversibly.
- Do not treat the electrode tips in ultrasonic baths!

Preparation

Rinse the electrode with ethanol and distilled water before use.

Regeneration

The RDE is a solid-state electrode with a stationary surface. During use, this surface becomes increasingly contaminated with the products of the electrode redox processes.

Proceed as follows to clean the electrode:

- Soak the electrode tip with its Pt surface for about 10 minutes in a sodium hydroxide solution $c(\text{NaOH}) = 0.5 \text{ mol/L}$ (at least p.A. quality).
- Thoroughly rinse the electrode tip with distilled water.

For CVS applications: as an alternative or in addition to cleaning with sodium hydroxide, the electrode can also be cleaned in pure Virgin Make-up Solution (VMS) applying cyclic voltage sweeps.

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