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
#### 6.0223.100 Solitrode HF

##### General

Immediately after receiving the electrode, check to make sure that it works properly. Electrodes that do not work properly must be sent back for warranty processing within two months (starting from the day of delivery). If the defect is proven to be due to a material or manufacturing defect, the electrode will be replaced at no charge. The transport costs are to the customer's account.

##### Caution

Electrodes must not be treated in the ultrasonic bath.

##### General preparations

- Open the closure of the refill aperture (1).
- Add the reference electrolyte c(KCl) = 3 mol/L up to the refill aperture.

Ensure that the refill aperture remains open during use.

Always use plastic vessels for samples containing hydrofluoric acid.

##### Polishing the Sb pin (Sb electrode)

You can polish the Sb pin (3) either using a polishing agent, e. g.  $\text{Al}_2\text{O}_3$  (6.2802.000) or toothpaste. Proceed as follows:

- Apply the polishing agent to a cloth.
- Moisten the polishing agent with water.

Make sure not to contaminate the diaphragm (2) with polishing agent.

##### Scope of application of the Sb electrode

We recommend using the Sb electrode primarily for titration. If using the Sb electrode for direct measurements, the offset voltage of the pH meter must be set to  $-366$  mV. No borate-based buffer solution may be used for calibrating the Sb

electrode. In addition, make sure that the stirrer is switched off for the measured value acceptance.

##### Caution


##### Cross-sensitivity of the Sb electrode

The Sb electrode is not suitable for titrations in solutions that contain complexing agents for  $\text{Sb}^{3+}$  ions. This applies, for example, to tartrates, citrates, oxalates and borates. Sulfides act as an electrode poison at concentrations as low as 100 mg/L.

##### Scope of application of the Solitrode HF

We recommend the Solitrode HF for titration. Its advantage compared to the Sb electrode lies in lower cross-sensitivity to other ions and better reproducibility.

The recommended fields of application of the electrodes are illustrated in the following graph.

- A Sb electrode
- B Solitrode HF
- C pH glass electrode

##### Storing electrodes

We recommend storing electrodes with c(KCl) = 3 mol/L as reference electrolyte in storage solution (6.2323.000).

### Troubleshooting

#### Sb electrode

- The measured value setting is sluggish.

You can polish the Sb pin (3) using either a polishing agent, e. g.  $\text{Al}_2\text{O}_3$  (6.2802.000) or toothpaste.

#### Sb electrode and Solitrode HF

- Solution containing protein.

Regularly soak the electrode for one hour in a 1% pepsin solution in c(HCl) = 0.1 mol/L. Then rinse the electrode well with distilled water.

- Sulfide-based solution: In the diaphragm, there is a black precipitate of  $\text{Ag}_2\text{S}$ .

Soak the electrode for one hour in a freshly prepared 7% solution of thiourea in c(HCl) = 0.1 mol/L. Then rinse the electrode well with distilled water.

- There are air bubbles in the electrode.

Eliminate the air bubbles in the electrode by applying a downwards centrifugal motion to it.

- The reference electrolyte is dirty or has dried out.

Aspirate and replace the electrolyte. If necessary, repeat this process several times.

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