
Rack base

Product information

Produktinformation

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EN

DE

CN

Metrohm AG
Ionenstrasse
CH-9100 Herisau
Switzerland
+41 71 353 85 85
info@metrohm.com
www.metrohm.com

EN

Rack base

6.02001.010

Product information

8.0108.8007EN /
2021-07-23

This documentation is protected by copyright. All rights reserved.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Disclaimer

Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc., are expressly excluded from the warranty. Unauthorized modifications to the product (e.g. conversions or attachments) exclude any liability on the part of the manufacturer for resulting damage and its consequences. Instructions and notes in the Metrohm product documentation must be strictly followed. Otherwise, Metrohm's liability is excluded.

Table of contents

1	Overview	1
1.1	Rack base – Product description	1
1.2	Rack base – Overview	1
1.3	OMNIS sample rack – Overview	2
1.4	Symbols and conventions	3
1.5	Further information	3
1.6	Accessories	4
2	Safety	5
2.1	Intended use	5
2.2	Responsibility of the operator	5
2.3	Requirements for operating personnel	6
2.4	Safety instructions	6
2.4.1	Danger from electrical potential	6
2.4.2	Danger from biological and chemical hazardous substances	6
2.4.3	Danger from highly flammable substances	7
2.4.4	Danger from leaking liquids	7
2.4.5	Danger during transport of the product	8
2.4.6	Dangers from hot surfaces and liquids	8
2.5	Design of warning messages	8
2.6	Meaning of warning signs	9
3	Technical specifications	11
3.1	Ambient conditions	11
3.2	Rack base / module base – Dimensions	11
3.3	Module base – Housing	11
3.4	Rack base – Housing	12
3.5	Rack base – Connectors specifications	12
3.6	OMNIS sample rack – Connectors specifications	12
3.7	OMNIS sample rack – Sample processing specifications	12

1 Overview

1.1 Rack base – Product description

The rack base is the skeletal structure for the sample robot. The OMNIS sample racks are placed on the rack holders of the rack base. Sensors in the rack holder let the system know if and what sample rack is present.

NOTICE

As a basic rule, the installation of the rack base is carried out by the regional Metrohm service representative.

1.2 Rack base – Overview


Figure 1 Rack base – Overview

1 Rack base

2 Rack holder

On each rack base (1-**1**), up to 2 OMNIS sample racks can be placed on the rack holders (1-**2**).

1.3 OMNIS sample rack – Overview


Figure 2 OMNIS sample rack – Overview

1	Sample vessels	2	Sample position
3	Transport handles	4	Spout hole
5	Label holder		

The sample vessels (2-1) are placed in sample positions (2-2) in the OMNIS sample rack.

The transport handles (2-3) allow for the OMNIS sample rack to be transported by hand so that it can be set down on or removed from the rack holder of the rack base. Several empty OMNIS sample racks can be stacked on top of one another on the transport handles.

The OMNIS sample rack contains spout holes (2-4) at the top for emptying liquids that seeped through due to overflowing, spilling or cleaning.

NOTICE

The OMNIS sample rack is not dishwasher-safe.

Table 1 Versions of the OMNIS sample rack

Vessel volume	Number of vessels	Article number
250 mL	9	6.02041.010

Vessel volume	Number of vessels	Article number
200 mL	9	6.02041.020
150 mL	9	6.02041.050
120 mL	16	6.02041.030
75 mL	25	6.02041.040

1.4 Symbols and conventions

The following formatting may appear in the documentation:

(5- 12)	Cross-reference to figure legend The first number refers to the figure number. The second number refers to the product part in the figure.
1	Instruction step Numbers indicate the order of the instructions steps.
Method	Names of parameters, menu items, tabs and dialog windows
File ► New	Menu path
[Continue]	Button or key

1.5 Further information

The Metrohm Knowledge Base <https://guide.metrohm.com> always provides the current version of this document. Further instructions, leaflets, release notes etc. may be available, depending on the product. You can directly access the required information or the associated PDF document using the full-text search function and filters.

Downloading the accessories list

- 

NOTICE

4 ■ ■ ■ ■ ■ ■ ■ ■

2 Safety

2.1 Intended use

Metrohm products are used for the analysis and handling of chemicals.

Usage therefore requires the user to have basic knowledge and experience in handling chemicals. Knowledge regarding the application of fire prevention measures prescribed for laboratories is also mandatory.

Adherence to this technical documentation and compliance with the maintenance specifications make up an important part of intended use.

Any utilization in excess of or deviating from the intended use is regarded as misuse.

Specifications regarding the operating values and limit values of individual products are contained in the "Technical specifications" section, if relevant.

Exceeding and/or not observing the mentioned limit values during operation puts people and components at risk. The manufacturer assumes no liability for damage due to non-observance of these limit values.

The EU declaration of conformity loses its validity as soon as modifications are carried out on the products and/or the components.

2.2 Responsibility of the operator

The operator must ensure that basic regulations on occupational safety and accident prevention in chemical laboratories are observed. The operator has the following responsibilities:

- Instruct personnel in the safe handling of the product.
- Train personnel in the use of the product according to the user documentation (e.g. install, operate, clean, eliminate faults).
- Train staff on basic occupational safety and accident prevention regulations.
- Provide personal protective equipment (e.g. protective glasses, gloves).
- Provide suitable tools and equipment to carry out the work safely.

The product may be used only when it is in perfect condition. The following measures are required to ensure the safe operation of the product:

- Check the condition of the product before use.
- Remedy defects and malfunctions immediately.
- Maintain and clean the product regularly.

- Basic regulations on occupational safety and accident prevention for chemical laboratories are known and complied with.
- Knowledge of handling hazardous chemicals is present. Personnel have the ability to recognize and avoid potential dangers.
- Knowledge regarding the application of fire prevention measures for laboratories is available.
- Safety-relevant information is communicated and understood. The personnel can operate the product safely.
- The user documentation has been read and understood. The personnel operate the product according to the instructions in the user documentation.

2.4 Safety instructions

2.4.1 Danger from electrical potential

- Operate the product only if it is in perfect condition. The housing must also be intact.
- Only use the product with the covers fitted. If covers are damaged or missing, disconnect the product from the energy supply and contact the regional Metrohm service representative.
- Protect live components (e.g. power supply unit, power cord, connection sockets) against moisture.
- Always have maintenance work and repairs on electrical components carried out by a regional Metrohm service representative.
- Disconnect the product from the energy supply immediately if at least one of the following cases occurs:
 - The housing is damaged or open.
 - Live parts are damaged.
 - Moisture penetrates.

2.4.2 Danger from biological and chemical hazardous substances

Contact with biological hazardous substances may cause poisoning from toxins or infections from microorganisms. Contact with aggressive chemical substances may cause poisoning or chemical burns. To avoid danger from biological or chemical hazardous substances, observe the following:

- Label the product according to regulations if it is used for substances that have a potential for chemical hazards and are generally subject to the Hazardous Substances Ordinance.
- Wear personal protective equipment (e.g. protective glasses, gloves).
- Use exhaust equipment when working with vaporizing hazardous substances.
- Dispose of hazardous substances in accordance with regulations.
- Clean and disinfect contaminated surfaces.
- Only use detergents that do not cause any unwanted side reactions with the materials to be cleaned.
- Dispose of chemically contaminated materials (e.g. cleaning material) in accordance with regulations.
- Proceed as follows in case of a return shipment to Metrohm AG or a regional Metrohm representative:
 - Decontaminate the product or product component.
 - Remove the labeling for hazardous substances.
 - Create a declaration of decontamination and enclose it with the product.

2.4.3 Danger from highly flammable substances

Using highly flammable substances or gases may cause fires or explosions. To avoid danger from highly flammable substances, observe the following:

- Avoid ignition sources.
- Use protective grounding.
- Use exhaust equipment.

2.4.4 Danger from leaking liquids

Leaking liquids may cause injuries and may damage the product. To avoid danger from leaking liquids, observe the following:

- Check the product and its accessories for leakages and loose connections.
- Replace leaking parts and connecting elements without delay.
- Tighten loose connecting elements.
- Do not loosen tubing connections under pressure.
- Do not remove aspiration tubing under pressure.
- Carefully pull the ends of the tubing out of the containers.
- Carefully let liquids from tubing drain into suitable containers.
- Insert the buret tips completely into the containers.
- Remove and dispose of leaked liquids in accordance with regulations.
- If you suspect that liquid has penetrated the instrument, disconnect the instrument from the energy supply. Then have the instrument checked by a regional Metrohm service representative.

2.4.5 Danger during transport of the product

Chemical or biological substances may be spilled during the transport of the product. Parts of the product may fall down or may be damaged. There is a risk of injury from chemical or biological substances and pieces of broken glass. To ensure safe transport, observe the following:

- Remove loose parts (e.g. sample racks, sample vessels, bottles) before transport.
- Remove liquids.
- Lift and transport the product with both hands on the base plate.
- Lift and transport heavy products only according to instructions.

2.4.6 Dangers from hot surfaces and liquids

Contact with hot surfaces or hot liquids may cause burns. To avoid the risk of injury, observe the following:

- Install and use the protective devices enclosed with the product.
- Allow hot surfaces to cool down before working on the product.
- Wear heat-resistant protective gloves.
- Clean up spilled liquids and solids immediately.

2.5 Design of warning messages

There are 4 hazard levels for warning messages. The following signal words are used for classifying the hazard levels in warning messages:

- **DANGER** indicates a hazardous situation which, if not avoided, will result in serious injury or death.
- **WARNING** indicates a hazardous situation which, if not avoided, could result in serious injury or death.
- **CAUTION** indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.
- **NOTICE** indicates a hazardous situation which, if not avoided, could result in property damage.

Warning messages differ in design (color and warning sign) depending on the hazard level:

DANGER

Type and source of danger

Consequences when not observing the notice: An irreversible injury that may result in death is very probable.

- Measures to avoid the danger

WARNING

Type or source of danger

Consequences when not observing the notice: A serious injury that may result in death is probable.

- Measures to avoid the danger

CAUTION

Type or source of danger

Consequences when not observing the notice: A minor to moderate injury is probable.

- Measures to avoid the danger

2.6 Meaning of warning signs

This documentation uses the following warning signs:

Table 2 Warning sign according to ISO 7010

Warning sign	Meaning
	General warning sign
	Warning of electrical voltage
	Warning of hand injuries
	Warning of sharp object
	Warning of hot surface
	Warning of biological hazard
	Warning of toxic materials

10 ■ ■ ■ ■ ■ ■ ■ ■

3 Technical specifications

3.1 Ambient conditions

Nominal function range	+5 to +45 °C	at max. 80% relative humidity, non- condensing
Storage	+5 to +45 °C	

3.2 Rack base / module base – Dimensions

Measurements

Width	280 mm
Height	115 mm
Depth	564 mm

Weight

Rack base	3.9 kg
Module base M/L	4.6 kg

3.3 Module base – Housing

Materials

Cover	AW-5754 H12 / H22	aluminum, coated
Base	AlSi12Cu	aluminum, coated

IP degree of protection	IP 40
-------------------------	-------

3.4 Rack base – Housing

Materials

<i>Cover</i>	AW-5754 H12 / H22	aluminum, coated
<i>Base</i>	AlSi12Cu	aluminum, coated

IP degree of protection	IP 20
-------------------------	-------

3.5 Rack base – Connectors specifications

Contacts	4	Contact surfaces for rack base
----------	---	--------------------------------

3.6 OMNIS sample rack – Connectors specifications

Contacts	4	Spring contacts
----------	---	-----------------

3.7 OMNIS sample rack – Sample processing specifications

Versions

<i>Number of beaker positions</i>	9, 16, 25	from Metrohm accessories
-----------------------------------	-----------	--------------------------

Metrohm AG
Ionenstrasse
CH-9100 Herisau
Schweiz
+41 71 353 85 85
info@metrohm.com
www.metrohm.com

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