


# Liquid Handling Station



## Manual

8.108.8011EN / v4 / 2024-05-30





Metrohm AG  
CH-9100 Herisau  
Switzerland  
+41 71 353 85 85  
info@metrohm.com  
www.metrohm.com

# **Liquid Handling Station**

## **Manual**

This documentation is protected by copyright. All rights reserved.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

### **Disclaimer**

Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc., are expressly excluded from the warranty. Unauthorized modifications to the product (e.g. conversions or attachments) exclude any liability on the part of the manufacturer for resulting damage and its consequences. Instructions and notes in the Metrohm product documentation must be strictly followed. Otherwise, Metrohm's liability is excluded.

# Table of contents

<b>1</b>	<b>Introduction</b>	<b>1</b>
1.1	Description .....	1
1.2	About the documentation .....	1
1.2.1	Symbols and conventions .....	2
1.3	Accessories and additional information .....	3
1.4	Safety instructions .....	3
1.4.1	General notes on safety .....	3
1.4.2	Electrical safety .....	3
1.4.3	Flammable solvents and chemicals .....	4
1.4.4	Recycling and disposal .....	5
<b>2</b>	<b>Overview of the instrument</b>	<b>6</b>
<b>3</b>	<b>Installation</b>	<b>8</b>
3.1	Prepare the Liquid Handling Station .....	8
3.2	Mounting the Liquid Handling Station .....	16
3.3	Defining work positions .....	20
<b>4</b>	<b>Operation and maintenance</b>	<b>24</b>
4.1	General notes .....	24
4.1.1	Care .....	24
4.1.2	Maintenance by Metrohm .....	24
<b>5</b>	<b>Technical specifications</b>	<b>25</b>
5.1	Magnetic stirrer .....	25
5.2	Energy supply .....	25
5.3	Interfaces and connectors .....	25
5.4	Ambient temperature .....	25
5.5	Reference conditions .....	25
5.6	Dimensions/material .....	26
5.7	Volume .....	26
	<b>Index</b>	<b>27</b>







## 1.3 Accessories and additional information

Additional information is available on the Metrohm website (<https://www.metrohm.com>):

- Product family
- Product versions
- Accessories
- Documents about the product

### Downloading the accessories list



#### NOTE

The accessories list is a part of the product documentation. Download the accessories list and store it as a reference.

1. Use the search function to search for the product.
2. Open the desired product version.
3. Download the accessories list.

## 1.4 Safety instructions

### 1.4.1 General notes on safety



#### WARNING

Operate this instrument only according to the information contained in this documentation.

This instrument left the factory in a flawless state in terms of technical safety. To maintain this state and ensure non-hazardous operation of the instrument, the following instructions must be observed carefully.

### 1.4.2 Electrical safety

The electrical safety when working with the instrument is ensured as part of the international standard IEC 61010.



#### WARNING

Only personnel qualified by Metrohm are authorized to carry out service work on electronic components.



## WARNING

Never open the housing of the instrument. The instrument could be damaged by this. There is also a risk of serious injury if live components are touched.

There are no parts inside the housing which can be serviced or replaced by the user.

## Supply voltage



## WARNING

An incorrect supply voltage can damage the instrument.

Only operate this instrument with a supply voltage specified for it (see rear panel of the instrument).

### Protection against electrostatic charges



## WARNING

Electronic components are sensitive to electrostatic charges and can be destroyed by discharges.

Do not fail to pull the power cord out of the power socket before you set up or disconnect electrical plug connections at the rear of the instrument.

### 1.4.3 Flammable solvents and chemicals




## WARNING

All relevant safety measures are to be observed when working with flammable solvents and chemicals.

- Set up the instrument in a well-ventilated location (e.g. fume cupboard).
- Keep all sources of flame far from the workplace.
- Clean up spilled liquids and solids immediately.
- Follow the safety instructions of the chemical manufacturer.

#### 1.4.4 Recycling and disposal



This product is covered by European Directive 2012/19/EU, WEEE – Waste Electrical and Electronic Equipment.

The correct disposal of your old instrument will help to prevent negative effects on the environment and public health.

More details about the disposal of your old instrument can be obtained from your local authorities, from waste disposal companies or from your local dealer.

## 2 Overview of the instrument

Liquid Handling Station in the left-handed version, complete with attachment holder for assembly on a sample changer.




Figure 1 Overview of the instrument Liquid Handling Station (left-handed version)

<b>1</b>	<b>Cover for the mixing vessel</b>	<b>2</b>	<b>Overflow</b> with connector
<b>3</b>	<b>Mixing vessel</b>	<b>4</b>	<b>Mixing vessel connector - UNF 10/32</b> sealed with threaded stopper
<b>5</b>	<b>Rinsing unit</b>	<b>6</b>	<b>Main body of the Liquid Handling Station</b> with magnetic stirrer dummy
<b>7</b>	<b>Waste connector</b>	<b>8</b>	<b>Rinsing connector - UNF 10/32</b> sealed with threaded stopper

<b>9</b>	<b>Cable clip</b>	<b>10</b>	<b>Support bracket</b>
<b>11</b>	<b>Base plate</b>	<b>12</b>	<b>Clamping fastener</b>

### 3 Installation

The installation of the Liquid Handling Station is accomplished in 3 main steps.

1. Prepare the Liquid Handling Station.
2. Mounting and aligning the Liquid Handling Station on the sample changer.
3. Defining the work positions of the sample needle for the Liquid Handling Station.

### 3.1 Prepare the Liquid Handling Station

If the Liquid Handling Station is to be used for dilutions, then it must be converted.

The following instruments are required for this purpose:

- Liquid Handling Station
  - Left-handed version, 6.2841.120  
or
  - Right-handed version, 6.2841.130
- Magnetic stirrer, 2.741.0010



## CAUTION

Use caution and do not apply excessive force when carrying out the following procedures.

The conversion is carried out through the following steps:

### Removing the magnetic stirrer dummy

- 1** Remove the mixing vessel from the main body by hand.



- 
- This exploded view diagram illustrates the assembly of the microscope. The components are arranged vertically from top to bottom: the eyepiece (yellow), the objective (yellow), the stage (purple), and the base (purple). The eyepiece is shown with a white lens and a green cap. The objective is shown with a green cap. The stage is shown with a white lens and a green cap. The base is shown with a white lens and a green cap. The diagram shows how these components fit together to form the complete microscope.

Figure 3 Removing the holder

- 3** Loosen the 2 screws on the underside of the main body and pull out the magnetic stirrer dummy.



## Installing the magnetic stirrer

- 1 Insert the magnetic stirrer from below into the main body and place the connection cable in the corresponding recess.




Figure 6 Inserting the magnetic stirrer

- 2** Insert the 2 screws and tighten them.



Figure 7 Tightening the magnetic stirrer

■■■■■■■■■■

- 3** Place the previously removed holder against the main body and connect the two parts with the three screws.




Figure 8 Screwing on the holder

- 4** Mount the mixing vessel and press it all the way down.



## NOTE

Observe the recess on the bottom rim of the mixing vessel for positioning.

It must snap into the main body.




Figure 9 Attaching the mixing vessel

- 5 Fixing the connection cable in place with the cable clip.



**NOTE**

Ensure that the connection cable is not pinched and damaged.




Figure 10 Fastening the connection cable in place

[illegible]

- 6** Connect the connection cable to the stirrer connector of the sample changer, see *Manual 858 Professional Sample Processor* (8.858.8002).




Figure 11 858 Professional Sample Processor – Stirrer connector

### Mounting the rinsing station drainage tubing

### Required accessories

- PVC tubing / 4 mm / 6 mm / 2 m (6.1801.120)

## 1 Mounting the drainage tubing

- Push one end of the PVC tubing over the waste connector (1-7) of the rinsing unit.
- Guide the other end of the PVC tubing to the waste container.

### Mounting the dilution station drainage tubing

### Required accessories

- Silicone tubing 8 mm internal diameter / 1 m (6.1816.070)
- Y tubing connector, PP, 3 x 7-8 mm (6.01808.010)




Figure 12 Drainage tubing

**1 Ion chromatograph drainage tubing**  
Bottle holder connection/detector chamber  
– Y connector

**2 Silicone tubing 6.1816.070**  
long tubing piece  
Dilution station overflow connection – Y  
connector

**3 Y tubing connector, PP, 3 x 7-8 mm**  
**6.01808.010**

**4 Silicone tubing 6.1816.070**  
short tubing piece  
Y connector connection – base tray drainage  
tubing connection

### 1 Converting the drainage tubing on the ion chromatograph

- Loosen the drainage tubing on the ion chromatograph (12-**1**) from the drainage tubing connection on the base tray.
- Attach the drainage tubing to one end of the Y connector (12-**3**).

## 2 Mounting the silicone tubing 6.1816.070

- Cut off an approx. 10 cm-long piece of the silicone tubing. Connect the Y connector and the drainage connection of the base tray with the short piece of tubing (12-4).
- Cut off an approx. 40 to 80 cm-long piece of the silicone tubing. The length of the tubing depends on the placement of the sample changer and the ion chromatograph. Cut the piece of tubing as short as possible. Connect the overflow of the dilution station and the Y connector with the long piece of tubing (12-2).



## NOTE

If necessary, shorten the long piece of tubing until the liquid always flows downwards. A U-shaped bend in the tubing or a horizontal position of the tubing could cause liquid to back up.

## 3.2 Mounting the Liquid Handling Station

The Liquid Handling Station requires precise mounting and alignment in order to function properly with the entire system.



## CAUTION

## Faulty, imprecise assembly

The components could become damaged during operation if assembled in a faulty and imprecise manner.


- Observe the assembly sequence.
- Observe instrument movements during setup and stop in a timely fashion if necessary.

The following description refers to the assembly of a left-handed Liquid Handling Station. The right-handed version is assembled in the same way but mirror-reversed.




- #### 4 Assembly rail





*Figure 15 Positioning the Liquid Handling Station*  
The precise position is shown in the following figure.



*Figure 16 Alignment with the retaining plate*

**1 Sample needle swing range**

- 3** Retighten the screw (13-2) on the clamping fastener.
- 4** Check whether the distance to the sample rack is consistent with *Figure 14* and whether the position in relation to the retaining plate (see *figure 16, page 19*) is suitable.

### 3.3 Defining work positions

The Liquid Handling Station has 3 positions. These work positions have to be defined in the software (e.g. **MagIC Net™**).

## Defining the work positions for the sample needle

The work positions are defined in the software as **external positions**.

You can find the value settings under the following path:

**Configuration ▶ Devices ▶ Properties ▶ Tower ▶ Swing Head ▶ External position**



## NOTE

### Value settings of the bracket **External position**

These value settings are guidelines that may need to be corrected when the Liquid Handling Station is set up.




Figure 17 External position – Guide values

<b>1</b>	<b>External position 1</b>	<b>2</b>	<b>External position 2</b>
<b>3</b>	<b>External position 3</b>	<b>4</b>	<b>External position 4</b>





Figure 18 External positions

<b>1</b>	<b>Rinsing unit - Rinse</b>	<b>2</b>	<b>Rinsing unit - Disposal</b>
<b>3</b>	<b>Dilution unit</b>	<b>4</b>	<b>Dilution unit</b>

**1** Move to **External position 1** (18-1).





- 10** If necessary, correct the angle value.





<i>Width</i>	131 mm
<i>Height</i>	94 mm
<i>Depth</i>	224 mm
<i>Weight</i>	660 g (without accessories)
<i>Material</i>	
<i>Mixing vessel/     rinsing unit</i>	PP
<i>Holder</i>	Aluminum sheet, stove-enameled

Mixing vessel, total	15 mL
Mixing vessel dead volume	5 mL
Rinsing unit - Rinse	3.1 mL
Rinsing unit - Disposal	5.1 mL

