Differential Amplifier

Manual
8.110.1998

Differenzverstärker Differential Amplifier Amplificateur différentiel Amplificador diferencial

**Gebrauchsanweisung
Instructions for Use
Mode d'emploi
Instrucciones para el uso**

Gebrauchsanweisung – Deutsch	1
Instructions for Use – English	9
Mode d'emploi – Français	17
Instrucciones para el uso – Español.....	25

Einleitung

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
Declaration of Conformity

This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.

CH-9101 Herisau, Switzerland
Tel. +41 71 353 85 85
Fax +41 71 353 89 01
www.metrohm.com

<i>Name of commodity</i>	Differenzverstärker
<i>System software</i>	Stored in ROMs
<i>Name of manufacturer</i>	Metrohm Ltd., Herisau, Switzerland
<i>Description</i> Device for titrators without built-in differential amplifier to make them usable for differential potentiometry.	
This Metrohm instrument has been built and has undergone final type testing according to the standards: <i>Electromagnetic compatibility: Emission</i> EN/IEC 61326, EN 55022 / CISPR 22 <i>Electromagnetic compatibility: Immunity</i> EN/IEC 61326, EN/IEC 61000-4-2, EN/IEC 61000-4-3, EN/IEC 61000-4-4, EN/IEC 61000-4-5, EN/IEC 61000-4-6, EN/IEC 61000-4-8, EN/IEC 61000-4-11, EN/IEC 61000-4-14 <i>Safety specifications</i> EN/IEC 61010-1, UL3101-1 <i>The instrument meets the requirements of the CE mark as contained in the EU directives 89/336/EEC and 73/23/EEC and fulfils the following specifications:</i>	
EN 61326	Electrical equipment for measurement, control and laboratory use – EMC requirements
EN 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use
Metrohm Ltd. Is holder of the SQS-certificate of the quality system ISO 9001 for quality assurance in design/development, production, installation and servicing. The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality and performance. The technical specifications are documented in the instruction manual.	
Herisau, December 20, 2002  Dr. J. Frank Vice President Head of R&D  Ch. Buchmann Vice President Head of Production Responsible for Quality Assurance	

Quality Management Principles

Metrohm Ltd., CH-9101 Herisau, Switzerland

 **Metrohm**
l o n a n a l y s i s
CH-9101 Herisau/Switzerland
E-Mail info@metrohm.com
Internet www.metrohm.com

Metrohm Ltd. Holds the ISO 9001 Certificate, registration number 10872-02, issued by SQS (Swiss Association for Quality and Management Systems). Internal and external audits are carried out periodically to assure that the standards defined by Metrohm's QM Manual are maintained.

The steps involved in the design, manufacture and servicing of instruments are fully documented and the resulting reports are archived for ten years. The development of software for PCs and instruments is also duly documented and the documents and source codes are archived. Both remain the possession of Metrohm. A non-disclosure agreement may be asked to be provided by those requiring access to them.

The implementation of the ISO 9001 quality system is described in Metrohm's QM Manual, which comprises detailed instructions on the following fields of activity:

Instrument development

The organisation of the instrument design, its planning and the intermediate controls are fully documented and traceable. Laboratory testing accompanies all phases of instrument development.

Software development

Software development occurs in terms of the software life cycle. Tests are performed to detect programming errors and to assess the program's functionality in a laboratory environment.

Components

All components used in the Metrohm instruments have to satisfy the quality standards that are defined and implemented for our products. Suppliers of components are audited by Metrohm as the need arises.

Manufacture

The measures put into practice in the production of our instruments guarantee a constant quality standard. Production planning and manufacturing procedures, maintenance of production means and testing of components, intermediate and finished products are prescribed.

Customer support and service

Customer support involves all phases of instrument acquisition and use by the customer, i.e. consulting to define the adequate equipment for the analytical problem at hand, delivery of the equipment, user manuals, training, after-sales service and processing of customer complaints. The Metrohm service organisation is equipped to support customers in implementing standards such as GLP, GMP, ISO 900X, in performing Operational Qualification and Performance Verification of the system components or in carrying out the System Validation for the quantitative determination of a substance in a given matrix.

Introduction

These Instructions for Use give you a overview of the applications of the Differential Amplifier.

Parts and controls


Fig. 2: Front view of the Differential Amplifier 6.2129.010

1 Measuring interface for reference electrode

2 Measuring interface for indicator electrode

3 Connecting cable
Connection to the measuring interface of the titrator

4 Measuring interface for auxiliary electrode

5 Connection for power supply unit

Functional principle

Potentiometric measurements in poorly conducting media, e.g. organic solvents, are seriously affected by static charges. This can go that far that evaluations of titration curves are complicated or even impossible. Especially strong field strengths are generated by friction on insulators such as plastic floors, clothing made of synthetic fibers etc.; conditions that occur in any laboratory. The interfering voltages are superimposed on the measuring signal and this can lead to imaginary endpoints in titrations which renders automatic evaluation almost impossible (see Application Notes T-41 and T-44).

Problems of this kind can be solved by differential potentiometry. In these measurements, indicator and reference electrode are each connected to a high impedance input. It is very important that both electrodes are shielded identically as far as possible and are therefore symmetrical with respect to their taking up of interfering voltages. An auxiliary electrode establishes the galvanic connection between the reference point of the amplifier circuit and the measuring solution.

Practical hints

- Glass electrodes should be conditioned for about 1 hour in the solvent used.
- For titrations: If the potential jump after the first volume increment is too high a small start volume may put things right.

Electrodes

Reference electrode

<i>Order-No.</i>	<i>Description</i>
6.0729.100	Ag/AgCl reference electrode, double junction, screened, with Metrohm plug-in head G; for measuring vessels
6.0729.110	Ag/AgCl reference electrode, double junction, screened, with Metrohm plug-in head G; for Sample Processors

Indicator electrode

<i>Order-No.</i>	<i>Description</i>
6.0431.100	Pt-Titrode with Metrohm plug-in head G; separate pH glass electrode and auxiliary electrode combined in one electrode
6.0150.100	Separate pH glass electrode with Metrohm plug-in head G

Auxiliary electrode

<i>Order-No.</i>	<i>Description</i>
6.0331.000	Separate Pt-rod electrode with Metrohm socket B
6.0301.100	Separate Pt-wire electrode with Metrohm socket B

Appendix

Technical data

Measuring amplifier

<i>Input voltage range</i>	- 4000 ... + 4000 mV
<i>Input impedance</i>	$> 10^{12} \Omega$
<i>Offset current</i>	$< 10^{-12} \text{ A}$
<i>Common mode rejection AC</i>	$> 30 \text{ dB (50 Hz)}$
<i>Common mode rejection DC</i>	60 dB (at $U_{\text{Ind/Ref}} \pm 2000 \text{ mV}$)

Power supply

<i>Power supply</i>	9 V AC
<i>Power consumption</i>	1 W

Ambient temperature

<i>Nominal working range</i>	+ 5 °C ... + 45 °C (bei max. 85 % Luftfeuchte)
<i>Storage</i>	- 20 °C ... + 60 °C
<i>Transport</i>	- 40 °C ... + 60 °C

Dimensions

<i>Housing material</i>	aluminium
<i>Width</i>	156 mm
<i>Height</i>	32 mm
<i>Depth</i>	54 mm
<i>Weight</i>	158 g

Standard equipment

Differential amplifier EU

No.	Order-No.	Description	
1	6.2117.000	Power supply unit EURO 230 V / 9 V AC	
1	6.2129.010	Differential amplifier	

Differential amplifier US


No.	Order-No.	Description	
1	6.2119.000	Power supply unit USA 115 V / 9 V AC	
1	6.2129.010	Differential amplifier	

Optional accessories

Electrodes for Sample Processors

No.	Order-No.	Description
1	6.1241.040	Shaft made of PP for electrode rods with socket B
1	6.1248.000	Pt-electrode rod for 6.1241.040
1	6.1248.040	Glassy-Carbon-electrode rod for 6.1241.040

Accessories for electrodes

No.	Order-No.	Description	
1	6.1236.020	SGJ sleeve for electrode 6.0301.100	

Cables for electrodes

<i>No.</i>	<i>Order-No.</i>	<i>Description</i>
1	6.2104.020	Connecting cable for Metrohm electrodes with plug-in head G, length 1 m
1	6.2104.030	Connecting cable for Metrohm electrodes with plug-in head G, length 2 m
1	6.2106.020	Connecting cable for Metrohm electrodes with plug B, length 1 m
1	6.2106.060	Connecting cable for Metrohm electrodes with plug B, length 2 m

Warranty and conformity

Warranty

The warranty on our products is limited to defects that are traceable to material, construction or manufacturing error which occur within 12 months from the day of delivery. In this case the defects will be rectified in our workshops free of charge. Transport costs are to be paid by the customer.

For day and night operation the warranty is limited to 6 months.

Glass breakage in the case of electrodes or other parts is not covered by the warranty. Checks which are not a result of material or manufacturing faults are also charged during the warranty period. For parts from outside manufacturers, insofar as these constitute an appreciable part of our instrument, the warranty stipulations of the manufacturer in question apply.

With the regard to the guarantee of accuracy the technical specifications in the instruction manual are authoritative.

Concerning defects in materials, construction or design as well as the absence of guaranteed features the purchaser has no rights or claims except those mentioned above.

If damage of the packaging is evident on receipt of a consignment or if the goods show signs of transport damage after unpacking, the carrier must be informed immediately and a written damage report demanded. Lack of an official damage report releases Metrohm from any liability to pay compensation.

If any instruments and parts have to be returned then the original packaging should be used if at all possible. This applies above all to instruments, electrodes, buret cylinders and PTFE pistons. Before embedding in wood shavings or similar material the parts must be packed in a dustproof package (for instruments the use of a plastic bag is essential). If open assemblies are included that are sensitive to electromagnetic voltages (e.g. data interfaces, etc.) then these must be returned in the associated original protective packaging (e.g. conductive protective bag). (Exception: assemblies with a built-in voltage source belong in non-conductive protective packaging).


For damage which arises as a result of non-compliance with these instructions no warranty responsibility whatsoever will be accepted by Metrohm.

Declaration of Conformity

This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.

CH-9101 Herisau, Switzerland
Tel. +41 71 353 85 85
Fax +41 71 353 89 01
www.metrohm.com

<i>Name of commodity</i>	Differential amplifier
<i>System software</i>	Stored in ROMs
<i>Name of manufacturer</i>	Metrohm Ltd., Herisau, Switzerland
<i>Description</i> <i>Device for titrators without built-in differential amplifier to make them usable for differential potentiometry.</i>	
This Metrohm instrument has been built and has undergone final type testing according to the standards: <i>Electromagnetic compatibility: Emission</i> EN/IEC 61326, EN 55022 / CISPR 22 <i>Electromagnetic compatibility: Immunity</i> EN/IEC 61326, EN/IEC 61000-4-2, EN/IEC 61000-4-3, EN/IEC 61000-4-4, EN/IEC 61000-4-5, EN/IEC 61000-4-6, EN/IEC 61000-4-8, EN/IEC 61000-4-11, EN/IEC 61000-4-14 <i>Safety specifications</i> EN/IEC 61010-1, UL3101-1 <i>The instrument meets the requirements of the CE mark as contained in the EU directives 89/336/EEC and 73/23/EEC and fulfils the following specifications:</i>	
EN 61326	Electrical equipment for measurement, control and laboratory use – EMC requirements
EN 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use
Metrohm Ltd. Is holder of the SQS-certificate of the quality system ISO 9001 for quality assurance in design/development, production, installation and servicing. The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality and performance. The technical specifications are documented in the instruction manual.	
Herisau, December 20, 2002  Dr. J. Frank Vice President Head of R&D  Ch. Buchmann Vice President Head of Production Responsible for Quality Assurance	

Quality Management Principles

Metrohm Ltd., CH-9101 Herisau, Switzerland

 **Metrohm**
l o n a n a l y s i s
CH-9101 Herisau/Switzerland
E-Mail info@metrohm.com
Internet www.metrohm.com

Metrohm Ltd. Holds the ISO 9001 Certificate, registration number 10872-02, issued by SQS (Swiss Association for Quality and Management Systems). Internal and external audits are carried out periodically to assure that the standards defined by Metrohm's QM Manual are maintained.

The steps involved in the design, manufacture and servicing of instruments are fully documented and the resulting reports are archived for ten years. The development of software for PCs and instruments is also duly documented and the documents and source codes are archived. Both remain the possession of Metrohm. A non-disclosure agreement may be asked to be provided by those requiring access to them.

The implementation of the ISO 9001 quality system is described in Metrohm's QM Manual, which comprises detailed instructions on the following fields of activity:

Instrument development

The organisation of the instrument design, its planning and the intermediate controls are fully documented and traceable. Laboratory testing accompanies all phases of instrument development.

Software development

Software development occurs in terms of the software life cycle. Tests are performed to detect programming errors and to assess the program's functionality in a laboratory environment.

Components

All components used in the Metrohm instruments have to satisfy the quality standards that are defined and implemented for our products. Suppliers of components are audited by Metrohm as the need arises.

Manufacture

The measures put into practice in the production of our instruments guarantee a constant quality standard. Production planning and manufacturing procedures, maintenance of production means and testing of components, intermediate and finished products are prescribed.

Customer support and service

Customer support involves all phases of instrument acquisition and use by the customer, i.e. consulting to define the adequate equipment for the analytical problem at hand, delivery of the equipment, user manuals, training, after-sales service and processing of customer complaints. The Metrohm service organisation is equipped to support customers in implementing standards such as GLP, GMP, ISO 900X, in performing Operational Qualification and Performance Verification of the system components or in carrying out the System Validation for the quantitative determination of a substance in a given matrix.

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
<i>Name of commodity</i>	Amplificateur différentiel
<i>System software</i>	Stored in ROMs
<i>Name of manufacturer</i>	Metrohm Ltd., Herisau, Switzerland
<i>Description</i>	<i>Device for titrators without built-in differential amplifier to make them usable for differential potentiometry.</i>
This Metrohm instrument has been built and has undergone final type testing according to the standards: <i>Electromagnetic compatibility: Emission</i> EN/IEC 61326, EN 55022 / CISPR 22 <i>Electromagnetic compatibility: Immunity</i> EN/IEC 61326, EN/IEC 61000-4-2, EN/IEC 61000-4-3, EN/IEC 61000-4-4, EN/IEC 61000-4-5, EN/IEC 61000-4-6, EN/IEC 61000-4-8, EN/IEC 61000-4-11, EN/IEC 61000-4-14 <i>Safety specifications</i> EN/IEC 61010-1, UL3101-1 <i>The instrument meets the requirements of the CE mark as contained in the EU directives 89/336/EEC and 73/23/EEC and fulfils the following specifications:</i>	
EN 61326	Electrical equipment for measurement, control and laboratory use – EMC requirements
EN 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use
Metrohm Ltd. is holder of the SQS-certificate of the quality system ISO 9001 for quality assurance in design/development, production, installation and servicing. The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality and performance. The technical specifications are documented in the instruction manual.	
Herisau, December 20, 2002  Dr. J. Frank Vice President Head of R&D  Ch. Buchmann Vice President Head of Production Responsible for Quality Assurance	

Quality Management Principles

Metrohm Ltd., CH-9101 Herisau, Switzerland

 **Metrohm**
l o n a n a l y s i s
CH-9101 Herisau/Switzerland
E-Mail info@metrohm.com
Internet www.metrohm.com

Metrohm Ltd. holds the ISO 9001 Certificate, registration number 10872-02, issued by SQS (Swiss Association for Quality and Management Systems). Internal and external audits are carried out periodically to assure that the standards defined by Metrohm's QM Manual are maintained.

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
<i>Name of commodity</i>	Amplificador diferencial
<i>System software</i>	Stored in ROMs
<i>Name of manufacturer</i>	Metrohm Ltd., Herisau, Switzerland
<i>Description</i>	<i>Device for titrators without built-in differential amplifier to make them usable for differential potentiometry.</i>
This Metrohm instrument has been built and has undergone final type testing according to the standards: <i>Electromagnetic compatibility: Emission</i> EN/IEC 61326, EN 55022 / CISPR 22 <i>Electromagnetic compatibility: Immunity</i> EN/IEC 61326, EN/IEC 61000-4-2, EN/IEC 61000-4-3, EN/IEC 61000-4-4, EN/IEC 61000-4-5, EN/IEC 61000-4-6, EN/IEC 61000-4-8, EN/IEC 61000-4-11, EN/IEC 61000-4-14 <i>Safety specifications</i> EN/IEC 61010-1, UL3101-1 <i>The instrument meets the requirements of the CE mark as contained in the EU directives 89/336/EEC and 73/23/EEC and fulfils the following specifications:</i>	
EN 61326	Electrical equipment for measurement, control and laboratory use – EMC requirements
EN 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use
Metrohm Ltd. is holder of the SQS-certificate of the quality system ISO 9001 for quality assurance in design/development, production, installation and servicing. The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality and performance. The technical specifications are documented in the instruction manual.	
Herisau, December 20, 2002  Dr. J. Frank Vice President Head of R&D  Ch. Buchmann Vice President Head of Production Responsible for Quality Assurance	

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