

IC equipment



IC equipment: Dose-in Gradient (6.5330.150)

Manual

8.110.8032EN / v2 / 2024-07-02



Metrohm AG
CH-9100 Herisau
Switzerland
+41 71 353 85 85
info@metrohm.com
www.metrohm.com

IC equipment

IC equipment: Dose-in Gradient (6.5330.150)

Manual

Technical Communication
Metrohm AG
CH-9100 Herisau

This documentation is protected by copyright. All rights reserved.

This documentation is an original document.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Disclaimer

Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc., are expressly excluded from the warranty. Unauthorized modifications to the product (e.g., conversions or attachments) exclude any liability on the part of the manufacturer for resulting damage and its consequences. Instructions and notes in the Metrohm product documentation must be strictly followed. Otherwise, Metrohm's liability is excluded.

[illegible][illegible]



Table of figures

Figure 1	Overview Dose-in Gradient	1
Figure 2	IC equipment: Dose-in Gradient – Parts	4



1 Introduction

1.1 Description of the IC equipment: Dose-in Gradient

You can upgrade an isocratic IC system to a gradient system very easily with an IC equipment: Dose-in Gradient. A gradient is created by mixing different eluents at a mixing ratio that varies over the analysis time.

An 800 Dosino doses a concentrated or diluted eluent into the main eluent stream at a defined ratio. This changes the total concentration of the mobile phase. The result is an optimum separation of all ions within an acceptable analysis time.

What sets the Dose-in Gradient apart, is that it can be built into any existent Metrohm IC system and can be used without restrictions. It can, of course, also be included in the system setup from the beginning.

Another function of the Dose-in Gradient allows you to add a rinsing solution without modifying the system. It rinses the system after the analysis or after a sample series and thereby increases the service life of the individual components.




Figure 1 Overview Dose-in Gradient



WARNING

Warning of optical radiation



CAUTION

This symbol draws attention to possible damage to instruments or instrument parts.



NOTICE

This symbol highlights additional information and tips.

2 Overview

2.1 Parts of the IC equipment: Dose-in Gradient




Figure 2 IC equipment: Dose-in Gradient – Parts

To create an eluent gradient, in addition to the IC equipment: Dose-in Gradient, you need:

- an 800 Dosino (2.800.0010)
- a suitable 807 Dosing Unit (6.1580.XX0)

3 Installation


Proceed as follows to install the IC equipment: Dose-in Gradient:

- Equip the eluent bottle.
- Assemble the 800 Dosino and connect all tubings.

3.1 Equipping the eluent bottle

Required accessories

- Bottle (6.1608.070) filled with the additional eluent
- Eluent bottle cap (6.1602.160)
- Adsorber tube (6.1619.000)
- Adapter for adsorber tube (6.1624.000)
- FEP aspiration tubing (6.1829.020)
- M8 stopper (6.1446.080), included in the accessories for the eluent bottle cap (6.1602.160)



1 Mounting the aspiration tubing

- Remove the clamping screw from the aspiration tubing.
- Insert the aspiration tubing into the M6 opening of the eluent bottle cap from above.
- Use the capillary cutter to cut the aspiration tubing to such a length that it touches the bottom of the bottle.

2 Inserting the stopper

- Tighten the M8 stopper in the M8 opening of the eluent bottle cap.

3 Mounting the adsorber tube

- Fill the adsorber tube with some cotton and adsorber material.
- Place the adsorber tube onto the adapter.
- Insert the adapter into the SGJ opening of the eluent bottle cap.

4 Mounting the eluent bottle cap

- Screw the eluent bottle cap onto the bottle filled with the additional eluent.

3.2 Mounting the Dosino

Attaching the Dosino to the 807 Dosing Unit

Required accessories

- 800 Dosino (2.800.0010)
- 807 Dosing Unit without accessories (6.1580.XX0)



CAUTION


Please read through the correct procedure in the Manual for the 800 Dosino before you attach the Dosino to the 807 Dosing Unit.

- 1 Attach the Dosino to the 807 Dosing Unit, (see *Manual for the 800 Dosino*).

Fastening the Dosino to the ion chromatograph

Required accessories

- Dosino (2.800.010) with 807 Dosing Unit without accessories (6.1580.XX0)
- Dosino holder (6.2057.210)
- Thread adapter (6.1618.020)
- 3 mm hex key (6.2621.100)



1 Fitting the Dosino holder onto the ion chromatograph

- Loosen the bottle holder on the ion chromatograph using the hex key.
- Clamp the Dosino holder under it sideways.
- Fasten the bottle holder again.

2 Attaching the Dosino to the holder

- Place the Dosino onto the Dosino holder.
- Fasten the Dosino to the Dosino holder by tightening the thread adapter from below.

3 Connecting the Dosino to the ion chromatograph



NOTE

The ion chromatograph **must** be **switched off** when the Dosino is being plugged into the MSB connector.

- Check whether the ion chromatograph is switched on.
If this is the case, switch off the ion chromatograph.
- Plug the Dosino cable into one of the ion chromatograph's MSB connectors.

3.3 Mounting the FEP tubing

Required accessories

- FEP tubing (6.1805.120)




- 1
 - Tighten one end of the FEP tubing in the M6 opening of the eluent bottle cap.
 - Tighten the other end of the FEP tubing in port 2 of the Dosino.

3.4 Mounting the PEEK capillary

Required accessories

- PEEK capillary, 0.75 mm ID (6.1831.030)
- Pressure screws (6.2744.070)
- Coupling M6 / UNF (6.2744.080)



- 1 ■ Tighten the coupling to Port 1 of the Dosino.

- 2 ▪ Tighten one end of the PEEK capillary to the coupling using a pressure screw.

3.5 Installing the 6-port adapter


The 6-port adapter functions as mixing valve for the eluent gradient.

Proceed as follows to connect the 6-port adapter:

Connecting the 6-port adapter

Required accessories

- 6-port adapter (6.2744.340)
- Tubing connection degasser – pump (6.1834.090)
- PEEK pressure screw, 14 mm (6.2744.350)
- PEEK capillary, 0.75 mm ID (6.1831.030), screwed onto the Dosino (see chapter 3.4, page 8)




- 1
 - Unscrew the tubing connection between the eluent degasser and the high-pressure pump at the inlet of the high-pressure pump.
 - Tighten the end of the tubing connection to a UNF ¼"-28 port of the 6-port adapter.

- 3

- ## Mounting the 6-port adapter

Self-adhesive hook-and-loop fastening tape (6.2053.030)

We recommend fastening the 6-port adapter to the high-pressure pump using the self-adhesive hook-and-loop fastening tape (6.2053.030).



- 1 Use scissors to cut the self-adhesive hook-and-loop fastening tape in half. (You need only one half to fasten the 6-port adapter.)
- 2 Tear the hook-and-loop fastener apart.
- 3 Stick one side onto the right half of the high-pressure pump.



NOTE

Make sure that the fastening screws of the high-pressure pump are not covered.

- 4 Stick the second side of the hook-and-loop fastening tape onto the 6-port adapter.
- 5 Fasten the 6-port adapter to the high-pressure pump.



NOTE

There is not much space between the high-pressure pump and the door. Therefore:

- mount the 6-port adapter on the far right.
- adjust the 6-port adapter so that the attached tubings and capillaries will not kink.


3.6 Installing the trap column

The high-capacity anion column is used to purify and mix the eluent flow. Even reagents of the highest quality level, e.g. "Ultrapur" or "puriss", can still contain minimal anion contaminations. These are reliably held back by the Metrosep A Trap 1 - 100/4.0. In the application described here, however, the trap column functions as mixer and ensures thorough mixing of the different media.

Install the trap column in place of the inline filter between the high-pressure pump and the pulsation absorber.

Removing the inline filter

1




Remove the inline filter.

Installing the trap column


Required accessories

- Metrosep A Trap 1 - 100/4.0 (6.1014.000)
- PEEK capillary, 0.25 mm ID (6.1831.120)
- Pressure screws (6.2744.070)

1




Remove the stoppers of the trap column.



The diagram shows a spray gun assembly. It includes a spray gun (a long, thin nozzle), a spray gun handle (a rectangular block with a trigger), and a spray gun trigger (a small, cylindrical component). An orange arrow points from the spray gun handle to the spray gun trigger, indicating the flow of material or the direction of use.

3




- Tighten one end of the PEEK capillary to the trap column inlet using a pressure screw.
- Tighten the other end of the capillary to the purge valve using a pressure screw.

4.1 807 Dosing Unit without accessories (6.1580.XX0)

5 Displaying accessories

Up-to-date information on the scope of delivery and on optional accessories can be found on the Metrohm website.

1 Searching for a product on the website

- Go to <https://www.metrohm.com>.
- Click on .
- Enter the article number of the product (e.g. **2.1001.0010**) into the search field and press **[Enter]**.

The search result is displayed.

2 Displaying product information

- To display the products matching the search term, click on **Product models**.
- Click on the desired product.

Detailed information regarding the product is displayed.

3 Displaying accessories and downloading the accessories list

- To display the accessories, scroll down to **Accessories and more**.
 - The **scope of delivery** is displayed.
 - Click on **[Optional parts]** for the optional accessories.
- To download the accessories list, click on **[Download accessories PDF]** under **Accessories and more**.



NOTE

Metrohm recommends keeping the accessories list for reference purposes.

E

Eluent
for gradient 5

6-port adapter	9
Bottle cap	5
Dosino	6
FEP tubing	8
Trap column	12

T	
Trap column	
Install	12