

Metrohm


Instructions for use

Sample Changer 676.0300

Series 06...

Control Unit 664.0020

Series 12...



Measurement in Chemistry

METROHM Ltd.
CH-9101 Herisau

Switzerland
Phone 071/53 11 33
Telefax 071/52 11 14
Telex 8827 12 Metro ch

8.676.1003

Sample Changer 676.0300 Series 06 . . .

Control Unit 664.0020 Series 12 . . .

Mains connection:

Mains voltage	U	= 100, 117, 220, 240 V \pm 10%
Mains frequency	f	= 50 ... 60 Hz
Power consumption	S	= 50 VA

8.676.1003 Instructions for Use

91.06 / dō

8.676.1003 Instructions for Use 676 Sample Changer

Table of Contents

	<u>Page</u>
1. Introduction	1
2. Control elements	2
2.1. 676 Sample Changer	2
2.2. 664.0020 Control Unit	4
3. Installation	7
3.1. Setting up the instruments	7
3.1.1. Packing	7
3.1.2. Check	7
3.1.3. Location	7
3.2. Equipping the titration head at the 676 Sample Changer	7
3.2.1. Equipping for the determination of hydroxyl and saponification numbers	7
3.2.2. Equipping for the COD determination	8
3.3. Tubing connections	9
3.3.1. Tubing connections for rinsing with 683 Pump	9
3.3.2. Tubing connections for rinsing with air pump and electromagnetic valve	9
3.3.3. Tubing connections for siphoning with the 683 Pump	9
3.3.4. Fastening the 6.1812.000 PTFE Tubing	10
3.3.5. Fastening tubing and cables	10
3.4. Mains connection and switching on the 664.0020 Control Unit	10
3.4.1. Setting the voltage	10
3.4.2. Mains cable	10
3.4.3. Mains connection	10
3.4.4. Switching on the instrument	10
3.5. Electrical interconnections	11
3.5.1. Interconnection with the 686 or 682 Titroprocessor	11
3.5.2. Interconnection with the 702 SET/MET Titrino	12
3.5.3. Interconnection with the 670 Titroprocessor	13
3.6. Settings at 664.0020 Control Unit	14
3.6.1. Program selection	14
3.6.2. Setting the delay and valve opening times	15
3.7. Filling the tubing	16
3.8. Setting the stirring speed	16
3.9. Conditioning	16
4. Operation	17
4.1. Operation scheme	17
4.2. Manual operation	18
4.2.1. Manual operation rinsing	18
4.2.2. Manual operation siphoning	18
4.2.3. Manual operation of the magnetic stirrer	18
4.2.4. Manual operation of the lifting arm	18
4.2.5. Manual operation of the sample transport	18
4.3. Notes on operation of the 676 Sample Changer	19
4.3.1. Sample rack	19
4.3.2. Stirring bars	19
4.3.3. Check on the electrolyte vessel	19
4.3.4. Start of the analysis program	19
4.3.5. Manual stop of the analysis program	19
4.3.6. Conditioning solution	19
4.3.7. Siphoning off with 683 Pump	19

	<u>Page</u>
5. Application examples	20
5.1. Determination of the hydroxyl number	20
5.1.1. Analysis procedure	20
5.1.2. Titration with the 686 or 682 Titroprocessor	20
5.2. Determination of the saponification number	21
5.2.1. Analysis procedure	21
5.2.2. Titration with the 686 or 682 Titroprocessor	21
5.3. COD determination	22
5.3.1. Analysis procedure	22
5.3.2. COD determination with the 686 Titroprocessor	22
5.3.3. COD determination with the 670 Titroprocessor	22
5.3.4. COD determination with the 702 SET/MET Titrino	24
6. Technical specifications	25
6.1. 676 Sample Changer	25
6.2. 664.0020 Control Unit	26
6.3. 6.2726.070 Electromagnetic Valve	27
6.4. 6.2817.000 Air Pump	27
6.5. Control inputs and outputs	28
7. Scope of delivery and ordering designations	33
7.1. 676 Sample Changer	33
7.2. 664.0020 Control Unit	33
7.3. Options	33
7.3.1. Reaction vessels	33
7.3.2. Accessories for the determination of hydroxyl and saponification numbers	33
7.3.3. Accessories for the COD determination	34
7.3.4. Accessories for 676 Sample Changer	34
7.3.5. Accessories for 664.0020 Control Unit	35
7.3.6. Auxiliary instruments	35
8. Warranty	36
9. Index	37

List of Figures

Fig. 1: 676 Sample Changer with 664 Control Unit, 686 Titroprocessor, 665 Dosimat and two 683 Pumps	1
Fig. 2: Front of 676 Sample Changer	2
Fig. 3: Rear of 676 Sample Changer	3
Fig. 4: Front of 664.0020 Control Unit	4
Fig. 5: Rear of 664.0020 Control Unit	4
Fig. 6: Top view of 664.0020 Control Unit with top panel removed	6
Fig. 7: Equipping the titration head	8
Fig. 8: Rinsing and siphoning the sample vessels	9
Fig. 9: Interconnection with 686/682 Titroprocessor and 665 Dosimat for the determination of hydroxyl and saponification numbers	11
Fig. 10: Interconnection with the 702 SET/MET Titrino for the COD determination	12
Fig. 11: Interconnection with the 670 Titroprocessor for the determination of hydroxyl and saponification numbers	13

676 Sample Changer with 664 Control Unit

Instructions for Use

1. Introduction

The **676 Sample Changer** allows series analyses in digestion tubes from the Behr, Barkey or Merck companies to be performed automatically. It is thus eminently suitable for all applications in which the sample has to be boiled or digested before the determination. The detachable sample rack holds 15 sample vessels and can be inserted directly in a suitable digestion block.

The 676 Sample Changer is normally controlled by the associated **644.0020 Control Unit**, which ensures not only coordinated sample transport, but also rational use of the required instruments and auxiliary equipment such as titrators and dispensers, valves and pumps in each particular case. To fulfill this task, the control unit has two special, built-in programs with the flow sequence of the analysis to be performed. The delay times incorporated in the analysis routine as well as the opening times of the valves can be entered directly on the control unit.

The **670 Titroprocessor** in conjunction with a **677 Drive Unit** can also be used to control the 676 Sample Changer in place of the 664 Control Unit. Special program commands in this case ensure sample transport, stirring, siphoning off and rinsing of the sample vessels.

The present Instructions for Use provide a detailed description of the above-mentioned possibilities for operation of the 676 Sample Changer. In addition, examples are used to illustrate the **determination of hydroxyl and saponification numbers** and the **COD determination** with various METROHM titrators.




Fig. 1: 676 Sample Changer with 664 Control Unit, 686 Titroprocessor, 665 Dosimat and two 683 Pumps

2. Control elements

2.1. 676 Sample Changer




Fig. 2: Front of 676 Sample Changer




Fig. 3: Rear of 676 Sample Changer

- | | |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>1 Indicator lamp of stirrer</p> <p>2 Rotary regulator for control of the stirring speed
rotational speed (stabilised):
200 ... 1900 /min</p> <p>3 Glide surface</p> <p>4 6.2041.030 Sample rack
for 15 sample vessels</p> <p>5 4.676.0700 Holding bracket for sample rack</p> <p>6 Sample vessel</p> <ul style="list-style-type: none"> • sample vessel Barkey/Behr:
6.1452.000 • sample vessel Merck:
6.1452.010 <p>7 Lift tower</p> <p>8 Burette tip</p> <ul style="list-style-type: none"> • for 6.1452.000 Sample vessels (Barkey/Behr):
6.1543.130 Burette tip • for 6.1452.010 Sample vessels (Merck):
6.1543.140 Burette tip <p>with connection for 6.1805.110 FEP Tubing
(supply of the titration solution from 665 Dosi-mat or 702 SET/MET Titrino)</p> <p>9 Connection for siphoning tubing
Connecting nipple for 6.1812.000 PTFE siphoning tubing
(connection to 683 Siphoning Pump)</p> <p>10 Lifting arm
with integral titration head</p> <p>11 Nipple with O-ring</p> <ul style="list-style-type: none"> • for mounting the 6.9902.048 pH Electrode (determination of hydroxyl and saponification numbers) • for mounting the 6.1240.050 or 6.1240.060 Electrode Vessels (COD determination) | <p>12 Connection for rinsing solution
Connecting nipple for 6.1812.000 PTFE Tubing
(connection to 683 Rinsing Pump or 6.2726.070 Electromagnetic Valve)</p> <p>13 Titration head</p> <p>14 6.2740.010 Rinse nozzle</p> <p>15 Siphoning tube</p> <ul style="list-style-type: none"> • for 6.1452.000 Sample vessels (Barkey/Behr):
6.1821.000 Siphoning tube • for 6.1452.010 Sample vessels (Merck):
6.1821.010 Siphoning tube <p>16 Turntable
with holder for conditioning vessel</p> <p>17 Connecting cable of electromagnetic valve
Connecting cable to 664.0020 Control Unit</p> <p>18 6.2726.070 Electromagnetic valve</p> <p>19 Electromagnetic valve outlet
Connecting nipple for 6.1812.000 PTFE tubing</p> <p>20 Electromagnetic valve inlet (not visible)
Connecting nipple for 6.1812.000 PTFE tubing</p> <p>21 Valve panel for 4 valves</p> <p>22 Connecting cable to 664.0020 Control Unit
integral</p> |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

2.2. 664.0020 Control Unit




Fig. 4: Front of 664.0020 Control Unit
(the parts not marked have no function)




Fig. 5: Rear of 664.0020 Control Unit
(the parts not marked have no function)

- 23 Power on lamp**
lights up when instrument switched on
- 24 Manual operation Siphoning**
Pressing this key *down* switches on the 683 Siphoning Pump (this must be attached to connection **34-2**).
- 25 Manual operation Rinsing**
Pressing this key *down* triggers the rinsing operation (683 Rinsing pump is switched on or electromagnetic valve **18** is opened). For this, the 683 Rinsing Pump or electromagnetic valve **18** must be attached to connection **34-3**.
- 26 Manual operation of the stirrer**
Pressing this key *up* switches on the magnetic stirrer at the sample changer (indicator lamp **1** lights up). The stirring speed can be set at the rotary regulator **2**.
- 27 Manual operation of the lifting arm**
By pressing this key *up*, lifting arm **10** is raised.

By pressing this key *down*, lifting arm **10** is lowered.
- 28 Manual operation of the sample transport**
By pressing this key *down*, operation of the turntable **16** is started.

If the sample vessels need to be moved only by 1 position, the key must be pressed only briefly (ca. 1 s).
- 29 Stop lamp**
lights up when the analysis program has been executed or stopped manually.
- 30 Mains switch**
On/off switching of the instrument

When instrument is in operation, the power-on lamp **23** is lit up
- 31 Thumbwheel switch to set the number of samples**
must be set to 00
- 32 Control key**
for starting and stopping the analysis program

By pressing the key *down*, the set analysis program is *started*.

By pressing the key *up*, the set analysis program is immediately *stopped* (an ongoing titration is first run to completion).
- 33 Connecting socket for stirrer**
- 34 Plug socket J with connecting sockets for siphoning pump (V2) and rinsing pump or electromagnetic valve (V3)**
- 35 Display of the set mains voltage**
Setting of the voltage, see section 3.4
- 36 Important note:**
"In case of service note the instructions inside this panel"
- 37 Instrument designation**
Instrument number, model number, serial number
- 38 3.540.2111 PCB Connector for connection of the 676 Sample Changer**
- 39 3.540.1640 PCB Connector for connection of peripheral units (e.g. 682, 686, 702)**
- 40 Important note:**
*"ATTENTION:
For continued protection replace only with the same type and rating of fuse!"*
- 41 Earth socket**
Earthing, see section 3.4
- 42 Mains socket**
Mains connection, see section 3.4
- 43 Fuse and power consumption data**
- 44 Fusible cutout**
Changing the fuse, see section 3.4




Fig. 6: Top view of 664.0020 Control Unit with top panel removed
(the parts not marked have no function)

45 Adjustment controls
for valve opening times V1, V2, V3
(procedure see section 3.6.2)

47 Adjustment controls
for delay times Δt_1 , Δt_2 , Δt_3
(procedure see section 3.6.2)

46 Designation of direction of rotation
of potentiometer

48 Selector switch for analysis program
4 or 5

49 Screw drivers / contact pins

3. Installation

3.1. Setting up the instruments

3.1.1. Packing

The 676 Sample Changer and 664 Control Unit are supplied together with the separately packed accessories in special packaging which provides exceptionally good protection. This comprises shock-absorbing foam liners in blue plastic film foamed to match the individual shape. The instruments themselves are packed in a dust-tight, evacuated polyethylene bag. It is essential that you store all these special packagings in a safe place as they and they alone assure indemnified transport of the instruments.

3.1.2. Check

On receipt of the goods, immediately check that the shipment is complete and is free from damage (compare with delivery note and accessories list in *section 8*). If transport damage is found, see instructions in *section 9 "Warranty"*.

3.1.3. Location


Set up the 676 Sample Changer and the 664 Control Unit at a place in the laboratory suitable for operation which is free from vibration, protected against any corrosive atmosphere and contamination by chemicals. For perfect functioning, the sample changer needs a horizontal surface. If none is available, the sample changer must be levelled with the aid of shims (check horizontal positioning with a level).

3.2. Equipping the titration head at the 676 Sample Changer

3.2.1. Equipping for the determination of hydroxyl and saponification numbers

Titration head **13** is adjusted in the factory for either Barkey/Behr or Merck sample vessels in accordance with your ordering instructions. The additional accessories required for each type of vessel (see *section 8.3.2*) are mounted in the titration head as follows (see *Figs. 2 and 7*):

- 1▶ Carefully (slowly) raise lifting arm **10** with permanently mounted titration head **13** by hand (*the 664 Control Unit must not be switched on!*).
- 2▶ Screw rinse nozzle **14** (6.2740.010) from below into opening **54** of the titration head.
- 3▶ Screw siphoning tube **15** (6.1821.000 for Barkey/Behr vessels; 6.1821.010 for Merck vessels) from below into opening **55** of the titration head so that the tip of the bevelled tube faces outward (towards the edge of the vessel).
- 4▶ Insert burette tip **8** (6.1543.130) from above into opening **50** of the titration head until the lower end is ca. 1 cm above the end of siphoning tube **15**.
- 5▶ Screw 6.1805.110 FEP Tubing to burette tip **8** (connection to 665 Dosimat).
- 6▶ Insert 6.9902.048 pH Microelectrode from above into opening **52** of the titration head until the lower end is ca. 1 cm above the end of siphoning tube **15**.
- 6▶ Fasten 6.9902.048 pH Microelectrode tightly with nipple **11** and fit with 6.2104.030 Electrode Cable (or 6.2104.060 when 670 Titroprocessor used).



Determination of hydroxyl and saponification number

50	6.1543.1X0 Burette tip
51	-
52	6.9902.048 pH Microelectrode
53	-
54	6.2740.010 Rinse nozzle
55	6.1821.0X0 Siphoning tube

COD determination

50	6.1543.1X0 Burette tip
51	6.1248.0X0 Gold rod electrode
52	6.1240.0X0 Electrolyte vessel with 6.0724.140 Reference system
53	-
54	6.2740.010 Rinse nozzle
55	6.1821.0X Siphoning tube

Fig. 7: Equipping the titration head

3.2.2. Equipping for the COD determination

Titration head **13** is adjusted in the factory for either Barkey/Behr or Merck sample vessels in accordance with your ordering instructions. The additional accessories required for each type of vessel (see section 8.3.3) are mounted in the titration head as follows (see Figs 2 and 8):

- 1▶ Carefully (slowly) raise lifting arm **10** with permanently mounted titration head **13** by hand (*the 664 Control Unit must not be switched on!*).
- 2▶ Screw rinse nozzle **14** (6.2740.010) from below into opening **54** of the titration head.
- 3▶ Screw siphoning tube **15** (6.1821.000 for Barkey/Behr vessels; 6.1821.010 for Merck vessels) from below into opening **55** of the titration head so that the tip of the bevelled tube faces outward (towards the edge of the vessel).
- 4▶ Insert 6.1248.0X0 Gold rod electrode (6.1248.060 for Barkey/Behr vessels; 6.1248.070 for Merck vessels) from above into opening **51** of the titration head until the lower end is ca. 1 cm above the stop and fit with 6.2120.030 Electrode cable.
- 5▶ Insert burette tip **8** (6.1543.130) from above into opening **50** of the titration head until the lower end is ca. 1 cm above the end of siphoning tube **15**.
- 6▶ Screw 6.1805.110 FEP Tubing to burette tip **8** (connection to 665 Dosimat or 702 SET/MET Titrino).
- 7▶ Fill 6.1240.0X0 Electrolyte Vessel (6.1240.050 for Barkey/Behr vessels; 6.1240.060 for Merck vessels) with $c(\text{KNO}_3) = \text{sat.}$ (6.2310.000 Electrolyte Solution), insert into opening **52** from above until the bottom end is ca. 2 cm above the end of siphoning tube **15** and then fasten tightly with nipple **11**.
- 8▶ Mount 6.0724.140 Reference System on 6.1240.0X0 Electrolyte Vessel and fit with 6.2106.060 Cable.

3.3. Tubing connections

3.3.1. Tubing connections for rinsing with 683 Pump

For rinsing of the sample vessels, use of a 683 Pump as shown in scheme **A** in Fig. 8 is recommended. The following connections must be made in such a case:

- 1▶ Connect connecting nipple **12** of the 676 Sample Changer (see Fig. 2) with the aid of the 6.1812.000 PTFE Tubing to connecting nipple **10** of the 683 Pump (see Instructions for Use of 683 Pump) (procedure, see section 3.3.4).
- 2▶ Attach 6.1812.000 PTFE Tubing to connecting nipple **6** of the 683 Pump (see Instructions for Use of 683 Pump) and insert in rinsing agent supply bottle (procedure, see section 3.3.4).

3.3.2. Tubing connections for rinsing with air pump and electromagnetic valve

The sample vessels can also be rinsed as shown in scheme **B** in Fig. 8 with the aid of a 6.2817.000 Air Pump and a 6.2726.070 Electromagnetic Valve. In this case, the supply bottle must be kept constantly under pressure by the air pump. When the electromagnetic valve opens, rinsing agent is delivered. For this type of rinsing, the following connections must be made:

- 1▶ Connect connecting nipple **12** of the 676 Sample Changer (see Fig. 2) with the aid of the 6.1812.000 PTFE Tubing to connecting nipple **19** of the electromagnetic valve **18** (see Fig. 3) (procedure, see section 3.3.4).
- 2▶ Attach 6.1812.000 PTFE Tubing to connecting nipple **20** of electromagnetic valve **18** (see Fig. 3) and insert in a compression-proof rinsing agent supply bottle (procedure, see section 3.3.4)
- 3▶ Connect 6.2817.000 Air Pump with the aid of the appropriate connecting tubing to the rinsing agent supply bottle.

3.3.3. Tubing connections for siphoning with the 683 Pump

For siphoning off of the sample vessel contents in accordance with scheme **C** in Fig. 8, the following connections must be made:

- 1▶ Connect connecting nipple **9** of the 676 Sample Changer (see Fig. 2) with the aid of the 6.1812.000 PTFE Tubing to connecting nipple **6** of the 683 Pump (see Instructions for Use of 683 Pump) (procedure, see section 3.3.4).
- 2▶ Attach 6.1812.000 PTFE Tubing to connecting nipple **10** of the 683 Pump (see Instructions for Use of 683 Pump) and insert in waste container (procedure, see section 3.3.4).

Note: *No solids can be siphoned off with the 683 Pump.
These could damage the pump.*




Fig. 8: Rinsing and siphoning the sample vessels

- | | |
|--------------------------|---------------------------------------|
| 1 Supply bottle | 5 6.2817.000 Air pump |
| 2 6.1812.000 PTFE Tubing | 6 6.2726.070 Electromagnetic valve V3 |
| 3 683 Pump | 7 Waste container |
| 4 Sample vessel | |

3.3.4. Fastening the 6.1812.000 PTFE Tubing

To fasten the 6.1812.000 PTFE Tubing to the appropriate connecting nipples, a suitable procedure is as follows:

- 1▶ Place the screw cap over the tubing end to be fastened.
- 2▶ Enlarge tubing using a pencil (or a similar object).
- 3▶ Grip tubing with fine abrasive paper and push over the connecting nipple.
- 4▶ Tighten screw cap firmly by hand (use no tools!)

3.3.5. Fastening tubing and cables

It is advisable to bind together all tubing as well as all electrode cables attached to titration head **13** using the optional 6.1815.010 Spiral Band.

3.4. Mains connection and switching on the 664.0020 Control Unit

3.4.1. Setting the voltage

Before switching on the 664 Control Unit for the first time, check whether the mains voltage set on the instrument (visible in display **36**) matches the available mains voltage. If this is not the case, you have to change the voltage as follows:

- 1▶ Disconnect mains cable and all other cables plugged into the control unit.
- 2▶ Unscrew instrument rear panel and take off.
- 3▶ Turn voltage selector switch with a screwdriver to the desired position.
- 4▶ Replace instrument rear panel and screw on.
- 5▶ Take out fuse **44** and compare its specifications with the data **43** on the rear panel.
Replace fuse if necessary.
220 ... 240 V: 0.25 A (slow-blow)
110 ... 117 V: 0.5 A (slow-blow)
- 6▶ Reconnect cables.

3.4.2. Mains cable

The requested mains cable supplied with the instrument:

- 6.2122.020 mit Stecker SEV 12 (Switzerland ...)
- 6.2122.040 mit Stecker CEE(7), VII (Germany ...)
- 6.2122.070 mit Stecker NEMA 5-15 (USA ...)

is three-core and fitted with a connector with an earthing pin. If a different connector has to be mounted, the yellow/green lead (IEC standard) must be connected to the protective earth (Class 1 instrument). If no socket with earthing is available, you have to connect the instrument to a perfect earthing conductor via its earth socket **41** (–).

3.4.3. Mains connection

Plug the mains cable into connection socket **42** of the 664 Control Unit.

3.4.4. Switching on the instrument

Switch on the 664 Control Unit with mains switch **30**.

When the power-on lamp **23** is lit up, the instrument is switched on.

3.5. Electrical interconnections

3.5.1. Interconnection with the 686 or 682 Titroprocessor

The 686 and 682 Titroprocessors can be used to determine not only the hydroxyl and saponification number, but also the COD.

Fig. 9 shows the interconnections with the 686 or 682 Titroprocessor for determination of the hydroxyl and saponification numbers. For the COD determination with the 686 Titroprocessor, a gold rod and reference electrode have to be attached in place of the pH electrode. Rinsing with a 683 Pump is also shown in Fig. 9 (for rinsing with an air pump and electromagnetic valve, see Fig. 10).

You will find examples of the programming of the 686/682 Titroprocessor in section 5.




Fig. 9: Interconnection with 686/682 Titroprocessor and 665 Dosimat for the determination of hydroxyl and saponification numbers

664 Control Unit:

Mode 4 or mode 5

V2: Siphoning off of the sample vessel contents with 683 Siphoning Pump

V3: Addition of the rinsing solution with 683 Rinsing Pump

686/682 Titroprocessor:

The 3.540.2191 Sample Changer Interface must be inserted in connector B.

3.5.2. Interconnection with the 702 SET/MET Titrino

Interconnection with the 702 SET/MET Titrino is especially suitable for COD determinations (see Fig. 10). Fig. 10 shows rinsing using an air pump and electromagnetic valve (for rinsing with a 683 Pump, see Fig. 9).

You will find an example of the programming of the 702 SET/MET Titrino in section 5.3.4.




Fig. 10: Interconnection with the 702 SET/MET Titrino for the COD determination

- | | |
|----------------------|------------------------------------------------------------------------------|
| 664 Control Unit: | Mode 4 or Mode 5 |
| V2: | Siphoning off of the sample vessel contents with 683 Siphoning Pump |
| V3: | Addition of the rinsing solution with the aid of Air Pump and Magnetic valve |
| 702 SET/MET Titrino: | Mode MET |

3.6. Settings at 664.0020 Control Unit

3.6.1. Program selection

Before the start of an analysis, the desired program must be set on the 664 Control Unit:

- 1► Switch on control unit (see *section 3.4*).
- 2► Loosen cover of the control unit by turning the fastening screw anti-clockwise and then remove. The adjustment controls of the 664 Control Unit are now accessible (see *Fig. 6*).
- 3► Switch in desired program by turning the selector switch **48**. Possible selections are the two programs **Mode 4** and **Mode 5** (see *Table 1*), which are distinguished only in the manner in which the sample vessels are rinsed (with Mode 4, rinsing is performed once, with Mode 5 twice).
- 4► Screw on cover of control unit.

	Mode 4	Mode 5
Function	Titration with 686, 682 or 702	Titration with 686, 682 or 702
Program sequence	1) Shift sample 2) Stirrer on 3) Lower lift 4) Stirrer off 5) Raise lift (Δt_2) 6) Stirrer on 7) Titration 8) Lower lift completely 9) Stirrer off 10) Open valve V2 11) Delay time (Δt_1) 12) Close valve V2 13) Rinsing with valve V3 ($2 \times V3 = \text{max. } 30\text{s}$) 14) Siphoning with valve V2 ($2 \times V2 = \text{max. } 30\text{s}$) 15) Raise lift completely 16) Delay time (Δt_3)	1) Shift sample 2) Stirrer on 3) Lower lift 4) Stirrer off 5) Raise lift (Δt_2) 6) Stirrer on 7) Titration 8) Lower lift completely 9) Stirrer off 10) Open valve V2 11) Delay time (Δt_1) 12) Close valve V2 13) Rinsing with valve V3 ($2 \times V3 = \text{max. } 30\text{s}$) 14) Siphoning with valve V2 ($2 \times V2 = \text{max. } 30\text{s}$) 15) Rinsing with valve V3 ($2 \times V3 = \text{max. } 30\text{s}$) 16) Siphoning with valve V2 ($2 \times V2 = \text{max. } 30\text{s}$) 17) Raise lift completely 18) Delay time (Δt_3)
664 Control Unit settings	V1: 0...15 s not programmed V2: 0...15 s Siphoning during $2 \times V2 = \text{max. } 30\text{s}$ V3: 0...15 s Rinsing during $2 \times V3 = \text{max. } 30\text{s}$ Δt_1 : 0...5 min Delay time (siphoning off the sample) Δt_2 : 0...15 s Raising the lift to the working position Δt_3 : 0...15 s Delay time	V1: 0...15 s not programmed V2: 0...15 s Siphoning during $2 \times V2 = \text{max. } 30\text{s}$ V3: 0...15 s Rinsing during $2 \times V3 = \text{max. } 30\text{s}$ Δt_1 : 0...5 min Delay time (siphoning off the sample) Δt_2 : 0...15 s Raising the lift to the working position Δt_3 : 0...15 s Delay time
Interconnections	686: <i>Fig. 9</i> 682: <i>Fig. 9</i> 702: <i>Fig. 10</i>	686: <i>Fig. 9</i> 682: <i>Fig. 9</i> 702: <i>Fig. 10</i>
Setting the delay and valve opening times	V1: not programmed V2: 3 s V3: 2 s Δt_1 : 15 s Δt_2 : 1-1.5 s Δt_3 : 4-5 s	V1: not programmed V2: 3 s V3: 2 s Δt_1 : 15 s Δt_2 : 1-1.5 s Δt_3 : 4-5 s

Table 1: Overview of the programs of 664.0020 Control Unit

3.6.2. Setting the delay and valve opening times

The delay and valve opening times are set in the factory to the values shown in *Table 1*. These values have been optimized for the two programs *Mode 4* and *Mode 5*, but if you wish to change these times proceed as follows:

- 1▶ Switch on control unit (see *section 3.4*).
- 2▶ Loosen cover of the control unit by turning the fastening screw anti-clockwise and then remove.
The adjustment controls **45** and **47** of the 664 Control Unit are now accessible (see *Fig. 6*).
- 3▶ Set valve opening times V2 and V3 (**45**) by turning the potentiometer screws with screwdriver **49**.

*Turning in a clockwise direction increases the time interval,
turning in an anti-clockwise direction decreases it.*

If the adjustment range of the potentiometers is exceeded, when the time is decreased a jump to infinite time occurs, whereas at the other end of the scale the time interval no longer increases despite further rotation in a clockwise direction.

The maximum time interval for the valve opening times is 15 s. To determine the current setting of the valve opening time, either briefly press the appropriate valve switches V2 (**24**) and V3 (**25**) or briefly plug the contact pin **49** into the appropriate contact socket of the adjustment controls **45**. This lights up the LED for a period corresponding to the set opening time.

- 4▶ Set delay times Δt_1 , Δt_2 , Δt_3 (**47**) by turning the potentiometer screws with screwdriver **49**:

*Turning in a clockwise direction increases the time interval,
turning in an anti-clockwise direction decreases it.*

If the adjustment range of the potentiometers is exceeded, when the time is decreased a jump to infinite time occurs, whereas at the other end of the scale the time interval no longer increases despite further rotation in a clockwise direction.

The maximum time interval for delay times Δt_2 and Δt_3 is 15 s. To ascertain the current valve opening time, the contact pin **49** is plugged briefly into the appropriate contact socket of the adjustment controls **47**. This lights up the LED for a period corresponding to the set opening time.

In the case of delay time Δt_1 , the adjustable time interval is 0...5 min (= 20 × 0...15 s). To ascertain the set time interval, plug in the contact pin and measure the on-period of the flashing LED (contact pin remains plugged in). This on-period is 1/20 of the set interval. In other words, should an interval of, e.g. 3 min be set, the on-period is regulated to $3 \times 60 \text{ s} / 20 = 9 \text{ s}$.

Note: The control key may not be switched to "start" until the manually initiated functions are complete and no LEDs are lit up!

- 5▶ Screw on cover of the control unit.

3.7. Filling the tubing

The operational procedure of the 676 Sample Changer with the 664 Control Unit is described. If the 676 Sample Changer is controlled by the 670 Titroprocessor, the corresponding program commands must be entered as described in the 670 Instructions for Use (see also section 3.5.3).

To fill the 6.1812.000 PTFE Tubing fastened to connecting nipple **12** with rinsing solution, first ensure that the 664 Control Unit is switched on (see section 3.4) and connected to the peripherals (see section 3.5), then proceed as follows:

- 1▶ Pull glass vessel detector (contact bridge in the titration position facing lift tower **7**) by hand in the direction of rotation and hold.
- 2▶ Press key **27** (↑) on 664 Control Unit *up* until the lifting arm **10** is right in the top rest position.
- 3▶ Press key **28** on 664 Control Unit *down*. As soon as the turntable **16** starts to rotate, the glass vessel detector can be released.
- 4▶ While the turntable **16** is rotating, place a sample vessel (6.1821.000 Bark-e/Behr or 6.1821.010 Merck) in the holder for the conditioning vessel permanently mounted on the turntable **16**. The turntable continues to turn until the conditioning vessel is positioned exactly under the titration head **13**.
- 5▶ Press *down* and hold key **27** (↓) on 664 Control Unit until the lifting arm **10** is right in the bottom rest position.
- 6▶ Press *down* and hold switch **25** (V3) on 664 Control Unit until only rinsing agent flows out of the rinse nozzle **14**.

The 6.1805.110 FEP Tubing attached to the burette tip **8** must then be filled with titrant. The procedure is described in the appropriate 665 *Dosimat* or 702 *SET/MET Titrimo Instructions for Use*.

3.8. Setting the stirring speed

The stirring speed of the magnetic stirrer built into the 676 Sample Changer is set at rotary regulator **2** to ca. "3".

If necessary, the stirring speed can be increased or decreased by turning switch **2**.

3.9. Conditioning

The operational procedure of the 676 Sample Changer with the 664 Control Unit is described. If the 676 Sample Changer is controlled by the 670 Titroprocessor, the corresponding program commands must be entered as described in the 670 Instructions for Use (see also section 3.5.3).

When not in use, the electrodes mounted in the titration head must always be immersed in a conditioning solution appropriate to the type of electrodes used.


- Determination of hydroxyl and saponification numbers: $c(\text{KCl}) = 3 \text{ mol/L}$ (6.2308.020 electrolyte solution)
- COD determination: $c(\text{KNO}_3) = \text{sat.}$ (6.2310.000 electrolyte solution)

For positioning of the conditioning solution, proceed as follows:

- 1▶ A sample vessel (6.1821.000 Barkey/Behr or 6.1821.010 Merck) is filled with 100 mL conditioning solution (see above).
- 2▶ Press key **27** (↑) on 664 Control Unit *up* until the lifting arm **10** is right in the top rest position.
- 3▶ Press key **28** on 664 Control Unit *down*. The turntable **16** starts to rotate.
- 4▶ While the turntable **16** rotates, remove the sample vessel used to fill the tubing from the holder.
- 5▶ Place sample vessel with conditioning solution in the permanently mounted holder on the turntable **16**. The turntable continues to rotate until the conditioning vessel is positioned exactly under the titration head **13**.
- 6▶ Press *down* and hold key **27** (↓) on 664 Control Unit until the lifting arm **10** is right in the bottom rest position.

4. Operation

4.1. Operation scheme



4.2. Manual operation

The operational procedure with the 676 Sample Changer and the 664 Control Unit is described. If the 676 Sample Changer is controlled by the 670 Titroprocessor, the appropriate program commands must be entered as described in the 670 Instructions for Use (see also section 3.5.3).

4.2.1. Manual operation rinsing

The rinsing operation is initiated by pressing key **25** (V3) down. Depending on the type of rinsing selected (see section 3.3), either the 683 Rinsing Pump is switched on or the 6.2726.070 Electro-magnetic Valve opened. The rinsing time is given by the set valve opening time V3 (see section 3.6).

4.2.2. Manual operation siphoning

The siphoning operation is initiated by pressing key **24** (V2) down. This switches on the 683 Siphoning Pump. The siphoning time is given by the set valve opening time V2 (see section 3.6).

4.2.3. Manual operation of the magnetic stirrer

By pressing the stirrer key **26** up, the magnetic stirrer built into the 676 Sample Changer is switched on.

The stirring speed (200 ... 1900/min) can be set at rotary regulator **2**.

To switch the stirrer off, the stirrer key **26** is switched to the middle position.

4.2.4. Manual operation of the lifting arm

By pressing key **27** (↑) up, the lifting arm **10** is raised.

By pressing key **27** (↓) down, the lifting arm **10** is lowered.

Note: *If there is no vessel in the titration position, for manual operation of the lifting arm the glass vessel detector (contact bridge at the titration position facing lift tower 7) must be moved by hand in the direction of rotation and held.*

*The lifting arm **10** can be moved with key **27** only when the turntable **16** is in the correct position (the glass vessel or the glass vessel holder must be directly under the titration head). It may thus be necessary to first move the turntable by hand into the correct position.*

4.2.5. Manual operation of the sample transport

By pressing key **28** down, movement of the turntable **16** is started. It rotates until the next sample vessel has reached the titration position.

Note: *If there is no vessel at the titration position, for manual operation of the sample transport lifting arm the glass vessel detector (contact bridge at the titration position facing lift tower 7) must be moved by hand in the direction of rotation and held.*

*The turntable **16** can be moved with key **28** only when the lifting arm **10** is in the top rest position.*

4.3. Notes on operation of the 676 Sample Changer

The operational procedure with the 676 Sample Changer and the 664 Control Unit is described. If the 676 Sample Changer is controlled by the 670 Titroprocessor, the appropriate program commands must be entered as described in the 670 Instructions for Use (see also section 3.5.3).

4.3.1. Sample rack

The sample rack **4** can hold 15 sample vessels of the Barkey/Behr type (6.1452.000) or Merck type (6.1452.010). The two holding brackets **5** are used to transport the sample rack (see Fig. 2).

The sample rack together with samples vessels can be directly inserted in the heating and cooling block of the digestion apparatus *csb-therm X* of the *Barkey Labor Technik* company. With other digestion apparatus, the sample vessels must be positioned individually in the sample holder **4** after digestion/heating.

4.3.2. Stirring bars

Before the start of the analysis program, a stirring bar must be added to every sample vessel:

- for Barkey/Behr sample vessels: 6.1906.010 Stirring Bars (triangular cross-section)
- for Merck sample vessels: 6.1906.020 Stirring Bars (oval)

4.3.3. Check on the electrolyte vessel

In COD determinations, before every start of the analysis program a check must be made to see if there is sufficient electrolyte solution in the 6.1240.0X0 Electrolyte Vessel. If need be, this vessel must be topped up with electrolyte solution (see section 3.2.2).

4.3.4. Start of the analysis program

Before the start of the analysis program, the sample counter **31** must be set to "00".

The set analysis program (see section 3.6.1) is started by pressing the control key **32** down.

The analysis program is then performed for all samples on the sample rack. After the last sample, the titration head is lowered into the conditioning beaker.

4.3.5. Manual stop of the analysis program

By pressing the control key **32** up, the current function is immediately terminated. An exception is the ongoing titration (incl. stirrer), which is allowed to run to completion before being stopped.

4.3.6. Conditioning solution

When not in use, the electrodes mounted in the titration head must always be immersed in a conditioning solution appropriate to the type of electrodes used.

- Determination of hydroxyl and saponification numbers: $c(\text{KCl}) = 3 \text{ mol/L}$
(6.2308.020 electrolyte solution)
- COD determination: $c(\text{KNO}_3) = \text{sat.}$
(6.2310.000 electrolyte solution)

On completion of the analysis program (that is, after max. 15 samples), the titration head together with the electrodes is automatically lowered into the vessel with conditioning solution (for first-time positioning of the conditioning solution, see section 3.9).

4.3.7. Siphoning off with 683 Pump

Note: The waste container must be **emptied** from time to time.

No solids can be siphoned off with the 683 Pump. These could damage the pump.

5. Application examples

5.1. Determination of the hydroxyl number

5.1.1. Analysis procedure

The analysis procedure depends on the nature of the samples.

Detailed instructions are described in *METROHM Application Bulletin 141* (for edible oils and fats) and *METROHM Application Bulletin 200* (for plastics).

5.1.2. Titration with the 686 or 682 Titroprocessor

The 686 and 682 Titroprocessors allow the automatic performance and evaluation of titrations. For operation and programming of the instruments, see the appropriate Instructions for Use.

Example: Determination of the hydroxyl number in edible oils and edible fats

- Sample preparation: see Appl. Bull. 141
- Titrant: $c(\text{KOH}) = 0.5 \text{ mol/L}$
- Calculation:

$$\text{Hydroxyl number (mg/g KOH)} = \frac{28.05 \cdot (B - A)}{E}$$

- A = mL KOH for the sample
- B = mL KOH for the blank value
- E = sample weight in g


- Instrument parameters, formulae and curve:

```
GET U ***** # 3
electr. input 1
titr.rate 1.00 ml/min
anticip. 5
stop V 24.00 mV
start V 15.00 mV
temp. 20.0 °C
EP crit. 8
=====
```

```
date 89-06-19 time 15:15
GET U ***** # 3
Id.#1 2305
Id.#2 1.1
C00= 7.908 g
U(init) 92 mV
V/ml U/mV
EP1 21.856 -152
RS1 1.631
stop V reached
=====
```

```
GET U ***** # 3
F1=(C01-EP1)*C02/C00+C03
;3;
C01= 22.062
C02= 28.05
C34= .90
=====
```

```
date 89-06-19 time 15:15
GET U ***** # 3
2.0ml/div ΔU=100mV/div
start V 15.000 mV
```



5.2. Determination of the saponification number

5.2.1. Analysis procedure

The analysis procedure depends on the nature of the samples.

Detailed instructions for edible oils and edible fats are described in *METROHM Application Bulletin 141*.

5.2.2. Titration with the 686 or 682 Titroprocessor

The 686 and 682 Titroprocessors allow the automatic performance and evaluation of titrations. For operation and programming of the instruments, see the appropriate Instructions for Use.

Example: Determination of the hydroxyl number in edible oils and edible fats

- Sample preparation: see *Appl. Bull. 141*
- Titrant: $c(\text{HCl}) = 0.5 \text{ mol/L}$
- Calculation:

$$\text{Saponification number (mg/g KOH)} = \frac{28.05 \cdot (B - A)}{E}$$

A = mL HCl for the sample
 B = mL HCl for the blank value
 E = sample weight in g


- Instrument parameters, formulae and curve:

```
GET U ***** # 0
electr. input 1
titr.rate .50 ml/min
anticip. 5
stop V 14.00 ml
temp. 20.0 °C
EP crit. 8
=====
```

```
date 89-06-15 time 07:36
GET U ***** # 0
Id.#1 2305
Id.#2 1.1
C00= 2.013 g
U(init) -450 mV
V/ml U/mV
EP1 13.132 -199
RS1 292.53
=====
```

```
GET U ***** # 0
F1=(C01-EP1)*C02/C00;2;
C01= 34.125
C02= 28.05
=====
```

```
date 89-06-15 time 07:36
GET U ***** # 0
2.0ml/div ΔU=100mV/div
start V 0.000 ml
```



5.3. COD determination

5.3.1. Analysis procedure

The analysis procedure for COD determinations is specified in *DIN 38409 Part 41*.

Detailed instructions can be found in *METROHM Application Bulletin 212*.

5.3.2. COD determination with the 686 Titroprocessor

The 686 and 682 Titroprocessors allow the automatic performance and evaluation of titrations. For operation and programming of the instruments, see the appropriate Instructions for Use.

- Sample preparation: see *DIN 38409 Part 41* and *Appl. Bull. 212*
- Titrant: $(\text{NH}_4)_2\text{Fe}(\text{SO}_4)_2$
- Instrument parameters

```
MET U ***** # 16
electr. input 1
vol.step      .08 ml
drift         70 mV/min
stop V        25.00 ml
stop EP#      1
temp.         25.0 °C
EP crit.U     40 mV
=====
```

```
date 87-12-09 time 16:04
MET U ***** # 16
C00=          10 ml
U(init)       -381 mV
              V/ml      U/mV
EP1           2.186      -270
RS1           616      ppm
stop EP# reached
=====
```

```
MET U ***** # 0
F1=(C30-EP1)*C01*C02/C00
;0;ppm
=====
```

- Calculation: COD value (mg/L = ppm) = $(C30-EP1) \cdot C01 \cdot C02 / C00$

- C01: Concentration of the titrant in mol/L
- C02: Equivalence factor $f = 8000$
- C00: Volume in mL of the water sample used in the investigation
- C30: Volume in mL of the titrant needed for the blank sample (twice dist. water)
- EP1: Volume in mL of the titrant needed for the analysis sample (calculated automatically)

5.3.3. COD determination with the 670 Titroprocessor

The 670 Titroprocessor allows the automatic performance and evaluation of titrations as well as control of the 676 Sample Changer. For the operation and programming of the instrument, refer to the appropriate Instructions for Use.

- Sample preparation: see *DIN 38409 Part 41* and *Appl. Bull. 212*
- Titrant: $(\text{NH}_4)_2\text{Fe}(\text{SO}_4)_2$

• Program page:

COD Determination with 676 Sample Changer

METHOD PAGE 3A
OPERATION SEQUENCE

OPERATIONS/PARAMETERS		OPERATIONS/PARAMETERS	
1	OMOVE1;	9c	M.VALUE OFF
2	STIR 1;	10	WPOS 1;
3	WPOS 1;	11	OSTIR1;
4	OSTIR1;	12	CTR 11; 20.0 s
5	RDIR 1;	13	OCTR11; 100 ms
5a	t.MOVE 1.0 s	14	CTR 12; 6.0 s
6	STIR 1; 3 s	15	OCTR12; 100 ms
7	MEAS 1;	16	CTR 11; 7.0 s
7a	QUANTITY U	17	OCTR11; 100 ms
7b	DRIFT /min 35 mV	18	CTR 12; 5.0 s
7c	M.DELAY 22 s	19	OCTR12; 100 ms
8	MONT 1;	20	CTR 11; 10.0 s
8a	V.MONOT 80 uL	21	OCTR11; 100 ms
8b	DOS.RATE MAX.	22	RPOS 1; 5 s
9	TSTOP ;	23	SHIFT1;
9a	N.EPs 1	24	CMOVE1;
9b	VOLUME 12.000 mL	25	BEEP ;END ;

• Calculations:

COD Determination with 676 Sample Changer

METHOD PAGE 6
CALC.STATEMENTS

a	b	c	d
RES.QUANTITY	R#, I#, Cx FORMULA	DPL	RES.UNIT
1 COD	R1=(CA-E1)*D1*C1/S0	0	mg/l

COD Determination with 676 Sample Changer

METHOD PAGE 7
ASSIGNMENTS

a	b	a	b
C# CALC.VAR	COMMENT	Cx COM.VAR	COMMENT
1 C1 8000.00	aequiv. factor	21 CA 9.54100 E- 3	Blank value

• Results:

COD Determination with 676 Sample Changer

METHOD
S0 SIZE(m,V) 20.0000 mL

TITR	EP V/mL	EP M.VAL	TEST M.VAL	FIN.V/mL	pK or HNP
#1 U					
E1	3.035	478.3 mV		3.280	

CALCULATIONS

COD = 312 mg/l

DATE 91-04-16

TIME 12:33

RUN # 5

• Curve example:

COD Determination

METHOD


V.BEGIN= 0 mL

V.DIV= 1.000 mL

U/mV

TITR # 1

300 360 420 480 540 600 660 720 780 840 900



FIN.VOLUME= 3.280 mL

DATE 91-04-16

TIME 12:33

RUN # 5

0022

5.3.4. COD determination with the 702 SET/MET Titrino

The 702 SET/MET Titrino allows the automatic performance and evaluation of titrations. For the operation and programming of the instrument, refer to the appropriate Instructions for Use.

- Sample preparation: see *DIN 38409 Part 41* and *Appl. Bull. 212*
- Titrant: $(\text{NH}_4)_2\text{Fe}(\text{SO}_4)_2$
- Instrument parameters

```

date 91-04-02      time 14:43      1
MET U      *****
parameters
>titration parameters
  V step          0.08 ml
  titr.rate       max. ml/min
  signal drift    70 mV/min
  equilibr.time   22 s
  start V:        OFF
  pause          0 s
  meas.input:     1
  temperature     25.0 °C
>stop conditions
  stop V:         abs.
  stop V          25 ml
  stop U          OFF mV
  stop EP         1
  filling rate    max. ml/min
>statistics
  status:         OFF
>evaluation
  EPC             40 mV
  EP recognition: all
  fix EP1 at U    OFF mV
>preselections
  req.ident:      OFF
  req,smpl size:  OFF
  
```

```

date 91-04-02      time 14:43      1
U(init)          21 mV MET U      *****
smpl size        10 ml
id#2             5
EP1              2.493 ml          427 mV
RS1              199 ppm
=====
  
```

```

date 91-04-02      time 14:43      1
MET U      *****
>calculations
RS1=(C30-EP1)*C01*C02/C00;0;ppm
  
```

- Calculation: COD value (mg/L = ppm) = $(C30-EP1) \cdot C01 \cdot C02 / C00$

C01: Concentration of the titrant in mol/L

C02: Equivalence factor $f = 8000$

C00: Volume in mL of the water sample used in the investigation

C30: Volume in mL of the titrant needed for the blank sample (twice dist. water)

EP1: Volume in mL of the titrant needed for the analysis sample (calculated automatically)

6. Technical specifications

6.1. 676 Sample Changer

- **Capacity** max. 15 sample vessels (+ 1 vessel for conditioning)

- **Drive mechanism**
 - 1 DC motor for raising and lowering the movable lifting arm with titration head
 - 1 DC motor for rotating the turntable
 - 1 DC motor for electromagnetic stirrer

- **Power supply** from 664.0020 Control Unit

- **Materials**

<i>Housing</i>	light alloy, multicoat stove-enamel
<i>Glide surface</i>	polyethylene terephthalate (PETP)
<i>Turntable</i>	light alloy, polyvinyl chloride (PVC)
<i>Lift tower</i>	polyvinyl chloride (PVC)
<i>Lifting arm</i>	light alloy, multicoat stove-enamel
<i>Titration head</i>	polytetrafluorethylene (PTFE)
<i>Sample rack</i>	light alloy, clear anodised

- **Dimensions**

<i>Width</i>	335 mm
<i>Height</i>	615 mm
<i>Depth</i>	490 mm

- **Weight** (incl. standard acc.) 9.5 kg

- **Sample vessels**

	<u>6.1452.000</u>	<u>6.1452.010</u>
<i>Application</i>	for digestion block Barkey or Behr	for digestion block Merck
<i>Material</i>	borosilicate glass 3.3	borosilicate glass 3.3
<i>Final volume</i>	170 mL	230 mL
<i>Diameter</i>	40.8 mm	40.8 mm
<i>Height</i>	200 mm	260 mm
<i>Weight</i>	112 g	90 g

6.2. 664.0020 Control Unit

- **Housing**
injection-moulded light alloy with several coats of stove enamel; top cover removable for easy access to valve and time settings
- **Number of samples**
15 (sample counter must be set to "00")
- **Analysis programs**
 - Number*
2 standard programs (Mode 4 + Mode 5), selected on thumbwheel switch under top cover
 - Program steps*
total 1024 (128 per program)
 - Program memory*
"Erasable Programmable Read Only Memory" (EPROM)
- **Functions**
 - Manual*
open valves 1...3
stir
raise and lower lifting arm
transport samples
start and stop program
 - Automatic*
the samples inserted into the sample rack are processed according to the program chosen
- **Time settings**
 - Delay times*
1 × 0...5 min ; 2 × 0...15 s
 - Valve opening times*
3 × 0...15 s
- **Ambient temperature**
 - Nominal operating range*
+5 ... +40 °C
 - Storage, transport*
-10 ... +70 °C
- **Mains connection**
 - Mains voltage*
 $U_{\sim} = 100, 117, 220, 240 \text{ V} \pm 10 \%$, switchable
 - Mains frequency*
 $f = 50...60 \text{ Hz}$
 - Power consumption*
ca. 50 VA, depending on extent of expansions
 - Fuse*
 $\varnothing = 5 \text{ mm}$, length $L = 20 \text{ mm}$
100 and 117 V: 0.5 A (slow-blow)
220 and 240 V: 0.25 A (slow-blow)
- **Safety specifications**
in accordance with IEC Publication 348 (class 1 instrument)
- **Dimensions**
 - Width*
350 mm
 - Height*
130 mm
 - Depth*
280 mm
- **Weight (incl. standard acc.)**
5 kg

6.3. 6.2726.070 Electromagnetic Valve



- **Application** rinsing valve for aqueous solutions
(not suitable for acetone or similar solvents)
- **Material** polytetrafluorethylene (PTFE)
polyvinylidene fluoride (PVDF)
- **Bore** 3 mm
- **Metering error** ca. 5 %
- **Temperature range** 0... + 100°C
- **Power supply** 24 V AC (from 664.0020 Control Unit)
- **Power consumption** 8 VA
- **Dimensions**
 - Width* 37 mm
 - Height* 95 mm
 - Depth* 106 mm
- **Weight** 250 g

6.4. 6.2817.000 Air Pump

- **Function principle** membrane
- **Pump capacity**
 - Delivery rate* ca. 2 ... 6 L/min, adjustable
 - Pressure head* ca. 1 ... 5 m water column (0.1 ... 0.5 bar), adjustable
- **Materials**
 - Membrane* ethylene polypropylene elastomer (EPDM)
 - Housing* aluminium
 - Tube nipple* aluminium
 - Inner tube* polyvinylchloride (PVC)
- **Ambient temperature**
 - Nominal operating range* +5 ... +40 °C
 - Storage, transport* -10 ... +70 °C
- **Mains connection**
 - Mains voltage* $U_{\sim} = 220 \dots 240 \text{ V}$
(other voltages on request)
 - Mains frequency* $f = 50 \dots 60 \text{ Hz}$
 - Power consumption* 13 VA
- **Dimensions**
 - Width* 190 mm
 - Height* 110 mm
 - Depth* 120 mm
- **Weight** 2.2 kg







6.5. Control inputs and outputs

Function		3.540.1640 (E)	extern
Start/Stop 1 titrate, dispense			<p>0.5 A = max. 72 V = max.</p>
Start/Stop 2 titrate, dispense			<p>1) as pulse: 650ms > t_{on} > 160ms</p> <p>2) as continuous pulse: (e.g. until stop by upper limit switch)</p>
Start/Stop 3 titrate, dispense	3.540.1640 (F)		
Program advance (command)			<p>$i_d \geq 3.5 \text{ mA}$ pulse length: > 20 ms</p>
Voltages			

Function	3.540.1620 (J)	extern
<p>Valves 1...3</p> <p>output</p>		<p>$U = 24\text{ V} \sim$</p> <p>maximum load of the 24 V \sim supply voltage:</p> <p>$I_{\text{max}} = 1.5\text{ A}$</p>
 Metrohm	664.0020 Control Unit	3.540.1640/664 4 E 4

Function	3.540.2111 (B)	extern
Top end switch at electrode lift L = stop raising		is closed when upper position is reached 0V
Bottom end switch at electrode lift L = stop lowering		is closed when lower position is reached 0V
Glass vessel detector and position switch at turntable L = turntable in position and start position or glass vessel detector L = turntable in position and start position		turntable position e 10 glass vessel detector sample counter e 1 start position position 0V
Motor for turntable		10 V I < 270 mA
Voltages	5V — e 3 0V — e 4/5	
	664.0020 Control Unit	3.540.2111/664 4 E 5

Function	3.540.2111 (B)	extern
Motor for lifting arm		
Non-occupied pins	e6, e7, e9, e16, e17, e18, e22, e25 eB1 ... eB26	
Contact location at the plug ordering number: 3.540.8440		
Any liability for damages due to wrong interconnections between instrumentation is refused.		
	664.0020 Control Unit	3.540.2111/664 4 E 6

Function	3.540.2420 (C)	extern																														
Start program		start (pulse)																														
Stop program		stop (pulse)																														
Program selection with Mode 8 : <table border="1" data-bbox="113 972 448 1099"> <thead> <tr> <th>flag 0</th><th>flag 1</th><th>program segment</th></tr> </thead> <tbody> <tr> <td>L</td><td>L</td><td>800-87F (Mode 0)</td></tr> <tr> <td>L</td><td>H</td><td>900-97F (Mode 2)</td></tr> <tr> <td>H</td><td>L</td><td>A00-A7F (Mode 4)</td></tr> <tr> <td>H</td><td>H</td><td>B00-B7F (Mode 6)</td></tr> </tbody> </table> with Mode 9: <table border="1" data-bbox="113 1151 448 1279"> <thead> <tr> <th>flag 0</th><th>flag 1</th><th>program segment</th></tr> </thead> <tbody> <tr> <td>L</td><td>L</td><td>880-8FF (Mode 1)</td></tr> <tr> <td>L</td><td>H</td><td>980-9FF (Mode 3)</td></tr> <tr> <td>H</td><td>L</td><td>A80-AFF (Mode 5)</td></tr> <tr> <td>H</td><td>H</td><td>B80-BFF (Mode 7)</td></tr> </tbody> </table>	flag 0	flag 1	program segment	L	L	800-87F (Mode 0)	L	H	900-97F (Mode 2)	H	L	A00-A7F (Mode 4)	H	H	B00-B7F (Mode 6)	flag 0	flag 1	program segment	L	L	880-8FF (Mode 1)	L	H	980-9FF (Mode 3)	H	L	A80-AFF (Mode 5)	H	H	B80-BFF (Mode 7)		flag 1 flag 0
flag 0	flag 1	program segment																														
L	L	800-87F (Mode 0)																														
L	H	900-97F (Mode 2)																														
H	L	A00-A7F (Mode 4)																														
H	H	B00-B7F (Mode 6)																														
flag 0	flag 1	program segment																														
L	L	880-8FF (Mode 1)																														
L	H	980-9FF (Mode 3)																														
H	L	A80-AFF (Mode 5)																														
H	H	B80-BFF (Mode 7)																														
Busy		busy = L																														
Voltages																																
	664.0020 Control Unit	3.540.2420/664 4 E 7																														

7. Scope of delivery and ordering designations

7.1. 676 Sample Changer

676.0300 Sample Changer **2.676.0300**

including the following accessories:

1 × Sample rack for 15 reaction vessels from Behr, Barkey or Merck incl. 2 holding brackets	6.2041.030
1 × Tubing connection of FEP (tetrafluoroethylene hexafluoropropylene copolymer) with 2 screw nipples with M6 thread, protected against light and buckling i.d. = 1.5 mm, length L = 80 cm	6.1805.110
1 × Rinse nozzle with M6 thread	6.2740.010
1 × PTFE/PVDF Electromagnetic valve (mounted at the instrument)	6.2726.070
2 × Connecting cable 664 Control Unit – 683 Pump, length L = 1 m	6.2121.010
1 × Instructions for Use (english)	8.676.1003

7.2. 664.0020 Control Unit

664.0020 Control Unit **2.664.0020**

for 676 Sample Changer,

with 3 adjustable valve opening times and 3 delay times,

including the following accessories:

1 × Plastic dust cover	6.2723.140
1 × Mains cable (to customer specification)	
<u>cable socket</u> <u>cable plug</u>	
type IEC 320/C 13 type SEV 12 (Suisse ...)	6.2122.020
type IEC 320/C 13 type CEE (7), VII (Germany ...)	6.2122.040
type CEE (22), V type NEMA 5-15 (USA ...)	6.2122.070

7.3. Options

7.3.1. Reaction vessels

Reaction vessel appropriate for Barkey or Behr digestion block	6.1452.000
Reaction vessel appropriate for Merck digestion block	6.1452.010

7.3.2. Accessories for the determination of hydroxyl and saponification numbers

For reaction vessels from Barkey or Behr

1 × Combined pH microelectrode with plug-in head; lower, tapered part of length L = 220 mm, Ø 6 mm	6.9902.048
1 × Cable with F plug, length L = 2 m, for 6.9902.048 pH electrode (for use with the 670 Titroprocessor, the 6.2104.060 Cable is needed instead)	6.2104.030
1 × 250 mL electrolyte solution c(KCl) = 3 mol/L, for 6.9902.048 pH electrode	6.2308.020
1 × Burette tip, length L = 280 mm	6.1543.130
1 × Siphoning tube, length L = 178 mm	6.1821.000
15 × Stirrer bars, length L = 25 mm, triangular cross-section	6.1906.010
1 × Spiral band for holding cables and tubing together Ø = 7 mm, length L = 50 cm	6.1815.010

Subject to modifications!

For reaction vessels from Merck

1 × Combined pH microelectrode with plug-in head; lower, tapered part of length $L = 220$ mm, \varnothing 6 mm	6.9902.048
1 × Cable with F plug, length $L = 2$ m, for 6.9902.048 pH electrode (for use with the 670 Titroprocessor, the 6.2104.060 Cable is needed instead)	6.2104.030
1 × 250 mL Electrolyte solution $c(\text{KCl}) = 3$ mol/L, for 6.9902.048 pH electrode	6.2308.020
1 × Burette tip, length $L = 312$ mm	6.1543.140
1 × Siphoning tube, length $L = 238$ mm	6.1821.010
15 × Stirrer bars, oval, 13 mm × 26 mm	6.1906.020
1 × Spiral band for holding cables and tubing together, $\varnothing = 7$ mm, length $L = 50$ cm	6.1815.010
1 × Set of spacers (2 pieces) for conversion of the standard version of the 676 Sample Changer to Merck reaction vessels (conversion in the factory)	6.2049.000

7.3.3. Accessories for the COD determination

For reaction vessels from Barkey or Behr

1 × Reference system Ag/AgCl	6.0724.140
1 × Cable with 2 B plugs, length $L = 2$ m; for reference system	6.2106.060
1 × Electrolyte vessel	6.1240.050
1 × 250 mL Electrolyte solution KNO_3 saturated, for reference system	6.2310.000
1 × Gold rod electrode (for use with the 670 Titroprocessor, an additional 6.2103.080 Adaptor is needed)	6.1248.060
1 × Cable for gold rod electrode; with F plug	6.2120.030
1 × Burette tip, length $L = 280$ mm	6.1543.130
1 × Siphoning tube, length $L = 178$ mm	6.1821.000
15 × Stirrer bars, length $L = 25$ mm, triangular cross-section	6.1906.010
1 × Spiral band for holding cables and tubing together, $\varnothing = 7$ mm, length $L = 50$ cm	6.1815.010

For reaction vessels from Merck

1 × Reference system Ag/AgCl	6.0724.140
1 × Cable with 2 B plugs, length $L = 2$ m; for reference system	6.2106.060
1 × Electrolyte vessel	6.1240.060
1 × 250 mL electrolyte solution KNO_3 saturated, for reference system	6.2310.000
1 × Gold rod electrode (for use with the 670 Titroprocessor, an additional 6.2103.080 Adaptor is needed)	6.1248.070
1 × Cable for gold rod electrode; with F plug	6.2120.030
1 × Burette tip, length $L = 312$ mm	6.1543.140
1 × Siphoning tube, length $L = 238$ mm	6.1821.010
15 × Stirrer bars, oval, 13 mm × 26 mm	6.1906.020
1 × Spiral band for holding cables and tubing together, $\varnothing = 7$ mm, length $L = 50$ cm	6.1815.010
1 × Set of spacers (2 pieces) for conversion of the standard version of the 676 Sample Changer to Merck reaction vessels (conversion in the factory)	6.2049.000

7.3.4. Accessories for 676 Sample Changer

Holding brackets for 6.2041.030 Sample rack	4.676.0700
PTFE Tubing, $\varnothing = 4/6$ mm, length $L = 400$ cm	6.1812.000
PTFE/PVDF Electromagnetic valve (<u>not suited</u> for acetone and similar solvents)	6.2726.070
PVDF Reducing valve for 6.2726.070 Electromagnetic valve	6.2726.040
PTFE/PVDF Electromagnetic valve (with reducing valve, <u>not suited</u> for acetone and similar solvents)	6.2726.080
PTFE Electromagnetic valve (specially suited for acetone and similar solvents)	6.2726.000

7.3.5. Accessories for 664.0020 Control Unit

6.2812.000	Contact pin for 664 Control Unit
3.980.2930	Connecting cable 664 Control Unit – 682/686/678
3.980.2990	Connecting cable 664 Control Unit – 682/686/678 + 586
3.980.3030	Connecting cable 664 Control Unit – 682/686/678 + 671
3.980.3560	Connecting cable 664 Control Unit – 702

7.3.6. Auxiliary Instruments

2.686.0100	686 Titroprocessor
	Potentiometric titrator for titrations with automatic endpoint recognition
	as well as titrations to a preset endpoint
3.540.2191	Sample changer interface for 686 Titroprocessor
3.980.3090	Connecting cable 686 Titroprocessor – 665 Dosimat
3.980.2930	Connecting cable 686 Titroprocessor – 664 Control Unit
2.665.0010	665 Dosimat
	for the precise dispensing of reagents
6.3012.213	Exchange Unit with 1 L reagent bottle of brown glass, micro outlet valve,
6.3012.223	burette volume 10 mL
6.3012.223	Exchange Unit with 1 L reagent bottle of brown glass, micro outlet valve,
	burette volume 20 mL
2.702.0010	702 SET/MET Titrimo
	Compact potentiometric Titrator for endpoint and monotonic titrations
3.980.3560	Connecting cable 702 SET/MET Titrimo – 664 Control Unit
6.3012.213	Exchange Unit with 1 L reagent bottle of brown glass, micro outlet valve,
6.3012.213	burette volume 10 mL
6.3012.223	Exchange Unit with 1 L reagent bottle of brown glass, micro outlet valve,
	burette volume 20 mL
2.670.0020	670 Titroprocessor
	Potentiometric titrator for dynamic, monotonic, endpoint and pH stat titrations
3.540.2360	Sample changer interface for 670 Titroprocessor
6.2124.010	Connecting cable 670 Titroprocessor – 665 Dosimat
3.980.3230	Connecting cable 670 Titroprocessor – 677 Drive Unit
2.677.0020	677 Drive Unit
	for control of 676 Sample Changer, stirrer, 683 Pumps, valves
	683 Diaphragm Pump
	Liquid pump with PTFE diaphragm, PTFE-coated valve plates, sealing rings
	and housing of PVDF (polyvinylidene fluoride);
2.683.0024	for 220 V supply
2.683.0025	for 110 V supply
2.683.0026	for 240 V supply
	For siphoning off the titrated solution and rinsing of the reaction vessel and internal fittings, two
	683 Diaphragm Pumps are needed.
	615 Power Distributor
	with built-in interference suppression filter
2.615.0010	with 6 sockets according to SEV (Switzerland)
2.615.0020	with 5 sockets according to VDE (Germany)

Digestion block, cooling block

These are available from the manufacturers of the digestion units, i.e. from Barkey, Behr or Merck.

8. Warranty

The warranty regarding our products is limited to rectification free of charge in our workshops of defects that can be proved to be due to material, design or manufacturing faults which appear within 12 months from the day of delivery. Transport costs are chargeable to the orderer.

For day and night operation, the warranty is valid for 6 months.

Glass breakage in the case of electrodes or other glass parts is not covered by the warranty. Checks which are not a result of material or manufacturing faults are also charged during the warranty period. For parts of outside manufacture insofar as these constitute an appreciable part of our instrument, the warranty stipulations of the manufacturer in question apply.

With regard to the guarantee of accuracy, the technical specifications in the Instructions for Use are authoritative.

concerning defects in material, construction or design as well as the absence of guaranteed features, the orderer has no rights or claims except those mentioned above.

If damage of the packaging is evident on receipt of a consignment or if the goods show signs of transport damage after unpacking, the carrier must be informed immediately and a written damage report demanded. Lack of an official damage report releases METROHM from any liability to pay compensation.

If any instruments and parts have to be returned, the original packaging should be used if at all possible. This applies above all to instruments, electrodes, burette cylinders and PTFE pistons. Before embedment in wood shavings or similar material, the parts must be packed in a dustproof package (for instruments, use of a plastic bag is imperative). If open assemblies are enclosed in the scope of delivery that are sensitive to electromagnetic voltages (e.g. data interfaces etc.) these must be returned in the associated original protective packaging (e.g. conductive protective bag). (Exception: assemblies with built-in voltage source belong in a non-conductive protective packaging). For damage which arises as a result of non-compliance with these instructions, no warranty responsibility whatsoever will be accepted by METROHM.

9. Index

A

Accessories	33-35
Adaptor (6.2103.080)	13,34
Adjustment controls 45	6,15
Adjustment controls 47	6,15
Air Pump (6.2817.000)	
Installation	9
Technical specifications	27
Analysis procedures	20-24
Analysis program	
Selector switch 48	6
Start	19
Stop	19
Application examples	
COD determination	22-24
Determination of hydroxyl number	20
Determination of saponification number	21
Auxiliary instruments	35

B

Burette tip	
Figure	3
Ordering designations	33-34
Burette tip 8	3,7-8,16

C

Cable (3.980.2930)	11-12,35
Cable (3.980.2990)	35
Cable (3.980.3030)	35
Cable (3.980.3090)	11,35
Cable (3.980.3230)	13,35
Cable (3.980.3560)	12,35
Cable (6.2104.030)	11,33-34
Cable (6.2104.060)	13,33-34
Cable (6.2106.060)	8,12,34
Cable (6.2120.030)	12,34
Cable (6.2121.010)	11-12,33
Cable (6.2124.010)	13,35
Capacity	25
Check	7
COD determination	
Accessories	34
Analysis procedure	22
Examples	22-24
with 702 SET/MET Titrimo	24
with 670 Titroprocessor	22-23
with 686 Titroprocessor	22
Conditioning	16
Conditioning solution	16,19
Connecting cable 17	3
Connecting cable 22	3
Connection 9	3
Connection 12	3,9,16
Connector 33	5
Contact pin 49	6,15,35
Control elements	2-6
Control inputs and outputs	28-32

Control key 32	5,15,19
Control Unit (664)	
Accessories	35
Control elements	4-6
Control inputs and outputs	28-32
Figures	4,6
Front	4
Mains connection	10
Ordering designations	33
Program selection	14
Program sequence	14
Rear	4
Scope of delivery	33
Settings	14
Setting the delay times	14-15
Setting the valve opening times	14-15
Switching on	10
Technical specifications	26
Top view	6
Cooling block	35

D

Data 43	5,10
Delay times	14-15,26
Designation 46	6
Determination of hydroxyl number	
Accessories	33
Analysis procedure	20
Example	20
with 686/682 Titroprocessor	20
Determination of saponification number	
Accessories	33
Analysis procedure	21
Example	21
with 686/682 Titroprocessor	21
Diaphragm pump → Pump (683)	
Digestion apparatus	19
Digestion block	35
Dimensions	25-26
Display 35 of set mains voltage	5
Display 36	10
Dosimat (665)	
Electrical interconnections	11,13
Ordering designations	35
Drive Unit (677)	
Electrical interconnections	13
General description	1
Ordering designations	35

E

Earth socket 41	5,10
Electrical interconnections	
with 702 SET/MET Titrimo	12
with 670 Titroprocessor	13
with 686/682 Titroprocessor	11
Electrode cable (6.2104.030)	7,11,33-34
Electrode cable (6.2104.060)	7,8,12,34
Electrode cable (6.2120.030)	8,12,34
Electrolyte solution (6.2308.020)	16,19,33-34
Electrolyte solution (6.2310.000)	8,16,19,34
Electrolyte vessel (6.1240.0X0)	8,19,34
Electromagnetic stirrer → Stirrer	

Electromagnetic valve 18	3,5,9
Electromagnetic valve (6.2726.000)	34
Electromagnetic valve (6.2726.070)	
Connecting cable 17	3
Figure	3
Inlet 20	3
Ordering designation	34
Outlet 19	3
Technical specifications	27
Electromagnetic valve (6.2726.080)	34
Equipping the titration head	7-8
Examples → Application examples	

F

Fastening the 6.1812.000	
PTFE Tubing	10
Fastening tubing and cables	10
FEP Tubing (6.1805.110)	3,7-8,16,33
Filling the tubing	16
Fuse	
Replacement	10
Technical specifications	26
Fuse 44	5,10

G

Glass vessel detector	16,18
Glide surface 3	3,25
Gold rod electrode (6.1248.0X0)	8,12,13,34

H

Holding bracket 5	3,19,33-34
Hydroxyl number → Determination of hydroxyl number	

I

Indicator lamp 1	3,5
Inlet 20	3,9
Installation	7-16
Instrument designation 37	5
Instrument number	5
Interface (3.540.2191)	11,35
Interface (3.540.2360)	13,35
Introduction	1

K

Key 24	5,15,18
Key 25	5,15,18
Key 26	5,18
Key 27	5,16,18
Key 28	5,16,18

L

Lifting arm	
Key 27	5
Manual operation	18
Lifting arm 10	3,5,7-8,16,18,25
Lifting tower 7	3,16,18,25

M

Magnetic stirrer → <i>Stirrer</i>	
Magnetic valve → <i>Electromagnetic valve</i>	
Mains cable	10,33
Mains connection	
Earthing	10
Fuse	10
Procedure	10
Mains cable	10
Mains socket 42	5
Mains voltage	10
Technical specifications	26
Mains socket 42	5
Mains switch 30	5,10
Mains voltage	
Changing the fuse	10
Display 35	5
Setting the voltage	10
Technical specifications	26
Manual operation	
Lifting arm	18
Rinsing	18
Sample transport	18
Siphoning	18
Stirrer	18
Mode 4	14
Mode 5	14
Model number	5

N

Nipple 11	3,7-8
Note 36	5
Note 40	5

O

Opening 50	7-8
Opening 51	8
Opening 52	7-8
Opening 53	8
Opening 54	7-8
Opening 55	7-8
Operation	
Manual operation	18
Notes	19
Operation scheme	17
Operation scheme	17
Options	33-35
Ordering designations	33-35
Outlet 19	3,9

P

Packing	7
PCB Connector (3.540.1640)	5,28
PCB Connector (3.540.2111)	5,30-31
pH Microelectrode (6.9902.048)	7-8,11,13,33-34
Plug socket 34	5
Plug socket 38	5
Plug socket 39	5
Power consumption	26
Power distributor (615)	35
Power on lamp 23	5,10
Program selection	14

Program sequence	14
PTFE Tubing (6.1812.000)	3,9-10,16,34
Pump (683)	
Electrical interconnections	11-13
Ordering designations	35

R

Reaction vessel → <i>Sample vessel</i>	
Reducing valve (6.2726.040)	34
Reference system (6.0724.140)	8,12,34
Rinse nozzle 14	3,7-8,16,33
Rinsing	
Key 25	5
Manual operation	18
Rinsing with air pump and electromagnetic valve	9,12
Rinsing with 683 Pump	9,11,13
Schematic view	9
Tubing connections	9
Rotary regulator 2	3,5,16,18

S

Safety specifications	26
Sample Changer (676)	
Accessories	34
Connecting cable to 664	3
Control elements	2-3
Figures	1-3
Front	2
General description	1
Ordering designations	33
Rear	3
Scope of delivery	33
Technical specifications	25
Titration head	8
Sample determination	17
Sample preparation	17
Sample rack 4	3,19,25,33
Sample transport	
Key 28	5
Manual operation	18
Sample vessel	
Figure	3
Ordering designations	33
Sample rack	19
Stirring bar	19
Technical specifications	25
Sample vessel 6	3
Saponification number → <i>Determination of saponification number</i>	
Scope of delivery	33
Screw driver 49	6,15
Selector switch 48	6,14
Serial number	5
SET/MET Titrimo (702)	
COD determination	24
Electrical interconnections	12
Ordering designations	35
Setting up the instruments	7
Siphoning	
Key 24	5
Manual operation	18
Schematic view	9
Siphoning with 683 Pump	9
Tubing connections	9
Siphoning tube	
Figure	3
Installation	7-8
Ordering designation	33-34

Siphoning tube 15	3,7-8
Spacer (6.2049.000)	34
Spiral band (6.1815.010)	10,33-34
Start of the analysis program	19
Startup procedure	17
Stirrer	
Connecting socket 33	5
Key 26	5
Manual operation	18
Rotary regulator 2	3
Rotational speed	3
Setting the stirring speed	16,18
Stirring bar (6.1906.010)	19,33-34
Stirring bar (6.1906.020)	19,34
Stop lamp 29	5
Stop of the analysis program	19
Supply bottle	9
Switching on the 664 Control Unit	10

T

Technical specifications	
Air pump (6.2817.000)	27
Control inputs and outputs	28-32
Control Unit (664)	26
Electromagnetic valve (6.2726.070)	27
Sample Changer (676)	25
Thumbwheel switch 31	5,19
Titration head	
Equipping for COD determination	8
Equipping for determination of hydroxyl and saponification numbers	7
Figure	8
Installation of accessories	7-8
Titration head 13	3,7-8,10,16,25
Titroprocessor (670)	
COD determination	22-23
Electrical interconnections	13
General description	1
Ordering designations	35
Titroprocessor (682)	
Determination of hydroxyl number	20-21
Electrical interconnections	11
Titroprocessor (686)	
COD determination	22
Determination of hydroxyl number	20-21
Electrical interconnections	11
Ordering designations	35
Transport damage	7,36
Tubing connections	
Filling	16
Installation	9
Turntable 16	3,5,16,18,25

V

Valve → <i>Electromagnetic valve</i>	
Valve opening times	14-15,18,26
Valve panel 21	3
Voltage selector switch	10

W

Warranty	36
Waste container	9
Weight	25-26