891 Professional Analog Out



2.891.0010

Manual 8.891.8001EN / 2019-12-16





Metrohm AG CH-9100 Herisau Switzerland Phone +41 71 353 85 85 Fax +41 71 353 89 01 info@metrohm.com www.metrohm.com

891 Professional Analog Out 2.891.0010

Manual

Technical Communication Metrohm AG CH-9100 Herisau techcom@metrohm.com

This documentation is protected by copyright. All rights reserved.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Table of contents

Table of contents

1	Instrument	description	1
2	About the	documentation	2
	2.1	Symbols and conventions	2
3	Safety inst	ructions	4
4	Overview o	of the instrument	5
	4.1	Front	5
	4.2	Rear	5
5	Installation	1	6
	5.1 5.1.1 5.1.2 5.1.3	Setting up the instrument Packaging Checks Location	6 6
	5.2	Connecting the 891 Professional Analog Out	7
6	Start-up		10
7	Operation		11
8	Technical s	pecifications	13
	8.1	Ambient temperature	13
	8.2	Housing	13
	8.3	Connectors	13
	8.4	Output channels	13
	8.5	Out 1 and Out 2 outputs	14
	8.6	Aux 1 and Aux 2 outputs	14
	8.7	Output and configuration possibilities	14
	8.8	Weight	15
9	Accessorie	s	16
	Index		17

Table of figures

Table of figures

Figure 1	Front	5
Figure 2	Rear	
Figure 3	Possible connections to the 891 Professional Analog Out	
Figure 4	Connecting the analog output	9

IV -----

1 Instrument description

1 Instrument description

With the **891 Professional Analog Out**, digital signals of Metrohm ion chromatographs can be converted to analog signals and processed in non-Metrohm software.

Like the ion chromatograph, the 891 Professional Analog Out is operated with the **MagIC Net** software. When the 891 Professional Analog Out is connected to an ion chromatograph or an extension module, MagIC Net recognizes it automatically. The software checks its functional readiness and controls and monitors all devices connected to one another.

The 891 Professional Analog Out is equipped with two high-resolution outputs, *Out 1* and *Out 2*, as well as with two low-resolution auxiliary outputs, *Aux 1* and *Aux 2*. All outputs can be used for transmitting detector data and information on pump flow, temperature, etc. The output of detector data on the low-resolution auxiliary outputs is not recommended.

2 About the documentation



CAUTION

Please read through this documentation carefully before putting the 891 Professional Analog Out into operation. The documentation contains information and warnings which the user must follow in order to ensure safe operation of the instrument.

2.1 Symbols and conventions

The following symbols and formatting may appear in this documentation:

(5- 12)	Cross-reference to figure legend		
	The first number refers to the figure number, the second to the instrument part in the figure.		
1	Instruction step		
	Carry out these steps in the sequence shown.		
Method	Dialog text, parameter in the software		
File ► New	Menu or menu item		
[Next]	Button or key		
	WARNING		
	This symbol draws attention to a possible life-threat- ening hazard or risk of injury.		
\wedge	WARNING		
	This symbol draws attention to a possible hazard due to electrical current.		
	WARNING		
<u></u>	This symbol draws attention to a possible hazard due to heat or hot instrument parts.		
	WARNING		
	This symbol draws attention to a possible biological hazard.		

 2 About the documentation



CAUTION

This symbol draws attention to possible damage to instruments or instrument parts.



NOTE

This symbol highlights additional information and tips.

3 Safety instructions



WARNING

Operate this instrument only according to the information contained in this documentation.



WARNING

Never open the housing of the instrument. The instrument could be damaged by this. There is also a risk of serious injury if live components are touched.

There are no parts inside the housing which can be serviced or replaced by the user.



WARNING

Only personnel qualified by Metrohm are authorized to carry out service work on electronic components.

4 Overview of the instrument

4 Overview of the instrument

4.1 Front




Figure 1 Front

1 Analog output 1

Main channel with high resolution, labeled *Out 1*.

3 Analog output 3

Auxiliary channel with low resolution, labeled *Aux 1*.

2 Analog output 2

Main channel with high resolution, labeled *Out 2*.

4 Analog output 4

Auxiliary channel with low resolution, labeled *Aux 2*.

4.2 Rear




Figure 2 Rear

1 Digital input

2 Type plate

5 Installation

5.1 Setting up the instrument

5.1.1 Packaging

The instrument is supplied in highly protective special packaging. Keep this packaging, as only this ensures safe transportation of the instrument.

5.1.2 Checks

Immediately after receipt, check whether the shipment has arrived complete and without damage by comparing it with the delivery note.

5.1.3 Location

The instrument has been developed for operation indoors and may not be used in explosive environments.

Place the instrument in a location of the laboratory which is suitable for operation and free of vibrations and which provides protection against corrosive atmosphere and contamination by chemicals.

The instrument should be protected against excessive temperature fluctuations and direct sunlight.

5 Installation

5.2 Connecting the 891 Professional Analog Out




Figure 3 Possible connections to the 891 Professional Analog Out



NOTICE

The total length of all 6.2156.0x0 connecting cables between the ion chromatograph and the 891 Professional Analog Out must not exceed 3.5 m.

Connecting the digital input

- Connect the cable (6.2156.070) to the connector *Digital Input* of the 891 Professional Analog Out.
- 2 Depending on the instrument, connect the other end of the cable to one of the following connectors:

On the instrument	to the connector
850 Professional IC / 940 Professional IC Vario	Extension Module
872 Extension Module / 942 Extension Module Vario	Out
881 Compact IC pro	Auxiliary
882 Compact IC plus	Auxiliary
930 Compact IC Flex	Auxiliary

Connecting the analog outputs

The cable (6.2128.200) has a plug on one end which fits to the analog outputs of the 891 Professional Analog Out. On the other end, a plug fitting to the target device must be mounted.

This is supported by the following color code:

- Green cable = positive voltage.
- Brown cable = negative voltage.
- White cable = negative voltage.
- Grounding

1 Connecting the analog cable to the 891 Professional Analog Out

- Connect the plug of the cable (6.2128.200) (4-1) to the required socket. Observe the alignment of the three contact pins.
- Tighten the knurled screw (4-2) at the front end of the plug by hand in clockwise direction to secure the plug.

5 Installation




Figure 4 Connecting the analog output

1 Cable 6.2128.200

2 Knurled screw

2 Connecting the analog cable to the target device

Plug the other end of the cable into the analog input of the target device.

6 Start-up

The 891 Professional Analog Out is put into operation together with the ion chromatograph it is connected to.

Putting ion chromatographs with 891 Professional Analog Out into operation

- 1 Start MagIC Net.
- Connect the ion chromatograph to the PC and switch it on.
 The 891 Professional Analog Out is recognized automatically by MagIC Net.

Additional information can be found in the *Start-up* chapter in the manual for the ion chromatograph as well as in the MagIC Net online help.

7 Operation

7 Operation

The **891 Professional Analog Out** is configured and the data sources are assigned to the output channels in the **MagIC Net** software. Additional information on operating the software can be found in the document "MagIC Net Tutorial" or in the online help.

Configuration

If the 891 Professional Analog Out is connected to an ion chromatograph or an Extension Module, it is recognized automatically when MagIC Net is started but it will not be entered into the device table as a stand-alone device.

Displaying configuration data

- 1 In the **Configuration** program part, open the properties of the IC instrument to which the 891 Professional Analog Out is connected.
- 2 Open the **Analog Out** tab.

The name and the serial number of the 891 Professional Analog Out are shown. These data cannot be modified.

The actual assignment of the data sources to the output channels takes place in the **Method** program part, on the four tabs for the analog outputs: **Out1 / Out2 / Aux1 / Aux2**. For this, the instrument to which the 891 Professional Analog Out is connected has to be selected in the device window. The four analog outputs can be found on the four tabs on the back.

Configuring output channels

- 1 Open the **Method** program part.
- 2 In the **Devices** subwindow, select the instrument to which the 891 Professional Analog Out is connected.
 - The data for the four analog outputs *Out 1 / Out 2 / Aux 1 / Aux 2* is to be found on the four tabs at the back.
- **3** The following indications can be made for each channel:

- Selection of the data source. The matching data channel is automatically assigned.
- Indication of the signal to be output as 1 volt (maximum output voltage).
- Selection of the sign of the output voltage.
- Optional: Indication of a value to be subtracted from every measured value before it is converted to the output voltage.

8 Technical specifications

8 Technical specifications

8.1 Ambient temperature

Operation $+5 - +45 \degree C$ Storage $-20 - +70 \degree C$ Transport $-40 - +70 \degree C$

8.2 Housing

Dimensions

Width 106 mm
Height 32 mm
Depth 123 mm

Material Aluminum

8.3 Connectors

Front 4 three-pin circular plugs for voltage output

Rear 1 D-sub plug, 15-pin

Power supply unit no power supply unit required

8.4 Output channels

Number 4 output channels

Absolute output 2,000 mV

range

Noise $< 2 \mu V \text{ rms}$

Channel assign- universal, without restrictions; e.g. detector signals, system pressure,

ment gradient percentages, column temperature, etc.

8.5 Out 1 and Out 2 outputs

Resolution 25 bit

Recommended as main channel for detector signals, e.g. conductivity 0 - 15,000

usage µS/cm

Output voltage -1,000 - +1,000 mV, bipolar

Accuracy typical: $< \pm 0.1\%$ deviation

8.6 Aux 1 and Aux 2 outputs

Resolution 16 bit

Recommended

usage

as auxiliary channel for system pressure, column temperature, etc.

Output voltage -1,000 - +1,000 mV, bipolar

Accuracy typical: $< \pm 0.2\%$ deviation

8.7 Output and configuration possibilities

Measured value absolute, in units of the source value, e.g. 10 μS/cm

offset

Voltage offset absolute, in mV

Full scale per volt rangeless for all channels in units of the source value per volt

Polarity positive, negative

Baseline reset baseline

Autozero programmable Autozero

8 Technical specifications

8.8 Weight

1.891.0010

300 g

9 Accessories

Up-to-date information on the scope of delivery and optional accessories for your product can be found on the Internet. You can download this information using the article number as follows:

Downloading the accessories list

- **1** Enter https://www.metrohm.com/ into your Internet browser.
- 2 Enter the article number (e.g. **2.891.0010**) into the search field. The search result is displayed.
- 3 Click on the product.

 Detailed information regarding the product is shown on various tabs.
- 4 On the **Included parts** tab, click on **Download the PDF**.

 The PDF file with the accessories data is created.



NOTICE

Once you have received your new product, we recommend downloading the accessories list from the Internet, printing it out and keeping it together with the manual for reference purposes.

Index

Index

F	Rear 5	R
Front 5		Rear 5
	0	
I	Overview of the instrument 5	S
Instrument		Service 4
Front 5		