

## Release Notes

EN

### 1 Purpose

This document describes new features as well as fixed bugs and problems in the following new firmware version:

- 5.917.0020

in comparison with the previous version:

- 5.917.0010

### 2 New features

- The document 8.840.4523 Titrando Compliance Guide 840/PC Control 4.1 has been updated for the 900 Touch Control and now includes 915 KF Ti-Touch, 916 Ti-Touch, and 917 Coulometer. New article number: 8.900.4007 Compliance Guide Touch Control.

### 3 Fixed Bugs and Problems

- If in the dialog **Dosing unit / PREP parameters** the volume was set to **cylinder volume**, the PREP run did not work correctly.
- If a determination was canceled in the dialog **Electronic signature**, the determination did not end correctly.
- If a report was manually created out of a determination by means of **Print ► Report ► Report seq. as def. in method**, under certain conditions the two chapters "Results" and "Messages" were printed on top of each other.
- When in a system a balance and a USB thermal printer Neo's were used, communication problems could occur.
- A bug in the Diagnosis area for the Metrohm Service was fixed.
- An attempt to load a Metrohm subsequence template in a language which does not support Metrohm method templates, lead to hidden Metrohm method templates upon returning to the dialog **Load method / New Method**.

The current firmware version of the 917 Coulometer complies with 21 CFR Part 11 and GAMP.

Metrohm AG, Ionenstrasse, CH-9100 Herisau/Switzerland

Metrohm AG holds the SQS certificate ISO 9001 – Quality management system for the development, production and sale of instruments and accessories for analyses.

Herisau, May 23rd, 2017

P. Hunzler

P. Hunziker  
Vice President,  
Head of Development

Dr. H. Loh

M. Saba  
Head of Quality Management