



Application Note AN-V-057

## 滴眼液中的硫柳汞

### Voltammetric determination of preservative in vaccines and various pharmaceutical and cosmetic solutions

Thiomersal (also called thimerosal) is a mercury containing organic molecule that has been widely used as preservative for vaccines and eye drops. Other applications of this substance have been for ear drops, storage and cleaning solutions for contact lenses, and in tattoo inks. It is very effective, even in very low concentrations, against a wide range of microorganisms and viruses.

To reduce the risk for consumers the maximum concentration of mercury in the products is limited by


the authorities. Typical limit values in the European Union are 0.007% (as Hg) in eye products (EU Commission Regulation 1223/2009) or up to 25 g thiomersal per vaccine injection dose if vaccines are distributed in multi-vaccine containers.

Polarography or voltammetry can be used to accurately determine the concentration of thiomersal in vaccines or other cosmetic and pharmaceutical solutions (such as eye drops). The method is simple to perform, specific, and free of interferences.

Vaccine against diphtheria, tetanus, pertussis

## EXPERIMENTAL

The sample and the supporting electrolyte are pipetted into the measuring vessel. The determination of thiomersal is carried out with an 884 Professional VA using the parameters specified in **Table 1**. The concentration is determined using external calibration with a calibration curve recorded from five standard concentrations.



**Figure 1.** 884 Professional VA.

**Table 1.** Parameters for voltammetric analysis of thiomersal in vaccine


Parameter	Setting
Working electrode	HMDE or DME
Mode	DP – Differential Pulse
Deposition	none
Start potential	-0.2 V
End potential	-0.6 V
Peak potential thiomersal	-0.45 V

## ELECTRODES

- Working electrode: Multi-Mode Electrode pro with standard glass capillaries
- Reference electrode: Ag/AgCl/KCl (3 mol/L) reference electrode with electrolyte vessel.  
Bridge electrolyte: KCl (3 mol/L)
- Auxiliary electrode: Platinum rod electrode

## RESULTS

The determination of thiomersal can be carried out in a simple and straightforward manner with voltammetry. The method is selective and free of interferences.



**Figure 2.** Determination of thiomersal in vaccine sample.

**Table 2.** Results of the thiomersal determination with the 884 Professional VA

Sample	Thiomersal [mg/L]
DTP vaccine	102

## RESULTS

Internal reference: AW BR4-0002-072011

## CONTACT

Metrohm AG Ionenstrasse 9100 Herisau	info@metrohm.com
--	------------------

## CONFIGURATION



(MME) 884 Professional VA manual  
(MME) 884 Professional VA manual pro  
scTRACE Gold / viva ,,

,«»(CVS)«»(CPVS)(CP),

**viva**

MME() 884 Professional VA manual , pro **viva**



VA pro Professional VA  
, pro