

1 Purpose

This document describes the differences between the software version **MagIC Net 2.3 (Build 87)** and the version **MagIC Net 2.3 (87 - 1)**.

2 Fixed bugs and problems

Program part Workplace

- When running multiple analyses in parallel, MagIC Net sometimes became inoperable. If live display was deactivated then the problem did not occur.

Program part Database

- An exception error occurred, if for the display of chromatograms the option **Scaling** in the dialog window **Properties - Graphics** was changed from **Relative** to **Absolute** or from **Absolute** to **Relative**, and if comma was used as decimal separator.
- In rare cases, when opening a determination in the dialog window **Reprocessing** only the baseline and the peak labels were shown. The chromatogram was missing.
- Changed default profile in Mode **PAD**.

3 Compliance

The current software version does not contain any modifications that affect conformity of **MagIC Net** regarding 21 CFR Part 11 and GAMP.

Herisau, February 28, 2012



Dr. J. Kleimann
Vice President, Head of R&D



U. Kürsteiner
Head of Quality Management