

Dosing Package pH Meter Package Conductivity Package



Complete packages for
professional liquid handling and
pH and conductivity measurement

Dosing Package for liquid handling ...

02

The Dosing Package is the practical and proven solution for precise, reliable liquid handling in accordance with directives, both in regulated and non-regulated environments.

You can rely on the Dosing Package for:

- Standard applications such as dosing, diluting, taking aliquots, pipetting
- The preparation of standard solutions from any given sample size
- The preparation of mixtures from various stock solutions



Your benefit

- Fulfills requirements of FDA regulation 21 CFR Part 11
- Each application can be programmed as a method and started at the press of a button
- Complete monitoring and traceability of applications and results
- Thanks to USB and Ethernet interface, a wide range of options for data backup (LIMS, PDF report, USB flash drive, local or network printer)
- Completely compatible with other Metrohm Titrandos instruments (pH and conductivity measurement, titration, automation)

Order number

- Dosing Package: 2.846.1010

Dosing Package


... and pH Meter Package for pH-value/ ion measurement and Conductivity Package for conductivity measurement

Both the 867 pH Module and the 856 Conductivity Module with Touch Control are measuring instruments that meet the highest standards.

03

They are particularly suitable for:

- Measurement of the pH value, the ORP, the temperature or an ion concentration
- pH measurement in accordance with USP 791
- Ion measurement using fully automated standard addition (with optional 800 Dosino)
- Conductivity measurement in accordance with USP <645>



Your benefit

- Fulfills requirements pursuant to FDA regulation 21 CFR Part 11
- Each application can be programmed as a method and started at the press of a button
- Complete monitoring and traceability of applications and results
- Thanks to USB and Ethernet interface, a wide range of opportunities for data backup (LIMS, PDF report, USB flash drive, local or network printer)
- Completely compatible with other Metrohm Titrandos instruments (titration, automation)

An optional 800 Dosino provides the 867 pH Module and the 856 Conductivity Module the same liquid handling functionality as the 846 Dosing Interface

Order numbers

- **pH Meter Package: 2.867.1410**
- **Conductivity Package: 2.856.0120**

Ordering options

2.846.1010 Dosing Package

	846 Dosing Interface
	800 Dosino
	900 Touch Control
6.2061.010	Bottle holder for Dosino
6.2065.000	Stacking frame for 846
6.3032.220	Dosing unit 20 mL
6.6064.010	USB flash drive for 900

2.867.1410 pH Meter Package

	867 pH Module
	900 Touch Control
6.6064.010	USB flash drive for 900
	including stand and electrode holder (sensor not included*)

2.856.0120 Conductivity Module with Touch Control and steel cell

	Conductivity Module
	900 Touch Control
6.0916.040	Conductivity measuring cell, stainless steel with Pt1000 m C=0.1
6.6064.010	USB flash drive for 900
	including stand and electrode holder

Optional accessories

2.800.0010	800 Dosino
2.801.0040	Magnetic stirrer MSB
6.3032.120	Dosing unit 2 mL
6.3032.150	Dosing unit 5 mL
6.3032.210	Dosing unit 10 mL
6.3032.220	Dosing unit 20 mL
6.3032.250	Dosing unit 50 mL

*You'll find the right electrode for your application under
www.metrohm.com/com/products/electrodes
or with our Electrode Finder App!



www.metrohm.com

