

Border security with Metrohm identification systems

Detect hazards fast
and accurately at
the touch of a screen



Mira – your material identification powerhouse

Customs plays a crucial role in border security. Challenges are diverse, ranging from counterfeits to illegal substances and imports all the way to explosives and toxic substances. With Mira handheld identification systems from Metrohm, you are prepared to detect them all – safely and accurately, at the touch of a screen.



Use Mira handheld identification systems to immediately...

- perform antifraud operations by determining the authenticity and origin of products
- identify illegal imports like narcotics and explosives
- protect consumers against dangerous goods or contaminated food
- safeguard the environment and identify endangered species
- authenticate documents
- combat terrorism

Instant answers for immediate action

Time can be everything: With Mira handheld identification systems, you benefit from on-the-spot-control and measurements in the field – no waiting for results from a laboratory.





A forensic lab in the field

Handheld analysis in every situation

Sampling is the first, crucial step in customs work – and can be difficult and dangerous. Not so with Mira handheld identification systems.

- Various sampling attachments keep you safe from direct contact with the sample
- Mere seconds from sampling to an accurate result
- Large libraries of substances
- Mixture analysis



1 Point-and-shoot Attachments 2 Contact Ball Probe 3 Tablet Holder 4 Vial Holder 5 Calibrate/Verify Accessory (CVA)
6 Universal Attachment 7 Standoff Attachment 8 Right Angle Attachment

Watch a demo on
www.metrohm.com/MIRA-DS



Ordering information

2.926.0020 Mira DS (Advanced Package, in carry case complete with):

- 1.926.0020 Mira DS Advanced
- 6.07501.010 Calibration Standard
- 6.07506.000 Right Angle Attachment
- 6.07506.010 Universal Attachment
- 6.06071.020 Mira Cal DS USB Stick
- 6.02707.000 Carrying case for Mira
- 6.2133.020 Lithium battery AA 1,5V/2900mAh
- 6.2151.110 Cable USB (OTG) Mini B PI-USB A PI 1.8m
- 6.2166.000 USB power adapter
- 6.7560.010 Laser goggles 785 nm
- 6.7430.060 NIR micro fiber cleaning cloth

Accessories

- 6.06071.020 Mira Cal DS USB Stick
- 6.07506.000 Right Angle Attachment
- 6.07506.010 Universal Attachment
- 6.07506.020 Standoff Attachment
- 6.07506.030 Contact Ball-Probe
- 6.06071.640 HazMasterG3 Software (1 license)

www.metrohm.com/Mira-DS

 **Metrohm**
Raman