

Solvotrode *easyClean* – cleaning at the touch of a button



NEW: The proven Solvotrode with the new *easyClean* diaphragm – clean and easy!



Easy cleaning of the diaphragm: Just press once on the electrode head and the electrolyte flows out. The diaphragm does not need to be touched any more for cleaning.

Highlights

- Easy contact-free cleaning of the diaphragm
- Greater accuracy and reproducibility of the electrolyte flow (glass on glass, the spring in the head returns to the defined starting position)
- Reduced immersion depth of the sensor

Ordering information

- 6.0229.010 Solvotrode *easyClean*, fixed cable, plug F
6.0229.020 Solvotrode *easyClean*, fixed cable 2 m, plug F

Visit us at www.metrohm.com

