

Combined Calcium ion-selective electrode for direct measurement and titration



Space-saving and robust



The advantages at a glance:

- **Space-saving** – the combined two-in-one sensor gives you the opportunity to use smaller sample vessels and less sample for analysis
- **Ideal for the titration of drinking water** because of a chloride-free electrolyte
- **Robust construction** with plastic shaft and protected membrane
- For titration, direct measurement or standard addition

Ordering information

6.0510.100	Combined Calcium Ion-Selective Electrode
6.2104.020	Electrode cable, plug-in head G, socket F, 1 m
6.2327.000	Electrolyte NH_4NO_3 , 1 mol/L

Visit us at www.metrohm.com