

1 Purpose

This document describes new functions, the changes, improvements and fixed bugs and problems in the following firmware version:

- 57.1008.0011 for OS version 005
- in comparison with the firmware version:
- 57.1008.0010 for OS version 004

2 New functions

The following new functions have been integrated:

- Support for PostScript USB printers
- Sample series
- Support for Sartorius balances via COM port (USB/RS-232 Converter required)
- COM port settings for balances

3 Improvements

The following improvements were made:

- Firmware optimizations for the Russian language.

4 Fixed bugs and problems

The following bugs and problems were fixed:

- Some translation mistakes were corrected.
- OHMTI-1087 – The fixed points FP1 and FP2 are now displayed.

5 Known bugs and problems

The following bugs or problems are known:

- Problems with saving or loading methods / sensors / solutions while a sample series is being processed.