

1 Purpose

This document describes the differences between the new software version **StabNet 1.1** and its predecessor **StabNet 1.0**.

2 New features

New instrument firmware

Instruments	Firmware version	Comment
892 Professional Rancimat	5.893.0014	The new firmware can only be installed if version 5.893.0012 is already running on the instrument. If the instrument still uses an older firmware version, please contact Metrohm Service.
893 Professional Biodiesel Rancimat	5.893.0014	
895 Professional PVC Thermomat	5.893.0014	

Operating systems

With version 1.1 and higher, StabNet runs on the following operating systems:

- **32-bit version of**
 - Windows 8/8.1 Professional/Enterprise
 - Windows Server 2012 R2
- **64-bit version of**
 - Windows 8/8.1 Professional/Enterprise
 - Windows Server 2012 R2

3 Compliance

The current software version does not contain any modifications that affect conformity of **StabNet** regarding 21 CFR Part 11 and GAMP.

Herisau, August 29, 2014



P. Hunziker
Vice President,
Head of Development



U. Kürsteiner
Head of Quality Management