



Polytron PT 1300 D

2.136.0100

Polytron PT 1300 D - Metrohm Version

Homogenisierer, welcher direkt durch die OMNIS Software, tiamo™ oder Touch Control angesteuert werden kann.

Der Polytron PT 1300 D besteht aus einem Kontrollgerät und einem Antrieb. Das Kupplungssystem am Antrieb erlaubt den einfachen und schnellen Austausch der Aggregate ohne zusätzliches Werkzeug.

Feste Proben können mühelos zerkleinert werden. Das Gerät ist aber auch hervorragend geeignet um in viskosen Proben für eine gute Durchmischung zu sorgen.

Lieferumfang 2.136.0100

Qt.	Order no.	Beschreibung
-----	-----------	--------------

1 PCS

1.136.0100

Polytron PT 1300 D

Polytron PT 1300 D - Metrohm version


Homogenizer that can be controlled directly by tiamo™ or Touch Control.




1 PCS

1.136.0110

Power supply unit for Polytron PT 1300 D



Optionales Zubehör

Order no.	Beschreibung	
6.02008.010	Polytron holder For fastening the Polytron in place, suitable for 601405040.	
6.02008.030	Holder for homogenizer for Eco Titrations instruments Holder for homogenizer for fastening the Polytron in place, suitable for 601405060 titration vessel lid	
6.1462.250	Swing arm for Polytron and 786 Swing Head, left swinging Swing arm for Polytron and 786 Swing Head, left swinging, with one buret tip for the addition of solvent or titrant and three spray nozzles for cleaning. The swing arm allows the usage of the Polytron for sample preparation on multi row racks.	
6.2047.110	Holder for Polytron Holder to fix the Polytron at the stand rod.	

6.9012.000

Dispersing aggregate for Polytron, 125 mm

Standard dispersing aggregate for Polytron 1300 D



6.9012.010

Dispersing aggregate for Polytron, 157 mm

Special dispersing aggregate (Metrohm version) with protruding knives

