




801 Stirrer mit Stativ

2.801.0040

Magnetrührer inklusive Bodenplatte, Stativstange und Elektrodenhalter zur Ergänzung von Titrino plus, Dosimat plus, Titrandos, Sample Processoren, 805 Dosimat und 780/781 pH Meter sowie der 856 und 867 Messmodule. Mit fest montiertem Kabel für MSB (Metrohm Serial Bus).

Lieferumfang 2.801.0040

Qt.	Order no.	Beschreibung
1 PCS	1.801.0010	801 Stirrer

Magnetic stirrer without stand for supplementing the Titrino plus, Dosimat plus, Titrandos, Sample Processors, 805 Dosimats and 780/781 pH Meters. With permanently attached cable for MSB (Metrohm Serial Bus).



1 PCS

6.1903.020

Stirring bar / 16 mm

Stirring Bar with magnetic core, PTFE coated, length 16 mm



1 PCS

6.1903.030

Stirring bar / 25 mm

Stirring bar with magnetic core, PTFE covering, length 25 mm.




1 PCS

6.2001.060

Support stand without support rod

To attach an 801, 804, 803 Stirrer to a Titrando, Titrino plus, Dosimat plus or 856 and 867.



1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2016.070

Support rod / 400 mm



1 PCS

6.2021.020

Electrode holder

Electrode holder for 4 electrodes and 2 buret tips



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors



1 PCS

6.2621.130

Hexagon key 2 mm

2 mm.



Optionales Zubehör

Order no.	Beschreibung	
6.2001.090	Base plate on both sides Double stand for attaching 2 x 801, 803, 804 to a Titrando	
6.2151.010	Cable MSB plug / socket Extension cable for MSB connections and for Touch Control - Titrando, Mini-DIN socket - mini-DIN plug	