



846 Dosing Interface

2.846.0010

USB-fähige Steuereinheit für den Anschluss von max. vier 800 Dosinos oder 805 Dosimaten für Dosier- und Liquid-Handling-Aufgaben. Zur Bedienung ist ein Touch Control oder der Anschluss an einen PC mit OMNIS Software, tiamo™, MagIC Net, viva oder 797 VA Computrace nötig.


Lieferumfang 2.846.0010

Qt.	Order no.	Beschreibung
1 PCS	1.846.0010	846 Dosing Interface

USB-capable control unit for connecting a maximum of four 800 Dosinos or 805 Dosimats for dosing and Liquid Handling tasks. A Touch Control or the connection to a PC with tiamo™, MagIC Net, viva or 797 VA Computrace is required for operation.



Optionales Zubehör

Order no.	Beschreibung
2.800.0010	800 Dosino The 800 Dosino is a drive with write/read hardware for intelligent Dosing Units. With fixed cable (length 150 cm). 
2.900.0010	900 Touch Control Operating unit for the Titrandos, USB Sample Processors, 856 Conductivity Module, 867 pH Module and 846 Dosing Interface. Touch-sensitive, high-resolution color display, simple and intuitive operation, thanks to Favorites for direct method access. With integrated Ethernet interface for direct connection to the Internet and USB interface for connecting USB printers or a USB memory stick. Dialog languages: German, English, Chinese, French, Spanish, Portuguese, Russian, Korean, Polish and Italian. 
6.3032.120	807 Dosing Unit 2 mL 807 Dosing Unit with integrated data chip with 2 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip. 
6.3032.150	807 Dosing Unit 5 mL 807 Dosing Unit with integrated data chip with 5 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip. 

6.3032.210

807 Dosing Unit 10 mL

807 Dosing Unit with integrated data chip with 10 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.



6.3032.220

807 Dosing Unit 20 mL

807 Dosing Unit with integrated data chip with 20 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.



6.3032.250

807 Dosing Unit 50 mL

807 Dosing Unit with integrated data chip with 50 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.

