



dPt-Titrode

6.00401.300

Digitale, kombinierte Platinringeletrode für OMNIS mit einer pH-Glasmembran als Referenzelektrode.

Diese wartungsfreie Elektrode eignet sich für Redox titrationen bei konstant bleibendem pH-Wert, z.B.

- Iodometrie
- Chromatometrie
- Cerimetrie
- Permanganometrie

Diese Elektrode wird in dest. Wasser aufbewahrt.

dTrodes können an OMNIS Titratoren verwendet werden.

Lieferumfang 6.00401.300

Qt.	Order no.	Beschreibung
-----	-----------	--------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Optionales Zubehör

Order no.	Beschreibung	
6.02104.300	Electrode cable plug-in head Q / plug P, 0.55 m Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	
6.02104.310	Electrode cable plug-in head Q / plug P, 1.5 m Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	