



dAg-Ringelektrode

6.00402.300

Digitale, kombinierte Silberringelektrode für OMNIS mit einem Festschliffdiaphragma.

Diese Elektrode eignet sich für Fällungstitrationsen bei variierendem pH-Wert (Titriermittel Silbernitrat) von z.B.

- Chlorid, Bromid, Iodid
- Sulfiden
- Schwefelwasserstoff
- Mercaptane
- Cyanide

Das Festschliffdiaphragma ist gegenüber Verschmutzung unempfindlich. Als Referenzelektrolyt und zur Aufbewahrung wird $c(\text{KNO}_3) = 1 \text{ mol/L}$ verwendet.

dTrodes können an OMNIS Titratoren verwendet werden.

Lieferumfang 6.00402.300

Qt.	Order no.	Beschreibung
-----	-----------	--------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Optionales Zubehör

Order no.	Beschreibung	
6.02104.300	Electrode cable plug-in head Q / plug P, 0.55 m Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	
6.02104.310	Electrode cable plug-in head Q / plug P, 1.5 m Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	