



Kombinierte dF-ISE mit Pt1000

6.00500.300

Kombinierte digitale Fluoridselektive Elektrode mit Kristallmembran für OMNIS mit integriertem Temperaturfühler Pt1000. Diese ISE eignet sich für:

- Ionenmessungen von F^- (10^{-6} mol/L bis sat.)
- Automatisierte Ionenmessungen
- Titrationsen

Referenzelektrolyt: $c(KCl) = 3$ mol/L

Die Elektrode wird im Referenzelektrolyt gelagert.

dTodes können an OMNIS Titratoren verwendet werden.

Lieferumfang 6.00500.300

Qt.	Order no.	Beschreibung
-----	-----------	--------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Optionales Zubehör

Order no.	Beschreibung	
6.02104.300	Electrode cable plug-in head Q / plug P, 0.55 m Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	
6.02104.310	Electrode cable plug-in head Q / plug P, 1.5 m Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	
6.2308.020	Electrolyte 3 mol/L KCl (250 mL) Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)	
6.2308.050	Electrolyte 3 mol/L KCl (50 mL) Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$ (for Ag/AgCl reference systems)	

6.2313.000

Electrolyte 3 mol/L KCl (1000 mL)

Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$ (for Ag/AgCl reference systems)

