



Messkopf (SPE) für den Messstand des 946 Portable VA Analyzer

6.01256.020

Messkopf mit eingebautem Rührer für Messungen mit Dickfilmelektroden (*screen-printed electrodes, SPE*).

Optionales Zubehör


Order no.	Beschreibung
6.01204.000	Stirrer tip

Stirrer tip for 946 Portable VA Analyzer.



6.01408.000 O-ring 46 x 2.5 mm

For the 946 Portable VA Analyzer measuring stand. Set of 2 pieces.



6.02621.000 4 mm wrench



6.1241.200 Nut for SPE electrode shaft

Nut for SPE electrode shaft 6.1241.090 and measuring head (SPE) 6.01256.020.



6.1241.210 Supporting ring for SPE electrode shaft

Supporting ring made of PTFE for SPE electrode shaft 6.1241.090 and measuring head (SPE) 6.01256.020.



6.1244.060

Silicone seal for SPE electrode shaft (10 pieces)

Seal made of silicone for SPE electrode shaft 6.1241.090 and measuring head (SPE) 6.01256.020, 10 pieces.



DRP-11L-U75

Screen-Printed Carbon Electrode (Aux.:C;Ref.:Ag/AgCl)

Screen-Printed Carbon Electrode (Aux.:C;Ref.:Ag/AgCl)

