



Biotrode

6.0224.100

Kombinierte pH-Elektrode für Messungen in sehr kleinen Probenvolumina (>50 uL) und biologische Proben.

Idrolyt (6.2308.040) wird als Referenzelektrolyt und Aufbewahrungslösung verwendet.

Lieferumfang 6.0224.100

Qt.	Order no.	Beschreibung
1 PCS	6.1236.050	Sleeve with SGJ 14/12 mm
		Sleeve with SGJ 14/12 mm, polyethylene.



Together with 6.2008.050 storage vessel holder, it provides a support for the electrode on 807 Dosing Units.



Optionales Zubehör

Order no.	Beschreibung	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	
6.2308.040	Idrolyte (250 mL) Idrolyte (for 6.0224.100 Electrode and Unitrode at T > 80°C)	
6.2325.100	Cleaning solution 3 x 50 mL Reliable measuring results over long periods of time can be guaranteed only if the pH glass membrane and the diaphragm receive preventive and regular care. Cleaning by etching with toxic chemicals or applying mechanical treatment to the diaphragm is not only complicated and expensive, it also accelerates the ageing of the pH glass electrode. The cleaning solution was developed for easy and gentle cleaning of pH glass electrodes. Regular use can considerably prolong their service life. This cleaning solution is also part of the pHit kit.	