




IC-Ausrüstung: Inline-Ultrafiltration 2 - push mode

6.05330.110

Zubehörset zur Inline-Ultrafiltration 2 - push mode. Zur Verwendung mit dem 858 Professional Sample Processor / 919 IC Autosampler plus.

Empfohlenes Set zu Kombination mit Inline-Verdünnung (MIDT).

Lieferumfang 6.05330.110

Qt.	Order no.	Beschreibung
1 PCS	6.02057.030	Ultrafiltration cell 2 holder

Holder for mounting the 6.02729.110 ultrafiltration cell 2 on the 858 Professional Sample Processor / 919 IC Autosampler plus.



For the Inline ultrafiltration of samples together with the 858 Professional Sample Processor / 919 IC Autosampler plus. As an alternative to the time-consuming manual filtration with a syringe filter.



Capillary for sample handling in IC.



Capillary for inline dialysis and filtration, for Dialysis Unit, IC Filtration Sample Processor



For suppressor solutions and for acceptor solutions for inline dialysis and for inline ultrafiltration.



1 PCS

6.1826.390

Pump tubing LFL (yellow/yellow), 3 stoppers

For the sample solution in Inline Ultrafiltration.



1 PCS

6.1831.160

PEEK capillary, 0.5 mm ID / 70 cm



1 PCS

6.2621.030

Hexagon key 4 mm



1 PCS

6.2714.020

Filtration membrane reg. cellulose, 0.2 µm, 50 pieces

Nominal pore size: 0.2 µm. Set of 50 pieces. For inline filtration.



1 PCS

6.2744.000

Pressure screw PVDF

With UNF 10/32 connection. For IC instruments, inline dialysis
(connection of PTFE capillaries)



1 PCS

6.2744.034

Coupling olive/UNF 10/32, 2x

Connection between pressure screw and pump tubing. 2 pieces. For IC
instruments with peristaltic pump.



1 PCS

6.2744.070

Pressure screw short

Short version. With UNF 10/32 connection. 5 pieces. For the connection
of PEEK capillaries




2 PCS

6.2744.160

Nozzle/UNF 10/32 coupling with security device

Security device for connecting pump tubing to a nozzle.



1 PCS

6.2831.010

Plastic tweezers

For dialysis and filtration membranes



Optionales Zubehör

Order no.	Beschreibung
6.5330.120	IC equipment: Inline dilution
	Accessory set for Inline Dilution. For use with the 858 Professional Sample Processor, 800 Dosino and 741 Magnetic Stirrer.

