



LL ISE Referenzelektrode

6.0750.100

Silber/Silberchlorid-Referenzelektrode mit Double-Junction System.

Diese Referenzelektrode eignet sich hervorragend für:

- automatisierte Anwendungen
- Ionenmessungen
- Tensidtitrationen

Das gegenüber Verschmutzung unempfindliche Schliffdiaphragma liefert einen konstanten und reproduzierbaren Elektrolytausfluss. Zudem wird der Referenzelektrolyt für noch bessere Signalstabilität geliert. Der Sensor wird mit $c(\text{KCl}) = 3 \text{ mol/L}$ als Zwischenelektrolyt ausgeliefert, dieser kann aber je nach Verwendung frei gewählt und ausgetauscht werden.

Lieferumfang 6.0750.100

Qt.	Order no.	Beschreibung
-----	-----------	--------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Optionales Zubehör

Order no.	Beschreibung	
6.2103.140	Adapter black 2 mm plug / B socket 4 mm For connecting plug B (4 mm) to 2 mm socket.	
6.2106.020	Strand / 1 m / 2 x B Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B	
6.2106.060	Strand / 2 m / 2 x B Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B	
6.2308.020	Electrolyte 3 mol/L KCl (250 mL) Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)	

6.2310.010

Electrolyte KNO₃-1 mol/L 250 mL

Electrolyte solution KNO₃ 1 M (reference electrolyte for combined silver electrodes)



6.2320.000

TEABr 0.4 mol/L in ethylene glycol (250 mL)

Electrolyte solution TEABr (tetraethylammonium bromide in ethylene glycol), c(TEABr) = 0.4 mol/L

