



Metrosep Carb 2 - 100/4.0

6.1090.410

Die IC-Säule Metrosep Carb 2 - 100/4.0 eignet sich speziell für die Bestimmung von Kohlenhydraten unter Verwendung alkalischer Eluenten und gepulster amperometrischer Detektion. Die hochkapazitive Anionenaustauschsäule basiert auf einem Styrol/Divinylbenzol-Copolymer. Sie ist im Bereich von pH = 0–14 stabil und trennt Glucose, Fructose und Sucrose. Darüber hinaus ist sie auch für die Analyse von einigen Zuckeralkoholen und Oligosacchariden geeignet. Auf der 100-mm-Variante der Metrosep Carb 2 Trennsäule lassen sich kurze Analysenzeiten erreichen.

Lieferumfang 6.1090.410

Qt.	Order no.	Beschreibung
2 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Optionales Zubehör

Order no.	Beschreibung	
6.1015.200	Metrosep BO3 Trap 1 - 100/4.0 Trap column for the removal of borate contaminants from the eluent. The Metrosep BO ₃ Trap 1 - 100/4.0 is used mainly in carbohydrate analysis with hydroxide eluents.	
6.1015.300	Metrosep CO3 Trap 1 - 100/4.0 Trap column for the removal of carbonate traces in hydroxide eluents.	
6.1090.500	Metrosep Carb 2 Guard/4.0 The Metrosep Carb 2 Guard/4.0 effectively removes contamination, thus protecting the analytical separation column. The design of the guard column has been selected in such a way that its influence on the chromatographic separation can be ignored.	
6.1090.510	Metrosep Carb 2 S-Guard/4.0 The Metrosep Carb 2 S-Guard/4.0 effectively removes contaminations, thus protecting the analytical separation column. The design of the guard column has been selected in such a way that its influence on the chromatographic separation can be ignored.	