



# Metrosep Carb 2 - 250/4.0

6.1090.430

Die IC-Säule Metrosep Carb 2 - 250/4.0 eignet sich speziell für die Bestimmung von Kohlenhydraten unter Verwendung alkalischer Eluenten und der gepulsten amperometrischen Detektion. Die hochkapazitive Anionenaustauschsäule basiert auf einem Styrol/Divinylbenzol-Copolymer. Sie ist im Bereich von pH = 0–14 stabil und trennt Mono- und Disaccharide. Darüber hinaus ist sie auch für die Analyse von Zuckeralkoholen, Anhydrozucker, Aminos Zucker usw. geeignet. Die 250-mm-Variante der Metrosep Carb 2 Trennsäule ist für komplexe Trennungen optimiert.

## Lieferumfang 6.1090.430

| Qt.   | Order no.  | Beschreibung     |
|-------|------------|------------------|
| 2 PCS | 6.2744.060 | Threaded stopper |

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



## Optionales Zubehör

| Order no.  | Beschreibung                                                                                                                                                                                                                                                                                 |                                                                                       |
|------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| 6.1015.200 | Metrosep BO3 Trap 1 - 100/4.0<br><br>Trap column for the removal of borate contaminants from the eluent. The Metrosep BO <sub>3</sub> Trap 1 - 100/4.0 is used mainly in carbohydrate analysis with hydroxide eluents.                                                                       |    |
| 6.1015.300 | Metrosep CO3 Trap 1 - 100/4.0<br><br>Trap column for the removal of carbonate traces in hydroxide eluents.                                                                                                                                                                                   |    |
| 6.1090.500 | Metrosep Carb 2 Guard/4.0<br><br>The Metrosep Carb 2 Guard/4.0 effectively removes contamination, thus protecting the analytical separation column. The design of the guard column has been selected in such a way that its influence on the chromatographic separation can be ignored.      |  |
| 6.1090.510 | Metrosep Carb 2 S-Guard/4.0<br><br>The Metrosep Carb 2 S-Guard/4.0 effectively removes contaminations, thus protecting the analytical separation column. The design of the guard column has been selected in such a way that its influence on the chromatographic separation can be ignored. |  |