

# Probenbecher aus Klarglas 250 mL, 1 Stück

6.1432.320

Probenbecher aus Klarglas 250 mL, passend zu Probenracks 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820, 6.2068.020 und 6.2068.030, 1 Stück

## Optionales Zubehör

Order no.	Beschreibung	
6.1458.010	Titration head, 6 x SGJ 14 and 3 x SGJ 9 openings.	
	With 6xSGJ 14 and 3x SGJ 9 openings.	



Propeller stirrer, length from lower edge of SGJ: 96 mm. For use with beakers with 722, 802 Rod Stirrer



6.2037.060

Lid for 250 mL sample beaker

Lid for 250 mL sample vessels in use with the Dis-Cover robotic arm on the Sample Processor.



6.2037.070

Lid with hole for 250 mL sample beaker


The lid with hole is used in the Soliprep system to pre-rinse applied syringe filters into the special beaker.



6.2041.310

Sample rack, 12 x 250 mL, (PVC)

Sample rack for 810 Sample Processor, 814 USB Sample Processor, 855 Robotic Titrosampler "Light", suitable for 12 sample beakers. The following sample beakers can be used: 6.1432.320, 6.1432.323, 6.1453.220, 6.1453.250.





Sample rack for 815 Robotic USB Sample Processor XL, 855 Robotic Titrosampler "Basic" for 126 sample vessels and 2 sample beakers. The following sample vessels can be used: 6.02743.050, 6.2743.057, 6.1432.320, 6.1432.323, 6.1453.220, 6.1453.250.

Plastic: Polyvinyl chloride (PVC)



#### 6.2041.510 Sample rack, 56 x 11 mL and 1 x 250 mL (PVC)

Sample rack for 919 Autosampler plus, 858 Professional Sample Processor suitable for 56 sample vessels and 1 rinsing bottle. The following sample vessels and rinsing bottle can be used: 100 mL digestion tubes (CSB).



Plastic: Polyvinyl chloride (PVC)

### 6.2041.820 Sample rack, 28 x 250 mL, (PVC)

Sample rack for 815 Robotic USB Sample Processor XL, 855 Robotic Titrosampler "Basic" suitable for 28 sample beakers. The following sample beakers can be used: 6.14332.320, 6.1432.323, 6.1453.220, 6.1453.250.



Plastic: Polyvinyl chloride (PVC)

#### 6.2068.020 Sample rack, 24 x 120 mL, (PVC)

Sample rack for 815 Robotic Soliprep, suitable for 24 sample beakers. The following sample beakers can be used: 6.1459.300, 6.1459.307.







Sample rack for 815 Robotic Soliprep, suitable for 10 sample beakers. The following sample beakers can be used: 6.1432.320, 6.1432.323, 6.1453.220, 6.1453.250.

