



Spritze 10 mL mit Luer-Anschluss

6.2816.020

Für verschiedene Anwendungen in IC und VA

Optionales Zubehör

Order no.	Beschreibung
6.2744.120	Luer F pressure screw adapter Connection between pressure screw and syringe.



6.2816.030

Needle with Luer connector

Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration.



6.2816.040

Purge tubing

With PTFE tubing and Luer connection. For syringes. For aspirating eluents.

