

Zubehör-Kit: Vario/Flex SeS

6.5000.020

Zubehörsatz zur sequentiellen Suppression (SeS) für Professional IC Vario und Compact IC Flex Geräte.Ohne Suppressor-Rotor.

Lieferumfang 6.5000.020

| 04 | Orderne | Doodyraibung |
|-------|------------|--|
| Qt. | Order no. | Beschreibung |
| 1 PCS | 6.1602.150 | Bottle attachment / GL 45 - 3 x UNF 10/32 |
| | | on of 1/16 in. capillaries. Used in conjunction with MSM and in Metrohm Inline Sample Preparation. |



Bottle for auxiliary solutions.



1 PCS 6.1803.030 PTFE capillary 0.5 mm i.d. / 3 m

Capillary for inline dialysis, for Dialysis Unit, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit



2 PCS 6.2744.040 2 x UNF 10/32 coupling

For connecting 1/16 in. capillaries. For IC instruments



1 PCS 6.2744.070 Pressure screw short

Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries





Long version. With UNF 10/32 connection. 2 pieces. For the connection of PEEK capillaries (MCS and Sample degasser)



1 PCS

6.2835.010

Connecting piece to MSM and SPM module

Universal connecting piece for MSM in Professional IC Vario and Compact Flex instruments and for MSM-HC and SPM module. With inlets and outlets.



1 PCS

6.2837.100

CO2 Absorber

Absorber cartridge for removing the ${\rm CO_2}$ from the air aspirated by the MCS. The cartridge is used without upstream ${\rm H_2O}$ absorber cartridge.

