



VA-Elektrodenausrüstung mit Bismuttropfenelektrode für Professional-VA- Geräte

6.5339.080

Kompletter Elektrodensatz für voltammetrische Bestimmungen von Schwermetallen. Enthält Bismuttropfenelektrode, Referenzelektrode, Glassy-Carbon-Hilfselektrode, Messgefäß, Rührer, Elektrolytlösung und weiteres Zubehör.

Lieferumfang 6.5339.080

Qt.	Order no.	Beschreibung
1 PCS	6.0346.000	Bismuth drop electrode
Bismuth drop electrode for the voltammetric determination of traces of metal ions in water samples. It can be used as a working electrode in the 884 Professional VA.		



1 PCS

6.0728.120

Ag/AgCl reference electrode

Inner reference system with electrolyte $c(\text{KCl}) = 3 \text{ mol/L}$. For use with electrolyte vessel made of plastic 6.1245.010.



1 PCS

6.1204.500

Stirrer for Professional VA instruments

Stirrer for Professional VA instruments, complete, metal-free.



1 PCS

6.1241.120

Electrode holder

Electrode holder for VA auxiliary electrodes with plug contact 5.5. mm. For use with glassy carbon rod (6.1247.000) or platinum pin (6.1247.010).



1 PCS

6.1245.010

Electrolyte vessel made of plastic

Electrolyte vessel made of plastic with ceramic diaphragm for double junction reference electrode.



1 PCS

6.1247.000

Glassy carbon rod (GC)

For electrode holder 6.1241.020 or 6.1241.120 for VA auxiliary electrode.



1 PCS

6.1415.210

Titration vessel / 10-90 mL



2 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.2308.020

Electrolyte 3 mol/L KCl (250 mL)

Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)



For 62753200 oil separator and for applications with the scTRACE Gold.

