

# **Eco Dosimat**

2.1007.0010

The Eco Dosimat is a dispensing unit with built-in magnetic stirrer and touch display for universal use.

Increase your efficiency and precision for classic laboratory tasks such as

- the precise dosing of any volume.
- manual titrations.
- the preparation of ion standards.
- controlled dosing with specified dosing rate.

A report in accordance with GLP can be stored as a PDF or printed out for all of the activities listed.

Including push-button cable for manual titration and dispensing control.


Freely selectable cylinder unit with 5, 10, 20, and 50 mL.



Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

## Scope of delivery 2.1007.0010

Qt.	Order no.	Description		
1 PCS	1.1007.0010	Eco Dosimat		
	The Eco Dosimat is a dispensing unit with built-in magnetic stirrer and touch display for universal use.			
	<ul> <li>Increase your efficiency and precision for classic laboratory tasks such as</li> <li>the precise measurement of volumes.</li> <li>of manual titration.</li> <li>the creation of ion standards.</li> <li>controlled dosing with specified dosing rate.</li> </ul>			
	A report in accordance with GLP can be stored as a PDF or printed out for all of the activities listed.			
	Including push-button cable for manual titration and dispensing control.			
	Freely selectable o	cylinder unit with 5, 10, 20, and 50 mL.		
1 PCS	6.1446.030 For the movable h	Ball stopper / B-14/(15)  nolding of buret tips in titration vessels		







1 PCS 6.1543.060 Tip / M6 thread

Tip with M6 thread. Together with the 6.1805.160 capillary and the 6.1446.030 link stopper, it makes up the entire 6.1537.010 buret tip.



1 PCS 6.1543.200 Titration tip / M6 thread

Titration tip with antidiffusion valve.



1 PCS 6.1546.040 Piston tool for Eco Dosimat/Eco Titrator

Tool for removing and mounting the piston.



1 PCS 6.1602.105 Bottle attachment / GL 45 / green

For direct dosing from reagent bottles with GL 45 thread.



1 PCS 6.1608.023 Amber glass bottle / 1000 mL / GL 45

For Exchange Units. Bottle for auxiliary solutions.





1 PCS

6.1619.010

Adsorber complete, for Exchange Unit



2 PCS

6.1805.090

FEP tubing / M6 / 31 cm

With light and kink protection



1 PCS

6.1805.100

FEP tubing / M6 / 40 cm

With light and kink protection



1 PCS

6.1819.020

Tube

For M6 threaded holes. For 6.1602.1XX bottle attachments





1 PCS

6.1903.020

Stirring bar / 16 mm

Stirring Bar with magnetic core, PTFE coated, length 16 mm



1 PCS

6.1903.030

Stirring bar / 25 mm

Stirring bar with magnetic core, PTFE covering, length 25 mm.



1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2016.070

Support rod / 400 mm



1 PCS

6.2021.020

Electrode holder

Electrode holder for 4 electrodes and 2 buret tips





1 PCS

6.2023.020

Clip for SGJ 14/15

Clip for SGJ 14/16



1 PCS

6.2107.100

Manual Dosing Controller

Handle for manual dosing and manual titration with the Eco Dosimat or Dosimat plus



1 PCS

6.2164.010

Power supply unit 100 - 240 V/24 V DC

Power supply unit for 915 KF Ti-Touch, 916 Ti-Touch, 917 Coulometer and all Eco Titrators / Eco Dosimat.



1 PCS

6.2709.070

Guiding sleeve

For fixation of buret and dosing tips in titration heads and electrode holders





## Optional accessories

Order no.	Description	
2.142.0100	Custom Q3X thermal printer	
	Compact printer with USB interface for	
	<ul> <li>900 Touch Control</li> <li>915 KF Ti-Touch</li> <li>916 Ti-Touch</li> <li>917 Coulometer</li> <li>877 / 848 Titrino plus</li> <li>865 / 876 Dosimat plus</li> <li>91X Meter (cable 6.2151.140)</li> <li>Eco Dosimat / Titrator</li> <li>862 Compact Titrosampler</li> <li>870 KF Titrino plus</li> <li>899 Coulometer</li> </ul>	
2.143.0200	Paper width 60 mm (40 characters). Including 6.2151.120 USB cable.  TX-900MH impact printer	
2.143.0200		
	<ul> <li>Compact impact printer with USB and RS232 interface for</li> <li>900 Touch Control (+)</li> <li>915 KF Ti-Touch (+)</li> <li>916 Ti-Touch (+)</li> <li>917 Coulometer (+)</li> <li>91X Meter (+)</li> <li>Eco Dosimat/Titrator (+)</li> <li>877/848 Titrino plus (additional 6.2151.100 required)</li> <li>865/876 Dosimat plus (additional 6.2151.100 required)</li> <li>862 Compact Titrosampler (additional 6.2151.100 required)</li> <li>870 KF Titrino plus (additional 6.2151.100 required)</li> <li>899 Coulometer (additional 6.2151.100 required)</li> </ul>	C C C C C C C C C C C C C C C C C C C
	(+) Cable included in scope of delivery	
6.03002.150	Cylinder unit Eco 5 mL	

Cylinder unit 5 mL for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.





### 6.03002.210 Cylinder unit Eco 10 mL

Cylinder unit 10 mL for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.



### 6.03002.220 Cylinder unit Eco 20 mL

Cylinder unit 20 mL for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.



### 6.03002.250 Cylinder unit Eco 50 mL

Cylinder unit 50 mL for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.

