



# Balance for 864 Robotic Balance Sample Processor, 230V

2.135.0206

The balance (Precisa) has a special construction for usage below the sample rack 6.2068.010 of the 864 Robotic Balance Sample Processor. The power supply is for 230V.

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Scope of delivery 2.135.0206

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

1.135.0206

Balance for 864 Robotic Balance Sample Processor,  
230V

The balance (Precisa) has a special construction for usage below the  
sample rack 6.2068.010 of the 864 Robotic Balance Sample Processor.  
The power supply is for 230V.



## Optional accessories

Order no.	Description	
6.2067.100	<p>Scale pan for sample rack 6.2068.010</p> <p>The scale pan is designed for the 120 mL PP beaker (6.1459.300) and is placed into the sample rack of the 864 Balance Sample Processor.</p>	
6.2068.010	<p>Sample rack, 20 x 120 mL, (PP)</p> <p>Sample rack for 864 Balance Sample Processor, suitable for 20 sample beakers. The following sample beakers can be used: 6.1459.300.</p> <p>Plastic: Polypropylene (PP)</p>	
6.2711.080	<p>Drip pan for 864 Robotic Balance Sample Processor</p> <p>The drip pan guards Robotic Balance Sample Processor and balance against contamination / corrosion by chemicals.</p>	