

## Polytron PT 1300 D

2.136.0100

Polytron PT 1300 D - Metrohm version

Homogenizer that can be controlled directly by OMNIS Software, tiamo™ or Touch Control.

The Polytron PT 1300 D is comprised of a control unit and a drive. The coupling system on the drive permits simple and rapid replacement of the aggregates without additional tools. Solid samples can be pulverized without difficulty. However the instrument is also highly suitable for ensuring thorough mixing in viscous samples.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

## Scope of delivery 2.136.0100

Qt. Order no. Description



Polytron PT 1300 D - Metrohm version

Homogenizer that can be controlled directly by tiamo $^{\mathsf{TM}}$  or Touch Control.



1 PCS

1.136.0110

Power supply unit for Polytron PT 1300 D



## Optional accessories

Order no.	Description
6.02008.010	Polytron holder
	For fastening the Polytron in place, suitable for 601405040.
6.02008.030	Holder for homogenizer for Eco Titrations instruments
	Holder for homogenizer for fastening the Polytron in place, suitable for 601405060 titration vessel lid

6.1462.250 Swing arm for Polytron and 786 Swing Head, left swinging

Swing arm for Polytron and 786 Swing Head, left swinging, with one buret tip for the additiona of solvent or titrant and three spray nozzles for cleaning. The swing arm allows the usage of the Polytron for sample preparation on multi row racks.



6.2047.110 Holder for Polytron

Holder to fix the Polytron at the stand rod.





6.9012.000

Dispersing aggregate for Polytron, 125 mm

Standard dispersing aggregate for Polytron 1300 D



6.9012.010

Dispersing aggregate for Polytron, 157 mm

Special dispersing aggregate (Metrohm version) with protruding knives

