



# Auto Boat Drive (AJ)

2.136.0710

The Auto Boat Drive (ABD) automatically guides the sample boats into the combustion oven. The sample can be put into the boat manually or using an appropriate sampler. This product made by Analytik Jena has been specially optimized for coupling with ion chromatographs.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

## Optional accessories

Order no.	Description
6.7302.100	Quartz boat 40 x 9 mm (AJ)
	Quartz glass boat for sample injection in Combustion Module, 10 pcs.

