

703 Ti Stand with stirrer and pump

2.703.0010

Compact titration stand with magnetic stirrer and built-in pump for addition of solvents and for aspirating off titration vessel contents.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 2.703.0010

Qt.	Order no.	Description
1 PCS	1.703.0010	703 Ti Stand with stirrer and pump
	•	stand with magnetic stirrer and built-in pump for ts and for aspirating off titration vessel contents.





1 PCS

6.1446.060

Stopper / B-14/15 / M10



1 PCS

6.1543.110

Dosing tip / M8 thread

Transfer tip with M8 thread



2 PCS

6.1602.105

Bottle attachment / GL 45 / green

For direct dosing from reagent bottles with GL 45 thread.



For Exchange Units. Bottle for auxiliary solutions.



1 PCS

6.1608.030

Round glass bottle / 1000 mL / GL 45



1 PCS

6.1609.010

Adsorber tube with tubing nipple

Lid with tubing nipple.



2 PCS

6.1801.000

PVC tubing / 4 mm / 7 mm / 1 m





1 PCS

6.1808.050

Tubing adapter olive / M8 outer

1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



1 PCS

6.1818.000

Aspiration tubing

For 681 Pump Unit and 703 Ti Stand



1 PCS

6.1819.030

Tube

For M8 threaded holes. For 6.1602.1XX bottle attachments and 6.1608.0 XX bottles





1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2023.020

Clip for SGJ 14/15

Clip for SGJ 14/16



1 PCS

6.2108.100

Connection cable for 622/649/703/728, 0.5 m

For supplying power to 622, 649, 728 Stirrers and 703 Ti Stand from control instruments and corresponding Dosimats



For tightening connectors



