



# 843 Pump Station (peristaltic) - rinse/aspirate for Sample Processors

2.843.0150

The 843 Pump Station (peristaltic) has two built-in peristaltic pumps. These can be controlled directly from the Sample Processors via pump commands. The rinse/aspirate version is provided with all the accessories needed for automatically emptying the titration beaker and rinsing the titration equipment. The accessories are intended for use with the 814 and 815 Sample Processors.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

## Scope of delivery 2.843.0150

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

1.843.0120

843 Pump Station (peristaltic)

The 843 Pump Station (peristaltic) has two built-in peristaltic pumps. These can be controlled directly via the interface using remote signals or manually by pressing a button.



---

1 PCS

6.1543.170

Aspiration tip / M8 thread

Aspiration tip for the Sample Processors.



---

2 PCS

6.1621.000

Container 10 L

As rinsing or waste container in automated systems.



---

2 PCS

6.1621.100

Cap for 10-L container

Cap for 10-L container. For easy and safe transportation of filled containers.



4 PCS

6.1801.120

PVC tubing / 4 mm / 6 mm / 2 m



---

3 PCS

6.1805.060

FEP tubing / M6 / 60 cm

With light and kink protection



---

3 PCS

6.1805.420

FEP tubing / M6 / 48 cm

With light and kink protection



---

1 PCS

6.1805.510

PTFE tubing / M8 / 60 cm

With protection against kinks.



1 PCS

6.1808.170

Distributor for Sample Processors for rinsing and aspirating

Distributor for installation on Sample Processors for rinsing and aspirating when using externally connected pumps.



2 PCS

6.1812.000

PTFE tubing 4 mm / 6 mm / 4 m



4 PCS

6.1820.050

M6 / Olive adapter

For 6.1826.100 pump tubing, to peristaltic pumps.



1 PCS

6.1826.100

Pump tubing 6.4 mm / 1 m

Replaces the pump tubing used with the 772 Pump Unit and the 843 Pump Station (peristaltic).



2 PCS

6.1828.000

PVDF connection nipple

For 6.1621.000 canister.



---

1 PCS

6.2141.300

Remote cable 843 Pump Station to Sample Processor

Remote cable for the direct connection of the 843 Pump Station to the pump connectors of the Sample Processors



---

3 PCS


6.2740.020

Spray nozzle

For the fine-spraying of the rinsing solution.



# Optional accessories

Order no.	Description	
6.2740.030	<p>Rinsing nozzle</p> <p>The rinsing solution is applied as a fine jet for the optimal removal of layers on the electrodes and titration accessories. For 6.5613.000 ECO titration equipment.</p>	

---