

## 849 Level Control

2.849.0010

Additional equipment for Sample Processors for monitoring the filling level of rinsing or waste canisters via Remote. Prevents pumps from running dry and/or canisters from overflowing and is suitable for use with aqueous solutions, solvents and suspensions.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

## Scope of delivery 2.849.0010

Qt. Order no. Description



Additional equipment for Metrohm instruments. For monitoring the filling level of bottles for eluent, auxiliary or waste solutions via remote control. The 849 Level Control prevents pumps from running dry and/or canisters from overflowing and is suitable for use with aqueous solutions, solvents and suspensions.



The metal-free sensor works in solutions with very low conductivity. Therefore it may be used in ion chromatography.

1 PCS 6.1113.000 Sensor for 849 Level Control

Sensor for 849 Level Control



1 PCS 6.2151.060 Cable for Level Control and Eluent Production Module

Connecting cable for sensor and Level Control or Eluent Production Module.

