



## 889 IC Sample Center – cool

2.889.0020

The 889 IC Sample Center – cool is the appropriate automation solution when you have only a very small amount of sample. In comparison with the 889 IC Sample Center, it possesses in addition a cooling function and is thus the ideal sample changer for biochemically relevant or thermally unstable samples.

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.889.0020

| Qt. | Order no. | Description |
|-----|-----------|-------------|
|-----|-----------|-------------|

1 PCS

1.889.0020

889 IC Sample Center – cool

The 889 IC Sample Center – cool is the appropriate automation solution when you have only a very small amount of sample. In comparison with the 889 IC Sample Center, it possesses in addition a cooling function and is thus the ideal sample changer for biochemically relevant or thermally unstable samples.



1 PCS

6.1608.110

Wash bottle 250 mL to IC Sample Center

Bottle for washing solution in the IC Sample Center.



2 PCS

6.1807.010

Y-connector for tubing 6-9 mm i.d.

Connector for waste tubing



2 PCS

6.1816.070

Silicone tubing 8/12 mm, 1 m

For the IC Sample Center



1 PCS

6.1831.010

PEEK capillary 0.25 mm i.d. / 3 m

For all IC components.



---

2 PCS

6.2041.200

Sample rack 48 x 11 mm to IC Sample Center

Sample rack for the IC Sample Center with 48 sample positions. For sample vials 6.2743.1xx.



---

1 PCS

6.2151.020

Cable USB A - USB B 1.8 m

USB connecting cable



---

1 PCS

6.2743.107

PP Vials 0.7 mL to IC Sample Center, 100 pieces

100 sample vials.



1 PCS

6.2743.127

Septum to 6.2743.1XX Si/PTFE

100 pieces



---

1 PCS

6.2744.010

Pressure screw 5x

With UNF 10/32 connection. For the connection of PEEK capillaries



---

1 PCS



6.2816.200

Air needle 80 mm, black

To IC Sample Center



## Optional accessories

| Order no.  | Description                                                                                                         |                                                                                       |
|------------|---------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| 6.1448.050 | Aluminum septum caps / 1000 pieces<br>Including Al sealing.                                                         |    |
| 6.1448.067 | Screw cap with septum, 100 pcs.<br><br>Screw cap with septum, suitable for sample glass 8 mL 6.2420.007, 100 pieces |    |
| 6.1562.170 | Tefzel transfer tubing 1000 µL<br><br>To IC Sample Center.                                                          |  |
| 6.1829.060 | Capillary set for syringe unit<br><br>For the IC Sample Center                                                      |  |

6.2041.210 Sample rack, 12 x 22 mm, for IC Sample Center

Sample rack for 12 vials.



---

6.2419.000 Sample glass, 6 mL, 1000 pieces

Sample glass for up to 6 mL sample, suitable for 6.1448.050 septum, 1,000 pieces



---

6.2420.007 Sample glass with thread, 8 mL, 100 pieces

Sample glass with thread, 8 mL, suitable for 6.1448.067 screw cap, 100 pieces



---

6.2621.110 Crimping tongs

Crimping tongs for sealing 6 mL sample vials (6.2419.000, 6.2419.007).



6.2743.100 PP Vials 0.7 mL to IC Sample Center, 1000 pieces  
1000 sample vials.



---

6.2743.110 PP Vials 0.3 mL to IC Sample Center, 1000 pieces  
1000 sample vials.




---

6.2743.117 PP Vials 0.3 mL to IC Sample Center, 100 pieces  
100 sample vials.



---

6.2743.120 Septum to 6.2743.1XX Si/PTFE  
1000 pieces



6.2816.210 Needle with capillary, 15  $\mu$ L  
For the IC Sample Center



---

6.5333.000 IQ/OQ Kit for IC  
The IQ/OQ Kit contains all parts and standard solutions required for IQ /OQ in ion chromatography.



---

6.7303.100 Glass vial, 1.5 mL, 100 pieces

For use with the following sample changers: Liquid Autosampler CIC (TEI), MMS 5000 Autosampler (AJ) and 889 Sample Center.



---

6.7303.110 Septum for 67303100, 100 pieces

For use with the 67303100 Glass vial on the following sample changers: Liquid Autosampler CIC (TEI), MMS 5000 Autosampler (AJ), and 889 Sample Center.

