



## Micro Au Titrode (installation length 30.8 cm)

6.00435.120

Combined gold ring electrode with a pH glass membrane as reference electrode, installation length = 30.8 cm, narrowed to a diameter of 6.4 mm.

This maintenance-free electrode is well suited for redox titrations when the pH value remains constant, e.g.:

- determination of the chemical oxygen demand (COD, ferrometry)
- determination of penicillin and ampicillin
- for titrations using  $\text{Hg}(\text{NO}_3)_2$
- for titrations in presence of chromium or iron



This electrode is stored in deion. water.

Below, the accessories are grouped into Scope of delivery and Optional accessories.


Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.00435.120

Qt.	Order no.	Description
1 PCS	6.1236.050	Sleeve with SGJ 14/12 mm Sleeve with SGJ 14/12 mm, polyethylene.
		
1 PCS	6.2008.040	Storage vessel Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.
		

# Optional accessories

Order no.	Description	
6.2104.020	Electrode cable / 1 m / F  For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	

---