

## Metrosep C Trap 1 - 30/4.0

6.01015.030

Short, high-capacity cation column, which is used to purify the eluent flow. Even reagents of the highest quality level, e.g. "Ultrapur" or "puriss", can still contain minimum amounts of cation contamination. These are reliably held back by the Metrosep C Trap 1 - 30/4.0.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

## Scope of delivery 6.01015.030

Qt.	Order no.	Description	
2 PCS	6.2744.060	Threaded stopper	
	For UNF 10/32. columns.	Stopper for IC, e.g. for the sealing of	

