



Titration head for KFT with reinforced tips (P&P)

6.01403.050

Titration head for Karl Fischer titration in Pick&Place modules of the OMNIS Sample Robot Pick&Place. In addition to 2 reinforced antidiffusion tips, the titration head offers an opening for using an homogenizer (Polytron 1300D) or rod stirrer.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.01403.050

Qt.	Order no.	Description
2 PCS	6.1808.000	Connection bushing M6



2 PCS



6.2726.090

Antidiffusion valve

For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140,
6.1543.150, 6.1543.160, 6.1543.200



Optional accessories

Order no.	Description	
6.01201.000	Rod stirrer SGJ sleeve for 6.01403.050	
	Rod stirrer SGJ sleeve for when using a rod stirrer instead of an homogenizer in the 6.01403.050 titration head	
<hr/>		
6.9012.010	Dispersing aggregate for Polytron, 157 mm	
	Special dispersing aggregate (Metrohm version) with protruding knives	
<hr/>		