



PRP-X100 guard column cartridge/4.0

6.1005.020

For protecting Hamilton PRP-X100 columns. The cartridge efficiently removes contamination in the form of particles due, for example, to the growth of bacteria or algae.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description
6.02821.000	Holder for the PRP-X100 Guard cartridge
	Holder for the Hamilton PRP-X 100 Guard cartridge for connection to the analytical separation column.

