

## Metrosep BO3 Trap 1 - 100/4.0

6.1015.200

Trap column for the removal of borate contaminants from the eluent. The Metrosep  $BO_3$  Trap 1 - 100/4.0 is used mainly in carbohydrate analysis with hydroxide eluents.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

## Scope of delivery 6.1015.200

Qt.	Order no.	Description
2 PCS	6.2744.060	Threaded stopper
	For UNF 10/32. columns.	Stopper for IC, e.g. for the sealing of

